SmartDrive

Power your push, your way

June 2021 enhancements sales team training





Curtis A. Merring, MOT, OTR

Director, Portfolio Marketing

Permobil Americas Region curtis.merring@permobil.com

Poll Question 1





SmartDrive overview & the why



Proactive vs. Reactive Approach

Current Reactive Approach New Proactive Approach



Proactive vs. Reactive Approach

Evidenced based practice - problems

> 50% of manual wheelchair (MWC) users will develop pain and injury

(Finley and Rodgers 2004)

Pain and injury can lead to significant decline in independence and quality of life for MWC users

(Gutierrez et al. 2007)

Only 55% of patients (USA) received recommended care (McGlynn, 2003)



Proactive vs. Reactive Approach

Evidenced based practice - solutions

Power assist wheels decrease heart rate and perceived exertion during MWC mobility

(Algood, 2005)

Power assist eases energy cost/perceptual responses to propulsion for person w/ shoulder pain and SCI

(Nash, 2008)

Power assist wheels increase distance traveled compared to conventional propulsion

(Levy, 2010)

MWC users push 80% less with SmartDrive vs. without measured by the Smart Evaluation Application

(Smart Evaluation Application, 2021)






SmartDrive overview & the why

Manual wheelchair users have UE challenges beyond propulsion





Poll Question 2



Value-Based Health Care Benefits

PATIENTS

Lower Costs & better outcomes

PROVIDERS

Higher Patient
Satisfaction
Rates &
Better Care
Efficiencies

PAYERS

Stronger Cost Controls & Reduced Risks

SUPPLIERS

Alignment of Prices with Patient Outcomes

SOCIETY


Reduced Healthcare Spending & Better Overall Health

NEJM Catalyst (catalyst.nejm.org) © Massachusetts Medical Society

Risk Factors:

Patient age,

years w/c dependence



Percentage with RCT

Gelmman, 1998 & Sie, 1992)



Evidenced based practice - problems

Rotator cuff tears are present in 49% of wheelchair users of which 70% were full thickness

(Akbar, 2011)

Carpal tunnel syndrome is a common side effect of long-term wheelchair use, & severity is associated with duration of use/age (Asheghan, 2015)



Direct costs

- Operation
 - > \$6,367 Medicare
 - > \$6,904 other insurers
- Rehab
 - Pre-op visits
 - Post-op visits (avg. 15)



Spectrumortho.com



Indrect costs

- No propulsion x 6 weeks
- NWB transfers x 4-6 months
- A w/ dressing, transfers, position changes
- PWC rental & transportation
- Extended time off work
- Caregiver fees



Gg120830570 GoGraph.com

SmartDrive (E0896):

Medicare Reimbursement ~ \$5750

- no configuration changes needed
- no costly remodeling needed to access same environment
- Removable for ease of transporting



TiLite ULMWC (K0005):

Medicare Reimbursement ~ \$2,100

- more durable
- more configurable
- lighter weight materials



Power your push, your way with SmartDrive

SmartDrive overview & the why

"It is the primary truth of complex rehab that no single solution works well for all wheelchair users. There's simply too much diversity in diagnoses, progression, function, goals and environments for one product to work well for everyone."

"Permobil's SmartDrive accomplishes this by giving manual wheelchair users several ways to access and use its world-expanding technology."

Laurie Watanabe – Mobility Management

https://mobilitymgmt.com/articles/2019/12/20/smartdrive-pushtracker-e2-switch-control.aspx?s=emob_060120&admgarea=emob&oly_enc_id=



SmartDrive overview & the why

Power your push, your way

	Stamina	Performance	Expert
Products available .	70		
Best for	Hemi-propeller, pediatric, geriatric, or caregivers	Pediatric/active users (smaller, simpler wearable)	Adult/active users (advanced wearable features)
On-demand power assist	0	0	0
Energy/shoulder preservation	0	0	0
Easy on/off buttons	0	0	0
Hand gesture activation		0	0
Advanced Bluetooth connection			0
Easy to read, touch screen settings control			0
Performance & maintenance notifications			0









SwitchControl

Is now included in with every SmartDrive order at no charge







SwitchControl programming

Separate setting than PushTracker/E2				
Set your speed	1-5.5mph/1.6-8.9kph	Separate setting than PushTracker /E2		
Set your mode	Momentary : Press to Go	Latched : Hold & Set		

Mobility considerations for SwitchControl				
Propulsion Technique	Mode options	Specific functional activity		
Bilateral Upper Extremity	Momentary alone	Ramps/Inclines		
Bilateral Lower Extremity	Latched alone	Indoor Mobility		
Hemi-Propulsion	Momentary with PushTracker/E2	Crowded Spaces		
Dependent/Caregiver operated	Latched with PushTracker/E2	Short Distances		
		Quick/Short Surface Change		
		Low/No PushTracker Battery		



SwitchControl (only) Set Up



PushTracker App



Original PushTracker & SwitchControl Set Up



PushTracker App





SmartDrive overview & the why

Power your push, your way


Set up in minutes
No phone or wi-fi required



Poll Question 3



SmartDrive enhancements, June 2021









Order form (updated)

SmartDrive enhancements, June 2021

permobil

SmartDrive Order & Request for Quote Form

Toll-Free: 800.736.0925 • Email: sales.smartdrive@permobil.com • Fax: 888.411.9027 Request for Quote Place an Order - P.O. #: Ouote #: (please complete entire form) (if you do not already have a quote, please complete entire form) **Provider Information** Shipping Ship to provider address on file Provider: Ship to alternate address: Provider account #: Alternate ship Provider location to recipient: (if applicable): Order contact name: Street: ATP / Therapist: City: State: E-mail: Mobility Products (all prices are in US dollars/USD) ■ SmartDrive MX2+ Power Assist System (HCPCS Code: E0986) MSRP: \$7131.00 (nowincludes standard SwitchControl buttons and your choice of PushTracker at no additional charge) Client/Mark for: PushTracker Options Included with Power Assist System at no charge Wheelchair make: Updated! PushTracker E2, MSRP: \$225.00, available in Wheelchair model: Shadow black · set up in minutes · no phone/Wi-Fi required Rear wheel size: · SmartDrive apps pre-installed 22" / 501 mm 24" / 540 mm 25" / 559 mm 26" / 590 mm Original PushTracker, MSRP: \$200,00. available in black with three band sizes provided: Is this an existing wheelchair or a new one being ordered at the Small, Medium and Large same time as this SmartDrive? Existing No PushTracker required Frame type: Extra Accessories/Parts (check box(es) if desired): Rigid (if the wheelchair is a Motion Composites New! Universal padded handle, MSRP: \$20.00 APEX, please provide the rear seat height): Carrying bag, MSRP: \$78.00 Folding (MUST provide the chair/frame width): USB power bank (2600 mAh), MSRP: \$20.00 Special frame, One-Arm drive, etc. Roller replacement kit, MSRP: \$205.00 (please specify): Order SwitchControl Button Options notes: SwitchControl buttons (two), MSRP: \$220.00 (standard option included at no charge) SwitchControl button (one) with mono jack option,


Send completed form via fax: 888.411.9027 or email: sales.smartdrive@permobil.com

Figure shown is for illustration purposes only and actual product may vary.

MSRP: \$370.00 (upgrade option with additional \$150.00 charge)

Brochure (updated)

SmartDrive enhancements, June 2021



Enhancements supplement

SmartDrive enhancements, June 2021





permobil

50% more

SmartDrive's updated electronics provide 50% more power to give your SmartDrive more power to handle more load over time.

In-field repair videos

SmartDrive enhancements, June 2021

SmartDrive in-field repair videos

- 1. Removing handle (for all applicable videos)
- 2. Circuit Board
- 3. Battery
- 4. Adapter bar set up with retrofit
- 5. Adapter bar set up without retrofit
- 6. Power switch
- 7. Charger receptacle
- 8. Receptacle cap

SmartDrive enhancements, June 2021







References

- Algood, S. D., Cooper, R. A., Fitzgerald, S. G., Cooper, R., & Boninger, M. L. (2005). Effect of a pushrim-activated power-assist wheelchair on the functional capabilities of persons with tetraplegia. Archives of Physical Medicine and Rehabilitation, 86(3), 380-386.
- Akbar, M., Brunner, M., Balean, G., Grieser, T., Bruckner, T., Loew, M., & Raiss, P. (2011). A cross-sectional study of demographic and morphologic features of rotator cuff disease in paraplegic patients. *Journal of Shoulder and Elbow Surgery*, 20(7), 1108–1113. https://doi.org/10.1016/j.jse.2011.03.021
- Asheghan, M., Hollisaz, M. T., Taheri, T., Kazemi, H., & Aghda, A. K. (2015). The prevalence of carpal tunnel syndrome among long-term manual wheelchair users with spinal cord injury: A cross-sectional study. *The Journal of Spinal Cord Medicine*, *39*(3), 265–271. https://doi.org/10.1179/2045772315y.0000000033
- Finley, M.A., Rodgers, M.M., 2004. Prevalence and identification of shoulder pathology in athletic and nonathletic wheelchair users with shoulder pain: a pilot study. J. Rehabil. Res. Dev. 41, 395–402
- Gutierrez, D.D., Thompson, L., Kemp, B., Mulroy, S.J., 2007. The relationship of shoulder pain intensity to quality of life, physical activity, and community participation in persons with paraplegia. J. Spinal Cord Med. 30, 251–255.
- Kim, S. S., Her, J. G., & Ko, T. S. (2015). Effect of different hand positions on trunk and shoulder kinematics and reaction forces in sitting pivot transfer. *Journal of physical therapy science*, 27(7), 2307–2311. https://doi.org/10.1589/jpts.27.2307
- H. Gellman, I. Sie, R.L. WatersLate complications of weight-bearing upper extremity in the paraplegic patient. Clin Orthop Relat Res, 233 (1998), pp. 132-135
- I.H. Sie, R.L. Waters, R.H. Adkins, H. Gellman. Upper extremity pain in the post rehabilitation spinal cord injured patient. Arch Phys Med Rehabil, 73 (1992), pp. 44-48



References

- Levy, C. E., Buman, M. P., Chow, J. W., Tillman, M. D., Fournier, K. A., & Giacobbi Jr, P. (2010). Use of Power Assist-Wheels Results in Increased Distance Traveled Compared to Conventional Manual Wheeling. American Journal of Physical Medicine & Rehabilitation/Association of Academic Physiatrists, 89(8), 625.
- McGlynn, E. A., Asch, S. M., Adams, J., Keesey, J., Hicks, J., DeCristofaro, A., & Kerr, E. A. (2003). Quality of Health Care Delivered to Adults in the United States. New England Journal of Medicine, 349(19), 1866-1868. doi:10.1056/nejm200311063491916.
- Narvy, S. et al (2016). Direct cost analysis of outpatient arthroscopic rotator cuff repair in medicare and non-medicare populations. Orthop J Sports Med. 4(10): 2325967116668829.
- Narvy, S. J., Didinger, T. C., Lehoang, D., Vangsness Jr, C. T., Tibone, J. E., Hatch III, G. F. R., ... & Gamradt, S. C. (2016). Direct cost analysis of outpatient arthroscopic rotator cuff repair in Medicare and non-Medicare populations. *Orthopaedic journal of sports medicine*, 4(10), 2325967116668829.
- Nash, M. S., Koppens, D., van Haaren, M., Sherman, A. L., Lippiatt, J. P., & Lewis, J. E. (2008). Power-Assisted Wheels Ease Energy Costs and Perceptual Responses to Wheelchair Propulsion in Persons With Shoulder Pain and Spinal Cord Injury. *Archives of Physical Medicine and Rehabilitation*, 89(11), 2080–2085. https://doi.org/10.1016/j.apmr.2008.05.018
- Paralyzed Veterans of America Consortium for Spinal Cord Medicine. (2005). Preservation of upper limb function following spinal cord injury: a clinical practice guideline for health-care professionals. The journal of spinal cord medicine, 28(5), 434.

Permobil webinar series

Up Next!

Permobil White Paper: A systematic review of the evidence for power standing wheelchairs

Tuesday June 22, 2021-**Same Webinar, Three Global Times!**

- 7am Central
- 1 pm Central
- 6 pm Central



Questions? education@permobil.com







