



TiLite X Series - Order Form



Dealer Name*: _____ Account #*: _____

City, Post Code, Country: _____

Phone: _____ Email: _____ P.O. Number _____

Clinician: _____ Client Reference*: _____ Client Weight: _____

Rep Name: _____

Special considerations that need to be addressed here (e.g., diagnosis): *Required field

Frame Style

AXFS1	Aero X Aluminum Folding Manual Wheelchair Full adjustability, 265lb (120kg) weight limit.	\$
G2FS1	2GX Titanium Folding Manual Wheelchair Full adjustability, 265lb (120kg) weight limit. Light weight and high strength. Not available with frame transit tie-down option.	\$

Frame Options

XFS1	Swing Away Front Frame	Standard
XFS2	Fixed Front Frame	No charge
XFS3	Hemi Height Frame Available with a front seat height of 13.5"-20". Not with depth adjustable back/frame, 12"-13" seat depth or heavy duty frame.	No charge
XHDR1	Heavy Duty Frame 350lb (158kg) weight limit. Includes tension adjustable by straps back upholstery. Available with 18"-22" seat widths and 16"-20" seat depths. Requires LiteSpeed billet aluminum front wheels and 0° camber. Not available with 90° front end or single sided forks.	\$
XPWR1	Power Adaptable Frame Requires LiteSpeed billet aluminum front wheels and 0° camber. Not available with single-sided forks.	\$
XRRF1	Reinforced Frame	\$

Transit Options

XTTDO1	Aluminum Frame Transit Tie-Down Consists of four tie-down brackets. Not with 2GX, heavy duty frame, hemi height frame, swing away tubular armrests or swing away flip back tubular armrests. Requires a folding stabilizer bar unless a depth adjustable backrest is selected. Requires minimum seat depth of 14". Please refer to the Aero X Transit Tie-Down Option Brochure for important information.	\$
XWTB	Transport Bracket Option Consists of four tie-down brackets. Unoccupied use only. Not available when swing away frame and heavy duty frame are both selected. Not with swing away tubular armrests or swing away flip back tubular armrests. Requires a minimum seat depth of 14", or 15" with hemi height frame.	\$

Frame Finishes

XFF1	Satin Finish Available with titanium frames only.	No charge																								
	XSCBT1 Buffing Pads for Satin Frame	\$																								
XFF	Powder Coated Frame Select below.	No charge																								
Select color:	<table border="0"> <tr> <td>XFF8 High Gloss Black</td> <td>XFF10 Red Pepper</td> <td>XFF29 Neon Pink</td> <td>XFF14 Ocean Blue Metallic</td> </tr> <tr> <td>XFF15 Tangerine Metallic</td> <td>XFF30 Acid Green Metallic</td> <td>XFF17 Midnight Blue Pearl</td> <td>XFF18 White Pearl</td> </tr> <tr> <td>XFF19 Black Lava</td> <td>XFF27 Titanium Grey</td> <td>XFF59 Pomegranate Red</td> <td>XFF62 Blue Mist</td> </tr> <tr> <td>XFF57 Indigo Purple</td> <td>XFF40 Hyper Blue</td> <td>XFF41 Flame Orange</td> <td>XFF42 OD Green</td> </tr> <tr> <td>XFF60 Silver Streak</td> <td>XFF63 Lavendar Dusk</td> <td>XFF58 Black Galaxy</td> <td>XFF61 Seafoam Green</td> </tr> <tr> <td>XFF64 Daffodil Yellow</td> <td></td> <td></td> <td></td> </tr> </table>	XFF8 High Gloss Black	XFF10 Red Pepper	XFF29 Neon Pink	XFF14 Ocean Blue Metallic	XFF15 Tangerine Metallic	XFF30 Acid Green Metallic	XFF17 Midnight Blue Pearl	XFF18 White Pearl	XFF19 Black Lava	XFF27 Titanium Grey	XFF59 Pomegranate Red	XFF62 Blue Mist	XFF57 Indigo Purple	XFF40 Hyper Blue	XFF41 Flame Orange	XFF42 OD Green	XFF60 Silver Streak	XFF63 Lavendar Dusk	XFF58 Black Galaxy	XFF61 Seafoam Green	XFF64 Daffodil Yellow				
XFF8 High Gloss Black	XFF10 Red Pepper	XFF29 Neon Pink	XFF14 Ocean Blue Metallic																							
XFF15 Tangerine Metallic	XFF30 Acid Green Metallic	XFF17 Midnight Blue Pearl	XFF18 White Pearl																							
XFF19 Black Lava	XFF27 Titanium Grey	XFF59 Pomegranate Red	XFF62 Blue Mist																							
XFF57 Indigo Purple	XFF40 Hyper Blue	XFF41 Flame Orange	XFF42 OD Green																							
XFF60 Silver Streak	XFF63 Lavendar Dusk	XFF58 Black Galaxy	XFF61 Seafoam Green																							
XFF64 Daffodil Yellow																										
XFF3	Tattoos Tattoos on front frame tubes. Finished in satin. Tattoos will be trimmed and adjusted to fit your frame.	\$																								
	XFF65 Liquid XFF66 Bloom XFF67 Rebel XFF68 Vector																									



Liquid



Bloom



Rebel



Vector

SBRS20	Comfort Unite 2 Back <i>Replaces Comfort Acta-Back Deep w/3" Wing. Inc 2-point quick release Compass hardware. Select shell height and cover.</i>	\$
	12" 14" 16" 18"	No charge
	20"	No charge
	Comfort-Tek Cover Diamond Spacer Cover	
	Privacy Shield	\$
SBRS22	Comfort Embrace Back <i>Replaces Comfort Acta-Embrace Back. Inc 2-point quick release Compass hardware. Select shell height.</i>	\$
	10" 12" 14" 16"	
	Comfort-Tek Cover Diamond Spacer Cover	
	Privacy Shield	\$
SBRS23	Comfort Relief Back <i>Replaces Comfort Acta-Relief Back. Inc 2-point quick release Compass hardware. Select shell height.</i>	\$
	10" 12" 14" 16"	
	Comfort-Tek Cover Diamond Spacer Cover	
	Privacy Shield	\$

Backrest

Back Type

XSB13	Fixed Angle, Height Adjustable Back	Standard
XLB1	Folding Stabilizer Bar <i>Required with transit tie-down option unless depth adjustable backrest is selected. Includes tension adjustable by strap back upholstery.</i>	\$
XSB14	Fixed Angle, Height Adjustable Back w/8° Bend <i>Not available with heavy duty frame.</i>	No charge
XSB18	Folding, Angle and Height Adjustable Back <i>Not available with heavy duty frame.</i>	\$
XSB16	Depth Adjustable Back Frame <i>Not available with hemi height frame, heavy duty frame, 19"-20" seat depths, swing away tubular or swing away flip back armrests. Seat depths 12"-17" allow 3" of rearward adjustment. Seat depth of 18" allows 2" of rearward adjustment. Requires tension adjustable seat upholstery.</i>	\$

Push Handles

XSB6	Integrated Push Handles	No charge
XSB7	Integrated Fold-Down Push Handles <i>Requires tension adjustable by straps back upholstery. Not available with heavy duty frame.</i>	\$
XSB9	Bolt-On Height Adjustable Push Handles <i>Requires tension adjustable by straps back upholstery. Not available with heavy duty frame, folding angle and height adjustable back, swing away tubular or swing away flip back tubular armrests.</i>	\$

Fixed Angle Adjustable Back Height 'E' | 11"-15" back height not with heavy duty frame, swing away tubular armrests 12"-15", bolt-on height adjustable push handles or folding stabilizer bar. Set at 1" increments.

Push Handles	None	11"-15"	16"-20"	Set at:
	Integrated	11"-15"	16"-20"	Set at:
	Fold-Down	11"-15"	16"-20"	Set at:
	Bolt-On Height Adj	N/A	16"-20"	Set at:
		Short	Tall	Specify
Fixed Angle Adj Back Height				

Folding Angle Adjustable Back Height 'E' | Not with folding stabilizer bar or heavy duty frame. Set at 1" increments.

Push Handles	None	11"-13"	13.5"-16.5"	16.5"-19.5"	Set at:
	Integrated	12.5"-14.5"	13.5"-16.5"	17"-19.5"	Set at:
	Fold-Down	15"-17"	16"-19"	18.5"-20.5"	Set at:
		Short	Medium	Tall	Specify
Folding Angle Adj Back Height					

Front Rigging

Front Angle 'H' | *Front angle is measured from inside of front frame tube to top of bottom frame tube.*

XFA3	Swing Away Front End	85° not with 6" front wheels. 90° requires aluminum footplate and not with heavy duty frame or 5"-6" front wheels.											No charge
		70° 80° 85° 90°	XFA0	Omit Hangers									
XFA6	Fixed Front End	85° not with 6" front wheels.											Standard
		70° 80° 85°											

Select one seat to footrest length and one footrest type.

Standard Footrests | *For front angles 70°-85°. Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than front seat height. 70° fixed front has a minimum seat to footrest of 14.5".*

Seat to Footrest 'E':	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	
	20"												
XFTR5	Flip-Up Footrest												Standard
XFTR9	Depth Adjustable Flip-Up Footrest												\$
XFTR10	Depth and Angle Adjustable Flip-Up Footrest	<i>Available for 15"-20" seat to footrest only.</i>											\$
XFTR11	Angle Adjustable Flip-Up Footrest	<i>Available for 15"-20" seat to footrest only.</i>											\$
XHLP1	Heel Loops	<i>Not with seat width ≤ 13" or one piece footrests.</i>											\$

High Mount Footrests | *Not with luggage carriers or impact guards.*

Seat to Footrest 'E':	7"	7.5"	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"	12"	12.5"	
	13"	13.5"	14"	14.5"									
XFTR20	High Mount Flip-Up												\$
XFTR15	Angle Adjustable High Mount Flip-Up Footrest												\$
XFTR16	Depth and Angle Adjustable High Mount Flip-Up Footrest												\$
XHLP1	Heel Loops	<i>Not with seat width ≤ 13" or one piece footrests.</i>											\$

90° Footrests

Seat to Footrest 'E':	7"	8"	9"	10"	11"	12"	13"						
XFTR5	Flip-Up Footrest												Standard
XFTR9	Depth Adjustable Flip-Up Footrest												\$
XFTR10	Depth and Angle Adjustable Flip-Up Footrest	<i>Available for 15"-20" seat to footrest only.</i>											\$
XHLP1	Heel Loops	<i>Not with seat width ≤ 13" or one piece footrests.</i>											\$

Elevating Legrests

Seat to Footrest 'E':	14"	15"	16"	17"	18"	19"	20"	21"					
XFA5	Depth Adjustable Flip-Up Footrest, Single												\$
		Left	Right										
XFTR17	Depth Adjustable Flip-Up Footrest, Pair												\$
XHLP1	Heel Loops	<i>Not with seat width ≤ 13" or one piece footrests.</i>											\$

Camber Options

Camber 'K'

XCBR8:	0°	2°	4°	6°									Standard
--------	----	----	----	----	--	--	--	--	--	--	--	--	----------

Front Wheels

See front seat height ranges chart for front seat height 'C' limitations.

XFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire		\$
	Silver Hub Black Hub		
XFW14	4"x0.75" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame or heavy duty frame.</i>		No charge
XFW3	4"x0.75" Light-Up Micro Caster <i>Not available with power adaptable frame or heavy duty frame.</i>		\$
XFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power assist, power add-on system or heavy duty frame.</i>		No charge
XFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire		\$
	Silver Hub Black Hub		
XFW6	5"x1" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame, heavy duty frame or with 90° front angle.</i>		No charge
XFW28	5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power adaptable frame, heavy duty frame or with 90° front angle.</i>		No charge
XFW30	5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 90° front angle.</i>		\$
	Silver Hub Black Hub		
XFW9	6"x1" Plastic Wheel w/Poly Tire <i>Not available with power adaptable frame, heavy duty frame or 85°-90° front angle.</i>		No charge
XFW32	6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not available with power adaptable frame, heavy duty frame or 85°-90° front angle.</i>		No charge
XFW33	6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not available with 85°-90° front angle.</i>		\$
	Silver Hub Black Hub		
XFW16	8"x1" Plastic Wheel w/Poly Tire <i>Not available with power assist frame, power add-on frame, heavy duty frame or 12"-14" seat depths. Only with 70° front angle.</i>		No charge

Front Forks

XFK1	TiLite Standard Fork		No charge
XGL1	Glide Suspension Fork <i>Black anodize only. Not available with 8" front wheels or heavy duty frame. Select rider weight.</i>		\$
	125lbs (57kg) 125-180lbs (57-82kg) 180-265lbs (82-120kg)		

Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

XOAD1	One Arm Drive System <i>Available with 24" mag wheels and standard aluminum handrim only. Rear seat height, camber or axle restrictions may apply.</i>		\$
	Install on Left Install on Right		
XRW14	TiLite SHADOW Rear Wheel <i>Select size. 25" and 26" wheels require 15" minimum seat depth.</i>		Standard
	22" (501) 24" (540) 25" (559) 26" (590)		
XRW7	Spinergy SPOX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth.</i>		\$
	22" (501) 24" (540) 25" (559) 26" (590)		
	Black Red Blue White Green		
XRW8	Spinergy LX Rear Wheel <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth. Only available with heavy duty frame when 0° or 2° camber is selected.</i>		\$
	22" (501) 24" (540) 25" (559) 26" (590)		
	Black Red Blue White Orange Pink Purple Green		
ZRW9	Spinergy Design Lab LX Rear Wheel <i>Black spokes only. Select size and anodizing color. 25" and 26" wheels require 15" minimum seat depth.</i>		\$
	22" (501) 24" (540) 25" (559) 26" (590)		
	Gold Green Blue Red Purple Orange		
XRW16	Spinergy Carbon Blade Rear Wheel <i>Select size. Not available with heavy duty frame. 25" wheels require 15" minimum seat depth.</i>		\$
	24" (540) 25" (559)		
XRW17	X-Core Mag Rear Wheel <i>Select size. Not available with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims. Requires aluminum silver anodized handrims.</i>		\$
	22" (501) 24" (540)		

Axles

XAXL1	Stainless Quick Release Axle		Standard
XAXL2	Titanium Quick Release Axle		\$
XAXL3	Quad Release Axle <i>Not available with Mag or 26" (590) Spinergy Spox or LX rear wheels.</i>		\$

Pneumatic Rear Tires

XRTR1	Treaded Tire 22"(501), 24"(540), 26"(590) <i>Not available with 25" rear wheels.</i>	Standard
XRTR2	Airless Inserts <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$
XRTR31	MBL SpeedLite Tire 22"(501), 24"(540), 25"(559), 26"(590)	No charge
XRTR30	Primo Sentinel Tire 24"(540), 25"(559), 26"(590)	\$
XRTR8	Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559), 26"(590)	\$
XRTR32	Primo C-764 Knobby Tire 24"(540)	\$

Solid Rear Tires

XRTR22	Primo Orion Tire 22"(501), 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$
XRTR23	Primo Active Tire 24"(540), 25"(559) <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$

Handrims

XRIM1	Aluminum Handrim, Silver Anodized	Standard
XRIM21	Aluminum Handrim, Black Anodized	No charge
XRIM12	Titanium Handrim, Satin Finish <i>Not with 22"(501) rear wheels or 24"(540) knobby tires.</i>	\$
XRIM20	Q-Grip Handrim <i>Exclusive high-friction, anti-slip coating.</i>	\$
XRIM5	Plastic Coated Handrim, Black	\$
XRIM10	The Surge Handrim <i>Large profile. All-In-One oval with Gription strip.</i>	\$
XRIM11	The Surge LT Handrim <i>Small profile. All-In-One oval with Gription strip.</i>	\$
XRIM14	The Natural-Fit Handrim, Standard Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
XRIM15	The Natural-Fit Handrim, Super Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
XRIM16	The Natural-Fit Handrim, Without Thumb Grip <i>Original, large profile.</i>	\$
XRIM17	The Natural-Fit LT Handrim, Standard Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
XRIM18	The Natural-Fit LT Handrim, Super Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
XRIM19	The Natural-Fit LT Handrim, Without Thumb Grip <i>Reduced weight, small profile.</i>	\$
XRIM22	Omit Handrims	No charge

Handrim Tab Length

XRIM7	Long Tabs	Standard
XRIM8	Short Tabs <i>Not available with Mag rear wheels, The Natural-Fit handrims with thumb grip, or Natural-Fit handrims on Shadow wheels.</i>	\$

Wheel Locks

XWLK7	Compact Composite Wheel Lock <i>Depending on configuration there may be compatibility issues with 12"-13" seat depths. Uni-Lock will be substituted at no charge.</i>	No charge
XWLK8	Compact Composite Pull to Lock <i>Depending on configuration there may be compatibility issues with 12"-13" seat depths. Uni-Lock will be substituted at no charge.</i>	No charge
XWLK10	Extension Handles for Compact Composite Push/Pull to Lock	\$
XWLK11	MBL Composite Scissor Lock <i>Not available with swing away frame.</i>	\$
XWLK12	MBL Performance Aluminum Scissor Lock <i>Not available with swing away frame.</i>	\$
XWLK3	Uni-Lock	\$
	Push to Lock Pull to Lock	
XWLK4	Extension Handles for Uni-Lock	\$
XWLK15	Grade Aids, Push to Lock <i>Not available with solid rear tires.</i>	\$

Back Upholstery

XBUP1	Padded Nylon Back Upholstery <i>Black only. Not available with heavy duty frame, folding back, folding stabilizer bar, depth adjustable frame, fold-down push handles or ultrasuede backrest stripe.</i>	Standard
XBUP3	Tension Adjustable Back Upholstery By Straps, Padded Nylon <i>Included with heavy duty frame. Required for aluminum fendered side guards and fold-down push handles.</i>	\$
	XUPC3 Ultrasuede Backrest Stripe, Black	\$
XBUP5	Tension Adjustable Back Upholstery By Straps, AIR Breathable <i>Included with heavy duty frame. Required for aluminum fendered side guards and fold-down push handles. Not with ultrasuede backrest stripe.</i>	\$
XBUP0	Omit Back Upholstery	No charge

Seat Upholstery

XSUP4	Tension Adjustable Seat Upholstery, Black	Standard
XSUP0	Omit Seat Upholstery	No charge

Seat Cushions

XCSH2	Comfort Curve Cushion <i>Stretch-Air cover.</i>	\$
XCSH23	Comfort Acta-Embrace Comfort-Tek Zero Elevation Cushion	\$
XCSH15	ROHO Hybrid Elite SR Single Compartment Cushion <i>Fits seat widths 14"-19"</i>	\$
XCSH14	ROHO Quadro Select Cushion <i>Fits seat widths 14"-19".</i>	\$
	Low Profile Mid Profile High Profile	
XCSH16	ROHO Hybrid Select Cushion <i>Fits seat widths 14"-18", and 20"</i>	\$

Side Guards

XSDG2	Aluminum Side Guards, Rigid Removable	\$
	Small Standard	
XSDG3	Aluminum Side Guards, Quad Removable w/Hole	\$
	Small Standard	
XSDG7	Carbon Side Guards, Rigid Removable <i>2x2 Twill Weave</i>	\$
	Small Standard	
XSDG8	Carbon Side Guards, Quad Removable w/Hole <i>2x2 Twill Weave</i>	\$
	Small Standard	
XSDG6	Aluminum Side Guards, Fendered Adjustable <i>Requires tension adjustable by straps back upholstery. Not available with angle adjustable fold-down back.</i>	\$
CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 0.5" and side guards type in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>	\$

Armrests

XARM1	Swing Away Tubular Armrests <i>Select range. Not available with heavy duty frame, folding stabilizer bar with short or medium back height ranges, or depth adjustable back.</i>	\$
	9"-12" 12"-15"	
XARM4	Swing Away Flip-Back Tubular Armrests <i>Not available with heavy duty frame, folding stabilizer bar or depth adjustable back.</i>	\$
XARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard <i>Select range and pad size. Requires 14" minimum seat depth.</i>	\$
	9"-12" 12"-15"	
	Desk Length (10") Full Length (14")	
XARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop <i>Select range and pad size. Requires 16" minimum seat depth.</i>	\$
	9"-12" 12"-15"	
	Desk Length (10") Full Length (14")	
XARM5	Armrests Covered in Black Ultrasuede <i>Not available with desk arms.</i>	\$

Positioning Belts

XBLT1	Velcro Adjustable Belt, Black Nylon Webbing	\$
XBLT2	Bodypoint 2" Belt, Push-Button Buckle (Auto Style)	\$
XBLT4	Bodypoint 1" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19"-22" rear seat widths.</i>	\$
XBLT5	Bodypoint 1.5" Padded Belt, Push-Button Buckle (Auto Style) <i>Not available with 19"-22" rear seat widths.</i>	\$
XBLT10	Bodypoint Optional Band Clamp Mounting Hardware	\$
XBLT16	Bodypoint EvoFlex Pelvic Stabilizer, Push Button	\$

Calf Straps

XCLF2	Velcro Adjustable	\$
	Small (9"-12" Footrest Widths) Medium (11"-15" Footrest Widths) Large (14"-17" Footrest Widths) Xtra Large (16"-19" Footrest Widths)	
XCLF3	Bodypoint® Padded Velcro Adjustable	\$
	Medium (12"-15" Footrest Widths) Large (16"-20" Footrest Widths)	
XCLF1	Buckle Adjustable, Black Nylon Webbing	\$

Residual Limb Supports

XRSL1	Residual Limb Support, 8"W x 10"L, Single	\$
	Left Right	
XRSL2	Residual Limb Support, 8"W x 10"L, Pair	\$

Anti-Tips

XTIP1	Rear Anti-Tips, Aluminum Black	\$
XTIP2	Rear Anti-Tips, User Friendly Flip-Up	\$
XTIP9	Tip Assist	\$

Accessories

XGRD1	Clear Spoke Guards <i>Not available with mag rear wheels.</i>	\$
XBCK3	Permobil Medical Necessities Backpack	\$
XPCH3	Permobil Medical Necessities Seat Pouch	\$
XPCH4	Bodypoint Mobility Bag	\$
XUEST1	Upper Extremity Support Tray <i>Requires armrests.</i>	\$
	12"-14" Seat Widths 13"-15" Seat Widths	
XLUG1	Luggage Carriers <i>Not available with high mount swing away footrests or 1-piece flip-up footrest.</i>	\$
XMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe <i>Not available with high mount footrests, 90° hangers or elevating legrests.</i>	\$
XCH1	Crutch Holder	\$

SmartDrive Power Assist System

XMAXSD1	Permobil SmartDrive MX2+ <i>Not available with rear seat heights less than 13.5".</i>	\$
XMAXSD18	Switch Control Buttons	No charge
XMAXSD8	SwitchControl Button with Mono Jack Option	No charge
XMAXSD19	Buddy Button <i>For use with SwitchControl button with mono jack option.</i>	No charge
XMAXSD4	Carrying Bag	\$
XMAXSD10	Universal Padded Handle	\$

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams

