



# TiLite Z Series - Order Form



Dealer Name\*: \_\_\_\_\_ Account #\*: \_\_\_\_\_

City, Post Code, Country: \_\_\_\_\_

Phone: \_\_\_\_\_ Email: \_\_\_\_\_ P.O. Number \_\_\_\_\_

Clinician: \_\_\_\_\_ Client Reference\*: \_\_\_\_\_ Client Weight: \_\_\_\_\_

Rep Name: \_\_\_\_\_

Special considerations that need to be addressed here (e.g., diagnosis): \_\_\_\_\_

\*Required field

## Frame Style

AZFS1	<b>Aero Z Aluminum Frame</b> Full adjustability. 265lb (120kg) weight limit. 150lb (68kg) weight limit for seat sizes 10"x10" to 12"x12".	\$
Z2FS2	<b>ZRA Titanium Frame</b> Full adjustability. 265lb (120kg) weight limit. 150lb (68kg) weight limit for seat sizes 10"x10" to 12"x12".	\$
Z2FS1	<b>ZR Titanium Frame</b> Minimal adjustability. 265lb (120kg) weight limit. 150lb (68kg) weight limit for seat sizes 10"x10" to 12"x12".	\$

## Frame Upgrades

ZHDR1	<b>Heavy Duty Frame</b> Not available with Aero Z. 300lb (136kg) weight limit. No charge when 20" rear seat width and 20" seat depth are selected. Includes billet aluminum front wheels. Requires tension adjustable by straps back upholstery, Shadow, Spinergy SPOX or LX rear wheels. Only available with 0° or 2° camber. Not available with Slipstream or Glide forks, carbon fiber camber tube, velcro adjustable upholstery, swing away tubular armrests or swing away flip back tubular armrests.	No charge
-------	--	-----------

## Transit Options

ZTTDO3	<b>Transit Tie-Downs</b> Not available with ZR. Tested to WC-19. Consists of four tie-down brackets, two pelvic belt mounting brackets and reinforced frame. Please refer to the Aero Z and ZRA Transit Tie-Down Option Brochure for important information regarding transit tie-down options.	\$
ZWTB	<b>Transport Brackets</b> Consists of four tie-down brackets. Unoccupied use only. Please refer to the ZR Transport Bracket Supplement for important information regarding transit options.	\$

## Frame Finishes

ZFF1	<b>Satin Finish</b> Available with titanium frames only.	No charge																								
	ZSCBT1 Buffing Pads for Satin Frame	\$																								
ZFF	<b>Powder Coated Frame</b> Select below.	No charge																								
<b>Frame Color:</b>	<table border="0"> <tr> <td>ZFF8 High Gloss Black</td> <td>ZFF10 Red Pepper</td> <td>ZFF29 Neon Pink</td> <td>ZFF14 Ocean Blue Metallic</td> </tr> <tr> <td>ZFF15 Tangerine Metallic</td> <td>ZFF30 Acid Green Metallic</td> <td>ZFF17 Midnight Blue Pearl</td> <td>ZFF18 White Pearl</td> </tr> <tr> <td>ZFF19 Black Lava</td> <td>ZFF27 Titanium Grey</td> <td>ZFF59 Pomegranate Red</td> <td>ZFF62 Blue Mist</td> </tr> <tr> <td>ZFF57 Indigo Purple</td> <td>ZFF40 Hyper Blue</td> <td>ZFF41 Flame Orange</td> <td>ZFF42 OD Green</td> </tr> <tr> <td>ZFF60 Silver Streak</td> <td>ZFF63 Lavendar Dusk</td> <td>ZFF58 Black Galaxy</td> <td>ZFF61 Seafoam Green</td> </tr> <tr> <td>ZFF64 Daffodil Yellow</td> <td></td> <td></td> <td></td> </tr> </table>	ZFF8 High Gloss Black	ZFF10 Red Pepper	ZFF29 Neon Pink	ZFF14 Ocean Blue Metallic	ZFF15 Tangerine Metallic	ZFF30 Acid Green Metallic	ZFF17 Midnight Blue Pearl	ZFF18 White Pearl	ZFF19 Black Lava	ZFF27 Titanium Grey	ZFF59 Pomegranate Red	ZFF62 Blue Mist	ZFF57 Indigo Purple	ZFF40 Hyper Blue	ZFF41 Flame Orange	ZFF42 OD Green	ZFF60 Silver Streak	ZFF63 Lavendar Dusk	ZFF58 Black Galaxy	ZFF61 Seafoam Green	ZFF64 Daffodil Yellow				
ZFF8 High Gloss Black	ZFF10 Red Pepper	ZFF29 Neon Pink	ZFF14 Ocean Blue Metallic																							
ZFF15 Tangerine Metallic	ZFF30 Acid Green Metallic	ZFF17 Midnight Blue Pearl	ZFF18 White Pearl																							
ZFF19 Black Lava	ZFF27 Titanium Grey	ZFF59 Pomegranate Red	ZFF62 Blue Mist																							
ZFF57 Indigo Purple	ZFF40 Hyper Blue	ZFF41 Flame Orange	ZFF42 OD Green																							
ZFF60 Silver Streak	ZFF63 Lavendar Dusk	ZFF58 Black Galaxy	ZFF61 Seafoam Green																							
ZFF64 Daffodil Yellow																										
ZFF3	<b>Tattoos</b> Tattoos on front frame tubes. Titanium frame finished in satin. Tattoos will be trimmed and adjusted to fit your frame.	\$																								
	ZFF65 Liquid   ZFF66 Bloom   ZFF67 Rebel   ZFF68 Vector																									



Liquid



Bloom



Rebel



Vector

## Color Anodize Package

<b>Standard:</b>	Black	Standard
<b>Anodized Package:</b>	Gold   Green   Purple   Orange   Blue   Red	\$

Includes caster mount, caster fork, caster caps and ZR camber tube bracket. If no anodize colour is selected, default will be black.

## TiFit: Tailored Individualized Frame

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or view our Creative Mobility Division (CMD) catalog [here](#) to view additional options. Please refer to Reference Diagrams at end of form for more information.

**Rear Seat Width 'A'** | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

ZSW1:	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	
ZSW2:	20"										\$

>20" sizes available through CMD.

**Front Seat Width (Seat Taper) 'T'** | Measure from inside of front frame tube to inside of opposite front frame tube. Required with V-front end, not with standard front end.

ZFSW1	<b>No Taper</b>	2.5" narrower than rear seat width.	Standard
ZFSW3	<b>1.5" Taper</b>	4" narrower than rear seat width. Includes tension adjustable by straps seat upholstery. Not with 10" and 11" seat widths. Bolt-on adjustable sling available through our Creative Mobility Division #20.	\$

**Seat Depth 'B'** | Measure from front of back post to the front edge of seat sling. Frame bend will start 0.5" in front of seat sling unless Custom Frame Depth is specified or the combination of specifications affect the rolling dynamics of the chair.

ZSD1:	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	
ZSD2:	20"										\$
ZSD3 Custom Frame Depth:	Specify:										No charge

Measure from front of seat sling to beginning of frame bend. Ranges from 0.5"-2" and can be given as specific as 0.25" increments. A custom frame depth outside the specified range may be available through Creative Mobility Division. A custom frame depth of 20" will have an additional charge.

**Ergonomic Seat 'U'** | Measure from front of back post to beginning of upward bend of seat tube.

ZERG1	<b>ZR Ergonomic Seat</b>	Comes standard with tension-adjustable by strap seat upholstery. Bolt-on adjustable sling available through our Creative Mobility Division #20. Solid seat available through our Creative Mobility Division CMD #16.	\$			
	5"	6"	7"	8"	9"	
ZERG2	<b>Aero Z and ZRA Ergonomic Seat</b>	Center of gravity adjustment may be affected by this option. Not with V-front or seat taper, seat back that locks when folded, desk armrests, depth adjustable backrests, fold down rigid sideguards or fendered adjustable sideguards. Bolt-on adjustable sling available through our Creative Mobility Division #20. Solid seat available through our Creative Mobility Division CMD #16. Minimum seat depth 10" for 6". Minimum seat depth 12" for 7". Minimum seat depth 13" for 8".	\$			
	6"	7"	8"			

**Front Seat Height 'C'** | Measure from floor to top of seat tube at front edge of seat sling.

ZFSH1:	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"	20"	20.5"	21"
--------	-----	-------	-----	-------	-----	-------	-----	-------	-----	-------	-----

**Rear Seat Height 'D'** | Measure from floor to top of seat tube at back post.

ZRSH1:	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"
	18.5"	19"	19.5"	20"								

Rear Wheel	26"	15.5"-20"
	25"	14.5"-20"
	24"	14.5"-20"
	22"	14"-19"
	20"	12.5"-18.5"

Rear Seat Height Range

**Front Angle or Overall Frame Length** | Select one. If Overall Frame Length is selected, front angle will be bent to achieve requested length.

ZFA1	<b>Frame Angle 'L'</b>	Measure from frame front to floor. 85° not available with 6" front wheels. 90° not available with 5" or 6" front wheels.	Standard			
	70°	75°	80°	85°	90°	
ZFA2	<b>Overall Frame Length 'M'</b>	Measure from front of backrest to front edge of footrest.	No charge			
	Specify:					

## Footrest

### Standard Footrests

<b>Seat to Footrest 'E':</b>	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"
	<i>Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height 'C'.</i>											
ZFTR3	<b>Angle Adjustable Footrest</b>											Standard
ZFTR1	<b>Titanium Footrest w/Flat ABS Cover</b>											\$
ZFTR2	<b>Titanium Open Loop Footrest</b>											\$
ZFTR12	<b>Angle Adjustable Flip-Back Footrest</b> <i>Not available with 10" or 11" rear seat width or heavy duty frame.</i>											\$

### High Mount Footrests

<b>Seat to Footrest 'E':</b>	7"-9"	9"-11"	11"-13"	Set at	
	<i>Measure from front edge of seat sling to top rear of footrest.</i>				
ZFTR4	<b>Angle Adjustable High Mount Footrest</b>				\$

### Footrest Width | Select standard front end 'I' if no seat taper 'T' has been selected. Select V-front end 'V' if a seat taper 'T' has been selected.

ZFE1 **No V-Front Taper 'I'** Find your rear seat width 'A' in the chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube. No charge

Rear Seat Width 'A'	10"	7.5"		
	11"	8.5"		
	12"	9.5"	8"	
	13"	10.5"	9"	
	14"	11.5"	10"	
	15"	12.5"	11"	9"
	16"	13.5"	12"	10"
	17"	14.5"	13"	11"
	18"	15.5"	14"	12"
	19"	16.5"	15"	13"
	20"	17.5"	16"	14"
		No Taper	1.5" Taper	3.5" Taper
		Footrest Width 'I'		

ZFE2 **V-Front Taper 'V'** Measure from inside of front frame tube to inside of opposite front frame tube 2.5" above footrest. Titanium footrests have 1" of adjustment up and down. Not with ergonomic seat frames.  
9" 10" 11" 12" 13" 14" 15"

## Back Supports

### ROHO Back Supports

SBRS14	<b>Agility Active Aluminum, Min Contour with Foam</b> Includes fixed hardware. Select shell height.	\$
SBRS15	<b>Agility Active Aluminum, Mid Contour with Foam</b> Includes fixed hardware. Select shell height.	\$
	8" 10"	
SBRS7	<b>Agility Carbon Min Contour Back with Foam/Air Insert</b> Includes fixed hardware. Select shell height.	\$
ACB1	<b>Agility Carbon Min Upgrade</b> Required with Agility Carbon Backs. Not available with bolt-on push handles.	\$ 0.00
SBRS9	<b>Agility Carbon Mid Contour Back with Foam/Air Insert</b> Includes fixed hardware. Select shell height.	\$
ACB2	<b>Agility Carbon Mid Upgrade</b> Required with Agility Carbon Backs. Not available with bolt-on push handles.	\$ 0.00
	10" 13" 16"	

### Comfort Back Supports

<b>SBRS22</b>	<b>Comfort Embrace Back</b> Replaces Comfort Acta-Embrace Back. Includes 2-point quick release Compass hardware. Select shell height.	\$
	10" 12" 14" 16"	
	Comfort-Tek Cover <b>Diamond Spacer Cover</b>	
	Privacy Shield	\$
<b>SBRS23</b>	<b>Comfort Relief Back</b> Replaces Comfort Acta-Relief Back. Includes 2-point quick release Compass hardware. Select shell height.	\$
	10" 12" 14" 16"	
	Comfort-Tek Cover <b>Diamond Spacer Cover</b>	
	Privacy Shield	\$

## Back Upholstery

ZBUP2	<b>Velcro Adjustable</b> <i>Not available with fold-down push handles.</i>	Standard
ZBUP3	<b>NEW Tension Adjustable By Straps, Padded Nylon</b> <i>Required for fold-down push handles.</i>	\$
ZBUP5	<b>Tension Adjustable By Straps, AIR Breathable</b> <i>Required for fold-down push handles.</i>	\$
ZUPC3	<b>Ultrasuede Backrest Stripe, Black</b> <i>Not available with AIR breathable upholstery.</i>	\$
ZBUP0	<b>Omit Back Upholstery</b>	No charge

## Backrest

### Back Type

ZSB1	<b>Folding Aluminum Adjustable Height Back</b> <i>Required for heavy duty frame.</i>	Standard
ZSB2	<b>Folding Titanium Adjustable Height Back, Rounded Rigidizer Bar</b> <i>Titanium back will match frame finish.</i>	\$
ZSB15	<b>Non-Folding Fixed Height Titanium Back, Rounded Rigidizer Bar</b> <i>Not available with Aero Z or ZRA, heavy duty frame, transport brackets or push handles. Titanium back will match frame finish.</i>	\$

### Back Options

ZSB11	<b>Backrest Locks When Folded</b> <i>Available on all backs. Not available with ergonomic seat.</i>	Standard
ZSB12	<b>4" Deep Backrest Rigidizer Bar</b> <i>Not a rounded style rigidizer bar. Required for heavy duty frame. Not available with non-folding fixed height titanium back type.</i>	No charge
ZSB16	<b>Depth Adjustable Back Frame</b> <i>Seat depths of 10"-18" allow 2" of rearward adjustment. Height adjustable push handles have a limited height range. Not with 19" or 20" seat depths, heavy duty frame, tension adjustable by straps seat upholstery or carbon fiber side guards.</i>	No charge
ZSB3	<b>TiShaft Back Release Bar, Aluminum</b> <i>Not available with non-folding fixed height titanium backs.</i>	Standard
ZSB5	<b>TiShaft Back Release Bar, 2x2 Twill Weave Carbon Fiber</b> <i>Not available with non-folding fixed height titanium backs.</i>	\$

### Push Handles

ZSB6	<b>Integrated Push Handles</b>	No charge
ZSB7	<b>Integrated Fold-Down Push Handles</b> <i>Requires tension adjustable by straps back upholstery.</i>	\$
ZSB8	<b>Bolt-On Push Handles</b> <i>Includes 4" deep backrest rigidizer bar. Not available with titanium back.</i>	\$
ZSB9	<b>Bolt-On Height-Adjustable Push Handles</b> <i>Includes 4" deep backrest rigidizer bar.</i>	\$

### Folding Back Height 'F' | ZSB1/ZSB2. Set at 0.5" increments.

Push Handles	None / Bolt-On	9.5"-11.5"	12"-16.5"	16"-20.5"	Set at:
	Integrated	11"-12.5"	13.5"-16.5"	17.5"-20.5"	Set at:
	Fold-Down	13.5"-15"	14.5"-18.5"	18.5"-21"	Set at:
		Short	Medium	Tall	Specify
Folding Back Height					

### Fixed Back Height 'F'

Select height:	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"
	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
	20"											

*Required for fixed side guards. Not with ZRA, swing-away tubular or swing-away flip-back tubular armrests, integrated push handles, integrated fold-down push handles or bolt-on height adjustable push handles.*

### Back Angle 'G' | ZSBA1 Measure angle from front of back post to floor.

Folding Back:	80°	83°	86°	89°	92°	95°	98°	101°				
Fixed Back:	80°	81°	82°	83°	84°	85°	86°	87°	88°	89°	90°	91°
	93°	94°	95°	96°	97°	98°	99°	100°	101°			

## Rear Wheel Positioning

**Center of Gravity 'H'** | Measure from front of back post at the seat tube to the center of the rear axle. If no center of gravity is selected the chair will ship with the center of gravity in the most rearward position. Ergo seat selection may limit available center of gravity. **0° not available with swing-away push handles.**

ZCOG1:	0"	0.25"	0.5"	0.75"	1"	1.25"	1.5"	1.75"	2"	2.25"	2.5"	2.75"
	3"	3.25"	3.5"	3.75"	4"	4.25"	4.5"	4.75"	5"			

**Rear Wheel Spacing 'J'** | Measure from outside of back post to inside of rear tire. **0.75" rear spacing not available with 0° camber, and not available with 2° camber for ZRA.** Ergonomic seat selection may limit available center of gravity.

ZRWS1	<b>Fixed Rear Wheel Spacing</b>											Standard
		0.75"	1"	1.25"								
ZCBR12	<b>Adjustable Rear Wheel Spacing Camber Plugs</b>											No charge
ZCBR9	<b>Additional Spacers to Increase Rear Wheel Spacing</b>											No charge
ZCBR11	<b>Zero Play Aluminum Camber Plugs</b>	Not available with power adaptable frame or more than 6° camber.										No charge

## Camber 'K'

ZCBR1:	0°	2°	4°	6°								
	More than 0° camber not with power adaptable frame. 8° camber may be available through our Creative Mobility Division.											
ZCBR15	<b>Aluminum Camber Tube</b>											Standard
ZCBR3	<b>Titanium Camber Tube</b>											\$
ZCBR4	<b>Carbon Fiber Camber Tube</b>	Only available on ZR and only with 0°, 2°, 4° zero play plugs. Not available with anti-tips, heavy duty frame.										\$

## Seat Upholstery

ZSUP1	<b>Tension Adjustable Bolt-On Upholstery, Black</b>	Not available with ergonomic seat or tapered seat.										Standard
ZSUP2	<b>Tension Adjustable By Straps, Black</b>	Standard with ergonomic seat and seat taper of 1"-3".										\$
ZSUP8	<b>Aluminum Solid Seat</b>	Ergonomic seat requires CMD #16. Not with 10" or 11" rear seat width.										\$
ZSUP0	<b>Omit Seat Upholstery</b>											No charge
CMD018	<b>Custom iSleeve</b>	Please provide sleeve width and depth (inches) in Notes section at the end of this order form. Add 1" to device depth and 0.25" to width. <b>Best when used with tension adjustable bolt-on upholstery.</b> Some limitations may apply.										\$

## Seat Cushions

ZCSH2	<b>Comfort Curve Cushion</b>	Stretch-Air cover.										\$
ZCSH23	<b>Comfort Acta-Embrace Comfort-Tek Zero Elevation Cushion</b>											\$
ZCSH11	<b>ROHO Single Compartment</b>	Fits seat widths 14"-19".										\$
		Low Profile	Mid Profile	High Profile								
ZCSH14	<b>ROHO QUADTRO SELECT Cushion</b>	Fits seat widths 14"-19".										\$
		Low Profile	Mid Profile	High Profile								
ZCSH16	<b>ROHO Hybrid Select Cushion</b>	Fits seat widths 14"-18", and 20"										\$

## Front Wheels

ZFW24	<b>3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire</b>											\$
		Silver Hub	Black Hub									
ZFW14	<b>4"x0.75" Plastic Wheel w/Poly Tire</b>	Not available with Slipstream fork.										No charge
ZFW3	<b>4"x0.75" Light-Up Micro Caster</b>	Not available with Slipstream fork.										\$
ZFW26	<b>4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire</b>	Not available with power assist or heavy duty frame.										No charge
ZFW27	<b>4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire</b>											\$
		Silver Hub	Black Hub									
ZFW6	<b>5"x1" Plastic Wheel w/Poly Tire</b>	Not available with 90° front angle or Slipstream fork.										No charge
ZFW28	<b>5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire</b>	Not available with 90° front angle.										No charge
ZFW30	<b>5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire</b>	Not available with 90° front angle.										\$
		Silver Hub	Black Hub									
ZFW9	<b>6"x1" Plastic Wheel w/Poly Tire</b>	Not available with 85° or 90° front angle or Slipstream fork.										No charge
ZFW32	<b>6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire</b>	Not available with 85° or 90° front angle or Slipstream fork.										No charge
ZFW33	<b>6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire</b>	Not available with 85° or 90° front angle.										\$
		Silver Hub	Black Hub									

## Precision Lock Front Forks

ZFK1	<b>TiLite Standard Fork</b>	Standard
ZFK2	<b>TiLite Slipstream Single-Sided Fork</b> <i>Not available with heavy duty frame.</i>	\$
ZGL1	<b>Glide Suspension Fork</b> <i>Black anodize only. Not available with heavy duty frame. Select rider weight.</i>	\$
	125lbs (57kg)    125-180lbs (57-82kg)    180-265lbs (82-120kg)	
ZTFS1	<b>Titanium Fork Stem</b> <i>Comes standard on ZR.</i>	\$

## Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

ZRW14	<b>TiLite SHADOW Rear Wheel</b> <i>Select size. 25" and 26" wheels require 15" minimum seat depth.</i>	Standard
	20" (451)    22" (501)    24" (540)    25" (559)    26" (590)	
ZRW7	<b>Spinergy SPOX Rear Wheel</b> <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth. Orange, pink, purple and green spokes only available in 20" size.</i>	\$
	20" (451)    22" (501)    24" (540)    25" (559)    26" (590)	
	Black   Red   Blue   White   Orange   Pink   Purple   Green	
ZRW8	<b>Spinergy LX Rear Wheel</b> <i>Select size and spoke color. 25" and 26" wheels require 15" minimum seat depth.</i>	\$
	22" (501)    24" (540)    25" (559)    26" (590)	
	Black   Red   Blue   White   Orange   Pink   Purple   Green	
ZRW9	<b>Spinergy Design Lab LX Rear Wheel</b> <i>Black spokes only. Select size and anodizing color. If same color is selected, anodized hubs and rims match colour anodized package. 25" and 26" wheels require 15" minimum seat depth.</i>	\$
	22" (501)    24" (540)    25" (559)    26" (590)	
	Gold   Green   Blue   Red   Purple   Orange	
ZRW16	<b>Spinergy Carbon Blade Rear Wheel</b> <i>2x2 twill weave carbon fiber. Not available with heavy duty frame. 25" wheels require 15" minimum seat depth.</i>	\$
	24" (540)    25" (559)	
ZRW17	<b>X-Core Mag Rear Wheel</b> <i>Select size. Not available with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i>	\$
	20" (451)    22" (501)    24" (540)	

## Axles

ZAXL1	<b>Stainless Quick Release Axle</b>	Standard
ZAXL2	<b>Titanium Quick Release Axle</b>	\$
ZAXL3	<b>Quad Release Axle</b> <i>Not available with Mag or 26" (590) Spinergy Spox or LX rear wheels.</i>	\$

## Pneumatic Rear Tires

ZRTR1	<b>Treaded Tire 22"(501), 24"(540), 26"(590)</b>	Standard
ZRTR2	<b>Airless Inserts</b> <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$
ZRTR31	<b>MBL SpeedLite Tire 20"(451), 22"(501), 24"(540), 25"(559), 26"(590)</b>	No charge
ZRTR10	<b>MBL Gazelle Tire 24"(540), 25"(559)</b> <i>Kevlar bead, lightest weight option. Not available on the Aero Z.</i>	\$
ZRTR30	<b>Primo Sentinel Tire 24"(540), 25"(559), 26"(590)</b>	\$
ZRTR8	<b>Schwalbe Marathon Plus Evolution Tire 24"(540), 25"(559), 26"(590)</b>	\$
ZRTR32	<b>Primo C-764 Knobby Tire 24"(540)</b>	\$

## Solid Rear Tires

ZRTR22	<b>Primo Orion Tire 20"(451), 22"(501), 24"(540), 25"(559)</b> <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$
ZRTR23	<b>Primo Active Tire 24"(540), 25"(559)</b> <i>Not available with Spinergy Carbon Blade rear wheels.</i>	\$

## Handrims

ZRIM1	Aluminum Handrim, Silver Anodized	Standard
ZRIM21	Aluminum Handrim, Black Anodized	No charge
ZRIM12	Titanium Handrim, Satin Finish <i>Not available with 20"(451) or 22"(501) rear wheels or 24"(540) knobby tires.</i>	\$
ZRIM20	Q-Grip Handrim <i>Exclusive high-friction, anti-slip coating. Not available with 20"(451) rear wheels.</i>	\$
ZRIM5	Plastic Coated Handrim, Black	\$
ZRIM10	The Surge Handrim <i>Large profile. All-In-One oval with Gription strip.</i>	\$
ZRIM11	The Surge LT Handrim <i>Small profile. All-In-One oval with Gription strip.</i>	\$
ZRIM14	The Natural-Fit Handrim, Standard Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM15	The Natural-Fit Handrim, Super Grip <i>Original, large profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM16	The Natural-Fit Handrim, Without Thumb Grip <i>Original, large profile.</i>	\$
ZRIM17	The Natural-Fit LT Handrim, Standard Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM18	The Natural-Fit LT Handrim, Super Grip <i>Reduced weight, small profile. Not available with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM19	The Natural-Fit LT Handrim, Without Thumb Grip <i>Reduced weight, small profile.</i>	\$
ZRIM22	Omit Handrims	No charge

## Handrim Tab Length

ZRIM7	Long Tabs	Standard
ZRIM8	Short Tabs <i>Not available with Mag rear wheels, The Natural-Fit handrims with thumb grip, or Natural-Fit handrims on Shadow wheels.</i>	\$

## Wheel Locks

ZWLK7	Compact Composite Push to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
ZWLK8	Compact Composite Pull to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
ZWLK10	Extension Handles for Compact Composite Push/Pull to Lock	\$
ZWLK11	MBL Composite Scissor Lock	\$
ZWLK12	MBL Performance Aluminum Scissor Lock	\$
ZWLK3	Uni-Lock Push Pull	\$
ZWLK4	Extension Handles for Uni-Lock	\$
ZWLK15	Grade Aids Push to Lock <i>Not available with solid rear tires.</i>	\$

## Side Guards

TSDG13	Side Guards Custom Cut to Profile of Rear Wheels <i>Available for aluminum and carbon non-fendered side guards.</i>	\$
CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 0.5" and side guards type (rigid or fold-down) in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>	\$
CMD	Fixed Side Guards <i>CMD33 (fendered carbon fiber) or CMD73 (flat panel aluminum or carbon fiber) available. ZRA also requires CMD19 (fixed back). Not available on the Aero Z. Please visit our CMD catalogue, pricing will vary.</i>	\$

## Aluminum

ZSDG2	Aluminum Side Guards, Rigid Removable Small Standard	\$
ZSDG14	Aluminum LT Side Guards, Rigid Removable <i>Not available on the Aero Z.</i> Small Standard	\$
ZSDG3	Aluminum Side Guards, Quad Removable w/Hole Small Standard	\$
ZSDG15	Aluminum LT Side Guards, Quad Removable w/Hole <i>Not available on the Aero Z.</i> Small Standard	\$
ZSDG4	Aluminum Fold Down Side Guards, Rigid <i>Not available with ergonomic seats.</i> Small Standard	\$
ZSDG6	Aluminum Fendered Side Guards, Adjustable <i>24"-25" rear wheels require rear seat height 14.5"-19.5". 26" rear wheels require rear seat height 15.5"-18.5".</i>	\$

## 2x2 Twill Weave Carbon Fiber

ZSDG7	<b>Carbon Side Guards, Rigid Removable</b>	\$
	Small Standard	
ZSDG16	<b>Carbon LT Side Guards, Rigid Removable</b> <i>Not available on the Aero Z.</i>	\$
	Small Standard	
ZSDG8	<b>Carbon Side Guards, Quad Removable w/Hole</b>	\$
	Small Standard	
ZSDG17	<b>Carbon LT Side Guards, Quad Removable w/Hole</b> <i>Not available on the Aero Z.</i>	\$
	Small Standard	
ZSDG9	<b>Carbon Fold Down Side Guards, Rigid</b> <i>Not available with ergonomic seats.</i>	\$
	Small Standard	
ZSDG12	<b>Carbon Fendered Side Guards, Removable</b>	\$

## Armrests

ZARM1	<b>Swing Away Tubular Armrests</b> <i>Select range. Not available with non-folding fixed height titanium back or height adjustable push handles.</i>	\$
	6.5"-9.5" 9.5"-12.5"	
ZARM6	<b>Swing Away Fixed Height Flip Back Tubular Armrests</b> <i>Not available with heavy duty frame.</i>	No charge
ZARM4	<b>Swing Away Flip Back Tubular Armrests, 8.5"-11.5"</b> <i>Not available with non-folding fixed height titanium back.</i>	\$
ZARM2	<b>Removable Height Adjustable Desk Arm with Rigid Side Guard</b> <i>Select range and pad size. Requires 14" minimum seat depth. Not available with ergonomic seat.</i>	\$
	9"-12" 12"-15"	
	Desk Length (10") Full Length (14")	
ZARM3	<b>Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop</b> <i>Select range and pad size. Requires 15" minimum seat depth. Not available with ergonomic seat.</i>	\$
	9"-12" 12"-15"	
	Desk Length (10") Full Length (14")	
ZARM5	<b>Armrests Covered in Black Ultrasuede</b> <i>Not available with desk arms.</i>	\$

## Positioning Belts

ZBLT1	<b>Velcro Adjustable Belt, Black Nylon Webbing</b>	\$
ZBLT2	<b>Bodypoint 2" Belt, Push-Button Buckle (Auto Style)</b>	\$
ZBLT4	<b>Bodypoint 1" Padded Belt, Push-Button Buckle (Auto Style)</b> <i>Not available with 19" or 20" rear seat widths.</i>	\$
ZBLT5	<b>Bodypoint 1.5" Padded Belt, Push-Button Buckle (Auto Style)</b> <i>Not available with 19" or 20" rear seat widths.</i>	\$
ZBLT8	<b>Q'Straint WC-19 Wheelchair-Anchored Pelvic Belt</b> <i>Recommended for use with Transit Tie Down. Not available with ZR.</i>	\$
ZBLT10	<b>Bodypoint Optional Band Clamp Mounting Hardware</b>	\$
ZBLT16	<b>Bodypoint EvoFlex Pelvic Stabilizer, Push Button</b>	\$

## Calf Straps

ZCLF2	<b>Velcro Adjustable</b>	\$
	Small (9"-12" footrest widths) Medium (11"-15" footrest widths) Large (14"-17" footrest widths) Xtra Large (16"-19" footrest widths)	
ZCLF3	<b>Bodypoint Padded Velcro Adjustable</b>	\$
	Small (8"-11" footrest widths) Medium (11"-15" footrest widths) Large (16"-19" footrest widths)	
ZCLF1	<b>Buckle Adjustable, Black Nylon Webbing</b>	\$

## Residual Limb Supports

ZRSL1	<b>Residual Limb Support, 8"W x 10"L, Single</b>	\$
	Left Right	
ZRSL2	<b>Residual Limb Support, 8"W x 10"L, Pair</b>	\$

## Anti-Tips

ZTIP1	<b>Rear Aluminum, Black</b>	\$
ZTIP2	<b>User Friendly Flip-Up</b>	\$

## Accessories

ZGRD1	<b>Clear Lexan Spoke Guards</b> <i>Not available with Mag rear wheels.</i>	\$
ZBCK3	<b>Permobil Medical Necessities Backpack</b>	\$
ZPCH3	<b>Permobil Medical Necessities Seat Pouch</b>	\$
ZPCH4	<b>Bodypoint Mobility Bag</b>	\$
ZUEST1	<b>Upper Extremity Support Tray</b> <i>Requires armrests.</i>	\$
	12"-14" seat widths    13"-15" seat widths	
ZLUG1	<b>Luggage Carriers</b>	\$
ZMPCT2	<b>Black Neoprene Impact Guards with Ultrasuede Stripe</b>	\$
ZCH1	<b>Crutch Holder</b>	\$

## SmartDrive Power Assist System

ZMAXSD1	<b>Permobil SmartDrive MX2+</b> <i>Not available with rear seat heights less than 13.5" or 20" rear wheels.</i>	\$
ZMAXSD9	<b>Battery Pack</b> <i>Enclosed in SmartDrive system.</i>	No charge
ZMAXSD18	<b>Switch Control Buttons</b>	No charge
ZMAXSD8	<b>SwitchControl Button with Mono Jack Option</b>	No charge
ZMAXSD19	<b>Buddy Button</b> <i>For use with SwitchControl button with mono jack option.</i>	No charge
ZMAXSD4	<b>Carrying Bag</b>	\$
ZMAXSD10	<b>Universal Padded Handle</b>	\$

## General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes

**TiFit Reference Diagrams**



**TiLite ZRA**



**TiLite ZR**

