

Panthera X3 - Order form



Dealer Name*:		Account #*:
City, Post Code, Country:		
Phone:	Email:	P.O. Number
Clinician:	Client Reference*:	Client Weight:
Rep Name:	*Required field	
Special considerations that need to be addressed here (e.g., diagnosis):		

Chassis

Standard configuration includes tapered frame with 100° front angle, 35-46cm depth-adjustable upholstery with 7° seat angle, backrest, front castors, rigid footplate, push to lock brakes, Spinergy rear wheels with 2° camber and aluminium handrims. Seat height A1 with 24" rear wheels: front 47cm, rear 43cm. Seat height can be adjusted with different size rear wheels. 125kg weight limit.

G355 Panthera X3 \$

Balance Position | Available in two different balance and height positions. The position cannot be changed later. Seat height 'B' is 15mm lower.

Standard Seat Height: A1 | Most Active A2 | Less Active
Low Seat Height: B1 | Most Active B2 | Less Active

Seat Measurements

Seat Width: 33cm 36cm 39cm 42cm 45cm

Backrest

Standard Backrest, Carbon Fibre Ovalised back tubes. Standard

Backrest Frame for Solid Backrests U3 Light style, round back tubes. 35cm height. Custom request, additional leadtime of 5 days. No charge

Backrest Height: 22cm 23cm 24cm 25cm 26cm 27cm 28cm 29cm 30cm 31cm 32cm 33cm 34cm
 35cm

Only for standard carbon fibre backrest.

Backrest Position: P1 | Most Active P2 | Active P3 | Less Active P4 | Least Active

Backrest position affects the balance. The position can be easily adjusted later.

Rear Wheels

Spinergy X Wheels: 24" 25" 26"
 Black Spokes White Spokes
 1100023-C Quick Release Axles, Steel Required for user weight >100kg. No charge

Tyres

Schwalbe Right Run 24"(540), 25" (559), 26"(590) Standard
 Schwalbe Marathon Plus 24"(540), 25" (559), 26"(590) \$
 Schwalbe One 24"(540), 25" (559) \$
 Rugged 24" (540) \$
 PU Solid 24" (540) \$

Handrims

190DD30-C Aluminium Standard
 190DD11-C Titanium \$
 19044DD-C Aluminium X \$
 19051DDBT-C Paragrip \$

19052	DDBT-C	Tetragrip	\$
19030	DD-C	Maxgrip 24" wheels only. Contains latex.	\$

Front Castors

2100016-C	X Style 87mm (3")	
2100204-C	S3 Style, 90mm (3.5")	
2100205-C	S3 Style, 120mm (4.7")	Not compatible with 24" rear wheels and B1/B2 balance position.

Footrest

36915	BB-C	U-Shape	Standard
36006	BB-C	U-Shape, Extended Reinforced	No charge
36916	BB-C	U-Shape, Small Heel Bar	No charge
36004	BB-C	Aluminium U-Shape	\$
65060	BB-C	Carbon Footplate	\$
36010	BB-C	Plastic Plate for Footrest	Not with Carbon Footrest. \$

Wheel Locks

4630000-C	Push to Lock S3	Standard
462BB02-C	One Hand Brake, Right Activation	No charge
462BB01-C	One Hand Brake, Left Activation	No charge

Side Guards

Standard Side Guards, S3 Style	\$
Carbon Side Guards, Without Fender	\$
Carbon Side Guards, Small Fender	\$
Carbon Side Guards, Large Fender	\$

Push Handles

4200016-C	Arm Hooks, X3	\$
4200017-C	Push Handles, Foldable	\$

Accessories

Seat Cushion

575	BB40-C	Seat Cushion 5cm	\$
575	BB35-C	Seat Cushion 5cm, Short	\$
572	BB40-C	Seat Cushion 2.5cm	\$
575	BB35-C	Seat Cushion 2.5cm, Short	\$
57710	BB-C	Wedge Cushion	Used to alter seat angle. \$
990021	X-C	Seat Bag	\$
		Mounted On Seat Upholstery	Mounted Under Cushion

Safety

4354104-C	Anti-Tips, Pair	\$
-----------	-----------------	----

Other

9024009-C	Transport Protection, Castor Housing	\$
9024020-C	Transport Protection, Rear	\$
9900256-C	Tube Protection, Front Chassis	\$
9900255-C	Tube Protection, Backrest Bar	\$
9900261-C	Backpack	\$
110030D-C	reTyre Skins	\$
	24" 25" 26"	
2600010X-C	FreeWheel Attachment <i>Requires extended reinforced footrest.</i>	\$
	Height 6-11.5cm Height 11.5-16.5cm	
905005X-C	Mototronik Quick Connector <i>Select seat width.</i>	\$
	33-36cm Width 39-42cm Width 45cm Width	


Set-Up Adjustments

Backrest Angle: Specify (mm) | A =
15mm standard.

Footrest Height: Specify (mm) | K =

Notes

Reference Diagrams




Balance Position



Backrest Placement



Backrest Angle Adjustment



Footrest Height