



PDG Mobility™

Comparison matrix

permobil

PDG Manual Wheelchairs



Fuze™ T50

Fuze T50n

Fuze M50

Fuze JR

Bentley™

Bentley LT

Bentley LT-R

Dynamic tilt range	Up to 50°	Up to 50°	Up to 50°	Up to 50°	Dynamic tilt range	Up to 20°	Up to 20°	Up to 20°
Seat width	13" - 24"	16" to 24"	15" - 24"	12" - 18"	Seat width	14" - 32"	14" - 22"	14" - 22"
Seat depth	15" - 22"	15" - 22"	15" - 22"	13" - 18"	Seat depth	16" - 22"	16" - 22"	16" - 22"
Seat to floor height	13" - 20"	16" - 20"	13" - 19"	13" - 20"	Seat to floor height	13" - 20"	13" - 20"	12" - 19"
Overall width	+7.5" - 8.5"	+4.0"	+8.5"	+7.5" - 8.5"	Overall width	+9.5" - 10.5"	+9.5" - 10.5"	+10.5"
Transport weight	20kg	20kg	24kg	19kg	Transport weight	23kg	15kg	15.5kg
Maximum user weight	181kg	181kg	181kg	113kg	Maximum user weight	204kg	113kg	113kg
Frame material	Aluminium	Aluminium	Aluminium	Aluminium	Frame material	Steel	Aluminium	Aluminium
Power tilt option	Yes, up to 159kg	Yes, up to 159kg	No	No	Power tilt option	Yes, up to 113kg	Yes	Yes
Take-Apart frame option	Yes	Yes	No	Yes	Take-Apart frame option	No	No	No
Self tilt	No	No	No	No	Self tilt	Yes	Yes	Yes
Transit tie-downs	Yes, up to 136kg	Yes, up to 136kg	Yes, up to 136kg	Yes, up to 136kg	Transit tie-downs	Yes, up to 136kg	Yes, up to 136kg	Yes, up to 136kg

Tilt Chairs

The Fuze is an adaptable tilt chair with a sleek aluminium frame and a small footprint. It offers up to 50° of effortless tilt. Set-up and adjustments of the Fuze are fast and easy and require only simple tools. The Bentley is a 20° no-knee-rise tilt chair that is ideal for foot propelling. It is available with a rugged steel frame or a with a lightweight aluminium frame weighing only 15kg.

Bariatric Chairs

The Fuze T50 and Bentley offer designs and features suitable for bariatric users. Please refer to the specifications above to find a chair that meets your requirements.

Take-Apart

The Take-Apart frame allows attendants to take the Fuze T50, Fuze T50n or Fuze JR apart in seconds without tools. This feature reduces the weight of the heaviest piece of the wheelchair to approximately 13kg, perfect for convenient transportation in the boot of a car or SUV.



Our innovative chair designs deliver an enhanced quality of life for both users and caregivers. Long product life, adjustability, and durability create value for all stakeholders, and our expertise in positioning wheelchairs ensures an excellent fit every time, no matter how basic or complex.

Seating and Positioning

Permobil provides a wide range of seating and positioning products to address the unique needs of individuals using PDG wheelchairs.

ROHO® and Comfort* offer a variety of cushions, back supports and accessories that are designed for pressure redistribution, stability, and comfort that can be tailored to meet the individual needs of users.

**Coming soon*