


Transportation Made Easy

The Take-Apart frame allows attendants to take the Fuze™ T50, Fuze T50n or Fuze JR apart in seconds without tools. This reduces the weight of the heaviest piece of the wheelchair to approximately 13kg, perfect for convenient transportation in the boot of a car or SUV.



Our innovative chair designs deliver an enhanced quality of life for both users and caregivers. Long product life, adjustability, and durability create value for all stakeholders, and our expertise in positioning wheelchairs ensures an excellent fit every time, no matter how basic or complex.

Seating and Positioning

Permobil provides a wide range of seating and positioning products to address the unique needs of individuals using PDG wheelchairs.

Comfort and ROHO® offer a variety of cushions, back supports and accessories that are designed for pressure redistribution, stability, and comfort that can be tailored to meet the individual needs of users.



Permobil Australia | 1300 845 483 | orders.au@permobil.com | www.permobil.com/en-au
Permobil New Zealand | 0800 115 222 | sales.nz@permobil.com | www.permobil.com/en-nz

BR-ANZ-PDGM-001.0

PDG Mobility™

Comparison matrix



PDG Manual Wheelchairs



Fuze™ T50

Fuze T50n

Fuze T50 Junior

Fuze M50

Stellar GL

Stellar GLT

Stellar Impact

Stellar HD

Bentley™

Bentley LT

Bentley LT-R

Eclipse

| | Fuze™ T50 | Fuze T50n | Fuze T50 Junior | Fuze M50 | Stellar GL | Stellar GLT | Stellar Impact | Stellar HD | Bentley™ | Bentley LT | Bentley LT-R | Eclipse |
|--------------------------------|------------------|------------------|-----------------|-------------|------------------|-------------|----------------|-----------------------------|------------------|-------------|--------------|-------------|
| Dynamic tilt range | 0° to 50° | 0° to 50° | 0° to 50° | 0° to 50° | -5° to 45° | -5° to 45° | 0° to 45° | 0° to 45° (25° at 272kg) | 0° to 20° | 0° to 20° | 0° to 20° | None |
| Seat width | 12" - 24" | 16" to 24" | 12" - 18" | 15" - 24" | 14" - 32" | 14" - 22" | 14" - 22" | 22" - 32" | 14" - 32" | 14" - 22" | 14" - 22" | 20" - 48" |
| Seat depth | 15" - 22" | 15" - 22" | 13" - 18" | 15" - 22" | 16" - 22" | 16" - 22" | 17" - 22" | 18" - 22" | 16" - 22" | 16" - 22" | 16" - 22" | 16" - 24" |
| Seat to floor height | 13" - 20" | 16" - 20" | 13" - 20" | 13" - 19" | 13" - 20" | 13" - 20" | 14" - 21" | 16" - 18" | 13" - 20" | 13" - 20" | 12" - 19" | 14" - 20" |
| Overall width* | +7.5" | +4.0" | +7.5" | +8.5" | +9.5" | +9.5" | +9.5" | +11" | +9.5" | +9.5" | +10.5" | +7.5" |
| Transport weight | 20kg | 20kg | 19kg | 24kg | 23kg | 15kg | 29kg | 28kg | 23kg | 15kg | 15kg | 20kg |
| Maximum user weight | 181kg | 181kg | 113kg | 181kg | 204kg | 113kg | 113kg | 272kg | 204kg | 113kg | 113kg | 454kg |
| Frame material | Aluminium | Aluminium | Aluminium | Aluminium | Steel | Aluminium | Steel | Steel | Steel | Aluminium | Aluminium | Steel |
| Power tilt option | Yes, up to 159kg | Yes, up to 159kg | Yes | No | Yes, up to 113kg | Yes | No | No | Yes, up to 113kg | Yes | Yes | No |
| Take-Apart frame option | Yes | Yes | Yes | No | No | No | No | No | No | No | No | No |
| Transit tiedowns | Up to 136kg | Up to 136kg | Up to 113kg | Up to 136kg | Up to 136kg | Yes | Yes | No | Up to 136kg | Up to 113kg | Up to 113kg | Up to 272kg |

*Overall width measurement is in addition to the seat width.