

Manual wheelchairs

Panthera S3 Swing



permobil

Adapts to the active life you lead.

Experience the Panthera S3 Swing, combining the comfort and lightweight design of the S3 series with the added convenience of swing-away footrests for effortless transfers.

Different frame options are available to suit both smaller and taller users, ensuring a tailored fit for every individual.

Key benefits:

- Made to manoeuvre
- Lightweight chrome molybdenum frame
- Integrated anti-tip device
- Easy to maintain



Specifications

	S3 Swing	S3 Swing Short
Seat width	36 - 39 - 42 - 45 - 50cm	33 - 36 - 39cm
Seat depth	40cm	35cm
Seat height	Front 45cm (24" rear wheels) Back 43cm (24" rear wheels)	Front 45cm (24" rear wheels) Back 43cm (24" rear wheels)
Seat angle	4°	4°
Total length	78 - 90cm	88cm
Total width	Seat width +21cm	Seat width +21cm
Front frame angle	75°	75°
Backrest height	20 - 45cm	20 - 45cm
Backrest angle	13° - (-8°)	13° - (-8°)
Footplate distance	38 - 46cm	38 - 46cm
Camber	2.2°, 4° (seat width 36 - 45cm) 8° (seat width 36 - 39cm)	2.2°, 4°, 6° (seat width 33cm), 8°
Transport weight	4.1 - 4.9kg	3.9 - 4.1kg
Wheelchair weight	9 - 9.9kg grams	8.8 - 9kg
Maximum user weight	100kg (150 kg for seat widths 45 - 50cm)	100kg

Permobil reserves the right to change specifications and design without prior notice.