

Manual wheelchairs

Panthera U3 Series



permobil

Designed for motion. Built for life.

The Panthera U3 is a lightweight, adjustable manual wheelchair designed for active users. It offers a secure, ergonomically optimised seating position and a sleek, slim profile. Enjoy enhanced manoeuvrability and personalise your comfort with a range of adjustable accessories for added support and peace of mind.

- Key benefits:
- Lightweight chrome molybdenum frame
 - Adjusts to individual needs
 - Integrated anti-tip device
 - Easy to manage



Specifications

	U3	U3 Y-front
Seat width	33 - 36 - 39 - 42 - 45 cm	33 - 36 - 39 - 42 - 45 - 50 cm
Seat depth	35 - 46 cm	35 - 46 cm
Seat height	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)	Front 47 cm (w. 24" rear wheels) Back 43 cm (w. 24" rear wheels)
Seat angle	7°	7°
Total length	84 cm	84 cm
Total width	Seat width +21 cm	Seat width +21 cm
Front frame angle	80°	80°
Backrest height	20 - 45 cm	20 - 45 cm
Backrest angle	13° - (-8°)	13° - (-8°)
Footplate distance	35 - 43 cm	35 - 43 cm
Camber	2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)	2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)
Transport weight	4.5 - 5.0 kg	4.5 - 5.4 kg
Wheelchair weight	8.0 - 8.5 kg	8.0 - 9.3 kg
Maximum user weight	100 kg, 150 kg (seat width 45 cm)	100 kg, 150 kg (seat width 45 - 50 cm)

	U3 Y-front Short	U3 Y-front Low
Seat width	33 - 36 - 39 - 42 - 45 - 50 cm	33 - 36 - 39 - 42 - 45 - 50 cm
Seat depth	35 - 46 cm	35 - 46 cm
Seat height	Front 47 cm (24" rear wheels) Back 43 cm (24" rear wheels)	Front 44.5 cm (24" rear wheels) Back 40.5 cm (24" rear wheels)
Seat angle	7°	7°
Total length	79 cm	84 cm
Total width	Seat width +21 cm	Seat width +21 cm
Front frame angle	80°	80°
Backrest height	20 - 45 cm	20 - 45 cm
Backrest angle	13° - (-8°)	13° - (-8°)
Footplate distance	35 - 43 cm	35 - 43 cm
Camber	2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)	2.2°, 4°, 6° (seat width 33 cm) 8° (seat width 33 - 39 cm)
Transport weight	4.5 - 5.4 kg	4.4 - 5.3 kg
Wheelchair weight	8.0 - 9.3 kg	8.0 - 9.3 kg
Maximum user weight	100 kg, 150 kg (seat width 45 - 50 cm)	100 kg, 150 kg (seat width 45 - 50 cm)

Innovating for individuals.