

Manual wheelchairs

Panthera U3 Light

permobil



Comfort that moves with you.

An ultra lightweight, manoeuvrable manual wheelchair, the Panthera U3 Light offers stability and comfortable posture support. Designed for optimal ergonomics when seated or mobilising, this wheelchair can also be easily lifted into a car thanks to its balanced, effortless-grip chassis.

- Key benefits:
- Ultra lightweight
 - Highly manoeuvrable
 - Easy to transport
 - Adaptable to individual needs



Specifications

	U3 Light	U3 Light Y-front
Seat width	33 - 36 - 39 - 42 - 45cm	33 - 36 - 39 - 42 - 45cm
Seat depth	35 - 46cm	35 - 46cm
Seat height	Front 47cm (w. 24" rear wheels) Back 43cm (w. 24" rear wheels)	Front 47cm (w. 24" rear wheels) Back 43cm (w. 24" rear wheels)
Seat angle	7°	7°
Total length	82cm	82cm
Total width	Seat width +21cm	Seat width +21cm
Front frame angle	100°	100°
Backrest height	25 - 35cm	25 - 35cm
Backrest angle	11.5° - (-4.5°)	11.5° - (-4.5°)
Footplate distance	35 - 43cm	35 - 43cm
Camber	2.2°	2.2°
Transport weight	3.4 - 4.0kg	3.7 - 3.9kg
Wheelchair weight	5.7 - 6.3kg	6.0 - 6.2kg
Maximum user weight	100kg	100kg

U3 Light Y-front 90°	U3 Light Y-front High 4°	U3 Light Y-front Low
33 - 36 - 39 - 42 - 45cm	33 - 36 - 39 - 42 - 45cm	33 - 36 - 39 - 42 - 45cm
35 - 46cm	35 - 46cm	35 - 46cm
Front 47cm (24" rear wheels) Back 43cm (24" rear wheels)	Front 51cm (24" rear wheels) Back 45cm (24" rear wheels)	Front 44.5cm (w. 24" rear wheels) Back 40.5cm (w. 24" rear wheels)
7°	4°	7°
73cm	82cm	82cm
Seat width +21cm	Seat width +21cm	Seat width +21cm
90°	100°	100°
25 - 35cm	25 - 35cm	25 - 35cm
11.5° - (-4.5°)	11.5° - (-4.5°)	11.5° - (-4.5°)
35 - 43cm	36 - 44cm	36 - 44cm
2.2°	2.2°	2.2°
4.0 - 4.2kg	3.7 - 3.9kg	3.7 - 3.9kg
6.3 - 6.5kg	6.0 - 6.2kg	6.0 - 6.2kg
100kg	100kg	100kg

Innovating for individuals.