

Comfort Embrace

Features

- Integrated wings attach to a planar aluminum back support shell and are adjustable for angle, height, and up to a total of 7.6 cm of width for standard wings and 3.8 cm for pediatric wings
- Moldable stays included for positioning, support and pressure redistribution for a variety of spinal curvatures and changes in conditions
- Updated cover design and 2 layer foam construction
- Comes with integrated chest harness strap guides
- Durable aluminum shell comes with pre-drilled 6 x 2.5 cm head support hole pattern for BodiLink Head Support
- Compatible with lateral trunk support hardware. A total of 28 cm of pad width adjustment (14 cm per side) & over 13 cm of height adjustment
- Cover options: Diamond Spacer Standard & Comfort-Tek Optional
- Compatible with BodiLink® Compass Hardware
- Warranty: 2 years back and hardware, 6 months for cover



Standard sizes

WIDTH

HEIGHT

	25.0 cm	30.5 cm	35.5 cm	40.5 cm
45.5 cm			35.5 x 45.5	
40.5 cm			35.5 x 40.5	
35.5 cm	25.0 x 35.5	30.5 x 35.5		

Compass 4 Hardware offer 5 cm width adjustments

Dimensions

Width	25.0–35.5 cm
Height	35.5–45.5 cm

Weight capacity

25.0 – 45.5 cm wide	113 kg
---------------------	--------