

Comfort Relief

Features

- The open back design gives the ability to off-load the pressure typically seen on spinal deformities through the back shell opening
- With the BOA® Fit System, equal and stable pressure distribution is accomplished in an “off the shelf” back design
- Updated cover design and 3-layer foam construction
- Cover options: Diamond Spacer Standard & Comfort-Tek Optional
- Compatible with BodiLink® Compass and Power Mount Hardware
- Warranty: 2 years back and hardware, 6 months for cover



Standard sizes

WIDTH

HEIGHT

	35,5 cm	40,5 cm	45,5 cm	51,0 cm
51,0 cm	35,5 x 51,0	40,5 x 51,0	45,5 x 51,0	51,0 x 51,0
45,5 cm	35,5 x 45,5	40,5 x 45,5	45,5 x 45,5	51,0 x 45,5
40,5 cm	35,5 x 40,5	40,5 x 40,5	45,5 x 40,5	51,0 x 40,5
35,5 cm	35,5 x 35,5	40,5 x 35,5	45,5 x 35,5	
30,5 cm	35,5 x 30,5	40,5 x 30,5	45,5 x 30,5	
25,5 cm	35,5 x 25,5	40,5 x 25,5	45,5 x 25,5	

Compass 4 Hardware offer 5 cm width adjustments

Dimensions

Width	35,5 – 51,0 cm
Height	25,5 – 51,0 cm

Weight capacity

35,5 – 45,5 cm wide	113 kg
51,0 cm wide	136 kg