

Comfort Unite 3

Features

- 17 cm of shell contour
- Adjustable moldable stays to conform and support the back
- Keyhole and stud strap securement for secure attachment of the cover to the shell
- Lightweight, durable aluminum shell providing pelvic and scapular relief
- ROHO insert option along the spine and in laterals
- Hardware quick install: easy adjustment for depth, height and angle while user is in chair
- Updated cover design and 2-layer foam construction
- Cover options: Diamond Spacer Standard, Comfort-Tek and GlideWear as Optional
- Compatible with BodiLink® Compass, Power Mount and Direct Mount Hardware
- Warranty: 2 years back and hardware, 6 months for cover



Standard sizes

WIDTH

HEIGHT

	35,5 cm	40,5 cm	45,5 cm	51,0 cm
51,0 cm	35,5 x 51,0	40,5 x 51,0	45,5 x 51,0	51,0 x 51,0
45,5 cm	35,5 x 45,5	40,5 x 45,5	45,5 x 45,5	51,0 x 45,5
40,5 cm	35,5 x 40,5	40,5 x 40,5	45,5 x 40,5	51,0 x 40,5
35,5 cm	35,5 x 35,5	40,5 x 35,5	45,5 x 35,5	51,0 x 35,5

Compass 4 Hardware offer 5 cm width adjustments

with ROHO insert

WIDTH

HEIGHT

	40,5 cm	45,5 cm	51,0 cm
51,0 cm	40,5 x 51,0	45,5 x 51,0	51,0 x 51,0
45,5 cm	40,5 x 45,5	45,5 x 45,5	51,0 x 45,5
40,5 cm	40,5 x 40,5	45,5 x 40,5	51,0 x 40,5

Dimensions*

Weight capacity

Width	35,5 – 51,0 cm	35,5 – 45,5 cm wide	113 kg
Height	35,5 – 51,0 cm	51,0 cm wide	136 kg
Contours	17,0 cm		

*Dimensions with ROHO insert: width 40,5–51,0 cm and height 40,5–51,0 cm