



PDG Bentley LT - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

***Required field**

Special considerations that need to be addressed here (e.g., diagnosis):

Wheelchair & Frame

Bentley LT 20° Manual Tilt Wheelchair *Lightweight aluminum frame. Aluminum frame not recommended for agitated users. Consider standard Bentley with steel frame option for agitated or heavier users. See separate order form.*

\$

Frame Type | All frames offer 300lbs (136kg) weight capacity.

Fixed Width Frame Select a fixed seat width below.

No charge

Adjustable Width Frame Select an adjustable seat width below.

\$

User Weight | Required for setup of tilt uplift force.

User Weight Range:

Up to 75lbs (34kg)

75lbs - 100lbs (34kg - 45kg)

100lbs - 115lbs (45kg - 52kg)

115lbs - 150lbs (52kg - 68kg)

150lbs - 200lbs (68kg - 90kg)

200lbs - 250lbs (90kg - 113kg)

250lbs - 300lbs (113kg - 136kg)

Frame Options

Transit Tie Down System (TTDS) *Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See owners manual for complete details.*

No charge

Seat Measurements

Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

Fixed:	15"	16"	17"	18"	19"	No charge
	20"	21"	22"	23"		\$
	24"					\$

Adjustable 12"-16":	12"	13"	14"	15"	16"	No charge
---------------------	-----	-----	-----	-----	-----	-----------

Adjustable 15"-20":	15"	16"	17"	18"	19"	No charge
---------------------	-----	-----	-----	-----	-----	-----------

	20"					\$
--	-----	--	--	--	--	----

Seat Depth

Select Depth:	16"	17"	18"	19"	No charge
	20"	21"			\$
	22"				\$

16" - 20" seat depths are adjustable.

Frame Finishes

Select Color:	High Gloss Black	Red Pepper	Neon Pink	Ocean Blue Metallic	Tangerine Metallic
	Acid Green Metallic	Midnight Blue Pearl	White Pearl	Black Lava	Titanium Grey
	Pomegranate Red	Blue Mist	Indigo Purple	Hyper Blue	Flame Orange
	OD Green	Silver Streak	Lavendar Dusk	Black Galaxy	Seafoam Green
	Daffodil Yellow				

Armrests

Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests \$

8"-11" Range 10-14" Range 12"-16" Range

Fixed Height Armrests 12" height, full length arm pads only. No charge

Armrest Position: Rear (Default) Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Armpads

Regular Armpads No charge

Desk Length 2" x 10" Full Length 2" x 14"

Waterfall Armpads \$

Desk Length 2" x 10" Full Length 2" x 14"

Gel Armpads May interfere with wheels in some tilt configurations. \$

Narrow 2" x 12" Wide 3.5" x 12"

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change. May cause tilt interference.

Armrest Wideners, 1" Total Additional Width \$

Armrest Wideners, 2" Total Additional Width \$

Clothing Guards

Clothing Guards, Single Post Armrests Not compatible with 8"-11" height range. \$

Armrest Extensions | Not compatible with flip-up armrests or desk length arm pads.

Armrest Extensions - 3" Forward \$

Armrest Extensions - 3" Rearward \$

Seat Height

Seat Height - 6" Casters

Low						High					
13"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	20"
Seat to Floor Height											

Seat Height - 8" Casters

Low				High			
15"	16"	17"	18"	17"	18"	19"	20"
Seat to Floor Height							

Seat Pan

No Seat Pan No charge

Solid Aluminium Seat Pan \$

Front Casters

Urethane Casters

6" x 1.25" Caster	No charge
8" x 1.25" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	No charge
6" x 1.50" Caster	\$
8" x 2.00" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Pneumatic Casters

6" x 1.25" Pneumatic Caster	\$
8" x 1.25" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$
8" x 2.00" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster <i>Urethane tire and silver anodized aluminum hub.</i>	\$
---	----

Frog Legs Suspension Forks

FrogLegs Suspension Forks <i>Only compatible with 16" and 24" rear wheels. Requires 6" soft-roll caster.</i>	\$
--	----

Rear Wheels

Composite Rear Wheel	No charge
12" 16" 20" 22" 24"	

Spoke Rear Wheel	\$
20" 22" 24"	

Rear Wheel Position: Rear Front (Default)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires | *Pneumatic tires may result in a finished seat height 1/2" greater than stated.*

Urethane Tires	No charge
Pneumatic Tires	\$
Pneumatic Tires with Flat-Free Inserts <i>Not available with spoke wheels.</i>	\$

Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$
Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Attendant Wheel Locks <i>Not compatible with wheel lock extensions.</i>	\$
Push to Lock Pull to Lock	
User and Attendant Wheel Locks <i>Not compatible with tall back posts (25") or drum brakes on 20" rear wheels.</i>	\$
Push to Lock Pull to Lock	
Drum Brakes <i>Available with 12" composite or 20" spoke rear wheels only. Urethane tires only. Not compatible with plastic coated handrims or two trigger tilt system.</i>	\$
User Activated Unilateral Wheel Lock <i>Pull to lock. Available up to 20" seat width only. Not compatible with hill holders.</i>	\$
Left Activation Right Activation	

Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	\$
Hill Holders <i>Requires pneumatic rear tires with flat-free inserts. Not compatible with unilateral wheel lock.</i>	\$

Anti-Tippers

Rear Anti-Tippers *Required safety option.*

\$

Backrest

Stroller Handle Back

Back Post Height: Short (22") Tall (25")

Measured from seat pan to handle.

Folding Style: Fold-Down Non-Folding

Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.

Seat to Back Angle: Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

Back angles >98° in combination with seat depths >19" may limit range of tilt.

Backrest Options

Angle/Height Adjustable Stroller Handle:

\$

Button on Outside Button on Inside

Endomorph Back *Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with dynamic rocker back or reclining backrest. Requires straight back post.*

\$

1" Offset Per Side 1.5" Offset Per Side 2" Offset Per Side

Spring Back Suspension *Not compatible with dynamic rocker back or reclining backrest.*

\$

Seating Dynamics Dynamic Rocker Back Interface *Not compatible with fold-down back, spring back suspension, endomorph back or 30° recline. Not available with 20", 21" and 22" seat depths. Finished back post height is 4" taller than selected back post height.*

\$

Recline

Reclining Backrest (Manual) *30° maximum recline. Not compatible with tall (25") back posts, fold-down back, dynamic rocker back or endomorph back. Requires longer frame. Maximum seat depth is 19". Finished back post height is 3" taller than selected back post height.*

\$

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post

No charge

Attendant Tilt Only, Left Back Post

No charge

Two Trigger Setup *A trigger mounted on each side of the stroller handle. Reduces grip force required to activate tilt. Not compatible with reclining backrest or drum brakes.*

No charge

User Tilt Only | Full length armpads only.

User Tilt Only, Right Armrest

\$

User Tilt Only, Left Armrest

\$

Attendant & User Tilt | Not compatible with tall (25") back posts or reclining backrest. Full length armpads only.

Attendant & User Tilt Operation

\$

Attendant Trigger Position: Right Back Post Left Back Post

User Trigger Position: Right Armrest Left Armrest

Power Tilt | 250lbs (113kg) weight capacity. 20° tilt range. Includes 2 batteries for off-chair charging. Power tilt is compatible with either O2 holder or IV pole, not both.

Standard Power Tilt

\$

Power Tilt with Manual Tilt Override *Allows for instant manual release of tilt in case of power tilt failure. Minimum seat width is 16".*

\$

Toggle Switch Control

No charge

Installed on Right Installed on Left

No charge

Button Switch Control

No charge

Installed on Right Installed on Left

No charge

Stealth Egg Switch Control *Not installed.*

No charge

Stealth Headrest Swing-Away Mounting Hardware *Only for Stealth egg switch. Not compatible with all headrest makes and models.*

\$

Front Rigging

No Front Rigging

No charge

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 70°

\$

Universal Adjustable Angle (0°-100°) Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.

\$

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle For shorter lower leg lengths of 8"-13".

\$

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

Swing Away 60°

No charge

Swing Away 70°

No charge

Swing Away Elevating Legrests (ELR) Includes calf pads.

\$

1" Offset Hangers for Swing Away Front Rigging

\$

Plug In Front Rigging | Legacy design. Not compatible with Easy-Out.

Plug In 70°

No charge

Footplates

Composite Footplates

No charge

Aluminium Angle Adjustable Footplates

\$

One Piece Fixed Footplate Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates.

\$

One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates.

\$

Interlocking Angle Adjustable Footplates Available for seat widths up to 21"

\$

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short

No charge

Standard Extension Tubes, Long

\$

Footplate Height:

Specify:

Ext	Short	8" - 13"	11" - 16"	12" - 16"	13" - 16"
	Long	N/A	14" - 22"	14" - 19"	14" - 17"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
Legrest					

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging Not compatible with gel padding. Includes reinforced premium calf pad mount.

\$

Calf Pads, Swing Away Front Rigging Not compatible with gel padding.

\$

Gel Padding Not compatible with calf pads or Easy-Out front rigging.

\$

Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging

\$

Gel Knee Pads, Swing Away ELR Front Rigging

\$

Length Adjustable Heel Loops

\$

Standard Calf Strap

\$

Positioning Belts

Q'Straint Lap Belt Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.

\$

Miscellaneous Options

One-Arm Drive Available with 22" and 24" composite rear wheels. Minimum 17" STF and 18" STF respectively. Large hand setup includes 24" rear wheels with 22" driving handrim. The maximum achievable tilt could be as low as 5°-10° and be further reduced to 0° when transit tie down system is installed. Check with PDG prior to ordering. Not compatible with user and attendant wheel locks, drum brakes, quick release rear axles or pneumatic tires.

\$

Installed On: Right Left

Handrim Setup: Standard Large Hand

Large hand setup only available with 24" rear wheels.

Spoke Guards Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.

\$

Black Frame color

Gel Hip Guides Reduces available seat width by 2"

\$

IV Pole Not compatible with attendant wheel locks.

\$

Installed on Right Installed on Left

O2 Holder Not compatible with attendant wheel locks.

\$

Installed on Right Installed on Left

Medical Necessities Backpack Black with Permobil logo.

\$

Rubber Shroud for Caster Journals

\$

Adhesive Back Hook Strip For cushion mounting, 2"x30".

\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes