



PDG Bentley LT-R - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

Wheelchair & Frame

Bentley LT-R 20° Manual Tilt Wheelchair Base Model *Lightweight aluminum frame and reverse wheel layout for improved self-propulsion. Choose standard Bentley with steel frame for agitated or heavier users.*

\$

Frame Type | All frames offer 300lbs (136kg) weight capacity.

Fixed Width Frame Select a fixed seat width below.

No charge

Adjustable Width Frame Select an adjustable seat width below.

\$

User Weight | Required for setup of tilt uplift force.

User Weight Range:

Up to 75lbs (34kg)

75lbs - 100lbs (34kg - 45kg)

100lbs - 115lbs (45kg - 52kg)

115lbs - 150lbs (52kg - 68kg)

150lbs - 200lbs (68kg - 90kg)

200lbs - 250lbs (90kg - 113kg)

250lbs - 300lbs (113kg - 136kg)

Frame Options

Transit Tie Down System (TTDS) *Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details.*

No charge

Seat Measurements

Seat Width | Overall chair width is +10.5".

Fixed: 15" 16" 17" 18" 19" No charge

20" 21" 22" 23"

\$

24"

\$

Adjustable 12"-16": 12" 13" 14" 15" 16" No charge

Adjustable 15"-20": 15" 16" 17" 18" 19" No charge

20"

\$

Seat Depth

Select Depth: 16" 17" 18" 19" No charge

20" 21"

\$

22"

\$

16" - 20" seat depths are adjustable.

Frame Finishes

Select Color:

High Gloss Black
Acid Green Metallic
Pomegranate Red
OD Green
Daffodil Yellow

Red Pepper
Midnight Blue Pearl
Blue Mist
Silver Streak

Neon Pink
White Pearl
Indigo Purple
Lavendar Dusk

Ocean Blue Metallic
Black Lava
Hyper Blue
Black Galaxy

Tangerine Metallic
Titanium Grey
Flame Orange
Seafoam Green

Armrests

Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests \$

8"-11" Range 10-14" Range 12"-16" Range

Fixed Height Armrests 12" height, full length armpads only. No charge

Armrest Position: Rear (Default) Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Armpads

Regular Armpads No charge

Desk Length 2" x 10" Full Length 2" x 14"

Waterfall Armpads \$

Desk Length 2" x 10" Full Length 2" x 14"

Gel Armpads May interfere with wheels in some tilt configurations. \$

Narrow 2" x 12" Wide 3.5" x 12"

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change. May cause tilt interference.

Armrest Wideners, 1" Total Additional Width \$

Armrest Wideners, 2" Total Additional Width \$

Clothing Guards

Clothing Guards, Single Post Armrests Not compatible with 8"-11" height range. \$

Armrest Extensions | Not compatible with flip-up armrests or desk length armpads.

Armrest Extensions - 3" Forward \$

Armrest Extensions - 3" Rearward \$

Seat Height

Seat Height - 6" Casters

Low						High					
12"	13"	14"	15"	16"	17"	14"	15"	16"	17"	18"	19"
Seat to Floor Height											

Seat Height - 8" Casters

Low				High			
14"	15"	16"	17"	16"	17"	18"	19"
Seat to Floor Height							

Seat Pan

No Seat Pan No charge

Solid Aluminium Seat Pan \$

Rear Casters

Urethane Casters

6" x 1.25" Caster	No charge
8" x 1.25" Caster <i>Only compatible with 24" rear wheels.</i>	No charge
6" x 1.50" Caster	\$
8" x 2.00" Caster <i>Only compatible with 24" rear wheels.</i>	\$

Pneumatic Casters

6" x 1.25" Pneumatic Caster	\$
8" x 1.25" Pneumatic Caster <i>Only compatible with 24" rear wheels.</i>	\$
8" x 2.00" Pneumatic Caster <i>Only compatible with 24" rear wheels.</i>	\$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster <i>Urethane tire and silver anodized aluminum hub.</i>	\$
---	----

Frog Legs Suspension Forks

FrogLegs Suspension Forks <i>Only compatible with 24" rear wheels. Requires 6" soft-roll caster.</i>	\$
--	----

Drive Wheels

Composite Drive Wheel	No charge
20" 22" 24"	

Spoke Drive Wheel	\$
20" 22" 24"	

Drive Wheel Position: Rear (Less Stable) Mid (Default) Front (More Stable)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference. Front position not compatible with 20"-22" seat depths.

Axles

Quick Release Axles	No charge
Fixed Axles	No charge

Tires

Urethane Tires	No charge
Pneumatic Tires <i>May result in a finished seat height 1/2" greater than stated.</i>	\$
Pneumatic Tires with Flat-Free Inserts <i>May result in a finished seat height 1/2" greater than stated. Not available with spoke wheels.</i>	\$

Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$
Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Drum Brakes <i>Available with 20" spoke rear wheels only. Urethane tires only. Not compatible with plastic coated handrims or two trigger tilt system.</i>	\$
User Activated Unilateral Wheel Lock <i>Pull to lock. Available up to 20" seat widths only. Not compatible with hill holders.</i>	\$
Left Activation Right Activation	\$

Wheel Lock Options

Wheel Lock Extensions	\$
Hill Holders <i>Requires pneumatic rear tires with flat-free inserts. Not compatible with unilateral wheel lock.</i>	\$

Anti-Tippers

Must select one option.

Front Anti-Tippers, Standard Required safety option. \$

Front Anti-Tippers, Spring Loaded Required safety option. Dampens the forward-tip movement of the chair during transfers and allows for improved ground clearance. Not removable. \$

Backrest

Stroller Handle Back

Back Post Height: Short (22") Tall (25")
Measured from seat pan to handle.

Folding Style: Fold-Down Non-Folding
Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.

Seat to Back Angle: Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°
Back angles >98° in combination with seat depths >19" may limit range of tilt.

Backrest Options

Angle/Height Adjustable Stroller Handle: Not compatible with headrest mounting plate. \$

Button on Outside Button on Inside

Endomorph Back Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with dynamic rocker back or reclining backrest. Requires straight back post. \$

1" Offset Per Side 1.5" Offset Per Side 2" Offset Per Side

Spring Back Suspension Not compatible with reclining backrest or dynamic rocker back. \$

Seating Dynamics Dynamic Rocker Back Interface Not compatible with fold-down back, spring back suspension, endomorph back or 30° recline. Not available with 20", 21" and 22" seat depths. Finished back post height is 4" taller than selected back post height. \$

Recline

Reclining Backrest (Manual) 30° maximum recline. Not compatible with fold-down back, dynamic rocker back or endomorph back. Requires longer frame. Maximum seat depth is 19". Finished back post height is 3" taller than selected back post height. \$

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post No charge

Attendant Tilt Only, Left Back Post No charge

Two Trigger Setup A trigger mounted on each side of the stroller handle. Reduces grip force required to activate tilt. Not compatible with reclining backrest or drum brakes. No charge

User Tilt Only | Full length armpads only.

User Tilt Only, Right Armrest \$

User Tilt Only, Left Armrest \$

Attendant & User Tilt | Not compatible with tall (25") back posts or reclining backrest. Full length armpads only.

Attendant & User Tilt Operation \$

Attendant Trigger Position: Right Back Post Left Back Post

User Trigger Position: Right Armrest Left Armrest

Power Tilt | 250lbs (113kg) weight capacity. 20° tilt range. Includes 2 batteries for off-chair charging. Power tilt is compatible with either O2 holder or IV pole, not both. Not available with 12" STF height.

Standard Power Tilt \$

Power Tilt with Manual Tilt Override Allows for instant manual release of tilt in case of power tilt failure. Minimum seat width is 16" \$

Toggle Switch Control No charge

Installed on Right Installed on Left No charge

Button Switch Control No charge

Installed on Right Installed on Left No charge

Stealth Egg Switch Control Not installed. No charge

Stealth Swing-Away Mounting Hardware Only for Stealth egg switch. Not compatible with all headrest makes and models. \$

Front Rigging

No Front Rigging

No charge

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 70°

\$

Universal Adjustable Angle (0°-100°) Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.

\$

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle For shorter lower leg lengths of 8"-12.5".

\$

Plug In Front Rigging | Legacy design. Not compatible with Easy-Out.

Plug In 70°

No charge

Footplates

Composite Footplates

No charge

Aluminium Angle Adjustable Footplates

\$

One Piece Fixed Footplate Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates.

\$

One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates.

\$

Interlocking Angle Adjustable Footplates Available for seat widths up to 21"

\$

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short

No charge

Standard Extension Tubes, Long

\$

Footplate Height:

Specify:

Ext	Short	8" - 13"	11" - 16"	12" - 16"	13" - 16"
	Long	N/A	14" - 22"	14" - 19"	14" - 17"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
Legrest					

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging Includes reinforced premium calf pad mount.

\$

Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging

\$

Length Adjustable Heel Loops

\$

Standard Calf Strap

\$

Positioning Belts

Q'Straint Lap Belt Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.

\$

Miscellaneous Options

Spoke Guards Available for 20"-24" composite drive wheels. Made from durable and lightweight aluminium.

\$

Black Frame Color

\$

Gel Hip Guides Reduces available seat width by 2"

\$

IV Pole Not compatible with attendant wheel locks.

\$

Installed on Right Installed on Left

\$

O2 Holder Not compatible with attendant wheel locks.

\$

Installed on Right Installed on Left

\$

Medical Necessities Backpack Black with Permobil logo.

\$

Adhesive Back Hook Strip For cushion mounting, 2"x30".

\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes