



PDG Bentley - Order Form



Dealer Name*:		Account #*:
City, Post Code, Country:		
Phone:	Email:	P.O. Number
Clinician:	Client Reference*:	Client Weight:
Rep Name:		
*Required field Special considerations that need to be addressed here (e.g., diagnosis):		

Wheelchair & Frame

Bentley 20° Manual Tilt Wheelchair Also available: Bentley LT with lightweight aluminium frame. See separate order form.

\$

Frame Type

250lbs (113kg) Weight Capacity	Available for 14" - 22" seat widths.	No charge
300lbs (136kg) Weight Capacity	Available for 14" - 22" seat widths.	\$
450lbs (204kg) Weight Capacity	Available for 20" - 32" seat widths.	\$

User Weight | Required for setup of tilt uplift force.

User Weight Range:	Up to 75lbs (34kg)	75lbs - 100lbs (34kg - 45kg)	100lbs - 115lbs (45kg - 52kg)
	115lbs - 150lbs (52kg - 68kg)	150lbs - 200lbs (68kg - 90kg)	200lbs - 250lbs (90kg - 113kg)
	250lbs - 300lbs (113kg - 136kg)	300lbs - 375lbs (136kg - 170kg)	375lbs - 450lbs (170kg - 204kg)

Frame Options

Extended Base with Rear Floor Jacks Lifts rear wheels off the ground for increased stability. Limited ground clearance. Not compatible with attendant wheel lock, HD anti-tipper, IV pole, O2 holder.

\$

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See owners manual for complete details. Not compatible with seat widths > 22". Tested to maximum weight capacity of 300lbs (136kg).

No charge

Seat Measurements

Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

Select Width:	14"	15"	16"	17"	18"	19"	No charge
	20"	21"	22"	23"			\$
	24"	25"	26"	27"			\$
	28"	29"	30"	31"	32"		\$

Seat widths 23"-32" require 450lb (204kg) Frame. Extended delivery lead times may apply.

Seat Depth

Select Depth:	16"	17"	18"	19"	No charge
	20"	21"			\$
	22"				\$

Frame Finishes

Select Color:	High Gloss Black	Red Pepper	Neon Pink	Ocean Blue Metallic	Tangerine Metallic
	Acid Green Metallic	Midnight Blue Pearl	White Pearl	Black Lava	Titanium Grey
	Pomegranate Red	Blue Mist	Indigo Purple	Hyper Blue	Flame Orange
	OD Green	Silver Streak	Lavendar Dusk	Black Galaxy	Seafoam Green
	Daffodil Yellow				

Armrests

Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests \$

8"-11" Range 10-14" Range 12"-16" Range

Dual Post, Height Adjustable Armrests \$

8"-11" Range 10-14" Range 12"-16" Range

Fixed Height Armrests 12" height, full length armpads only. No charge

Single Post Dual Post

Armrest Position: Rear Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Locking Flip-Up Armrests Not recommended for users with excessive strength or tone. Not compatible with fold-down back, or clothing guards. Only with regular armpads. \$

Set at 10" Set at 11" Set at 12" Set at 13" Set at 14"

Armpads

Regular Armpads No charge

Desk Length 2" x 10" Full Length 2" x 14"

Waterfall Armpads \$

Desk Length 2" x 10" Full Length 2" x 14"

Gel Armpads May interfere with wheels in some tilt configurations. \$

Narrow 2" x 12" Wide 3.5" x 12"

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change. May cause tilt interference.

Armrest Wideners, 1" Total Additional Width \$

Armrest Wideners, 2" Total Additional Width \$

Clothing Guards | Not compatible with flip-up armrests.

Clothing Guards, Single Post Armrests Not compatible with 8"-11" height range. \$

Clothing Guards, Dual Post Armrests Not compatible with 8"-11" height range. \$

Armrest Extensions | Not compatible with flip-up armrests or desk length armpads.

Armrest Extensions - 3" Forward \$

Armrest Extensions - 3" Rearward \$

Seat Height

Seat Height - 6" Casters

RW	Low					High					
	24"										
	22"										
	20"										
	16"										
	12"										
	13"	14"	15"	16"	17"	15"	16"	17"	18"	19"	20"
Seat to Floor Height											

Seat Height - 8" Casters

RW	Low			High			
	24"						
	16"						
	15"	16"	17"	17"	18"	19"	20"
Seat to Floor Height							

Seat Pan

No Seat Pan	No charge
Solid Aluminium Seat Pan	\$

Front Casters

Urethane Casters

6" x 1.25" Caster	No charge
8" x 1.25" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	No charge
6" x 1.50" Caster	\$
8" x 2.00" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Pneumatic Casters

6" x 1.25" Pneumatic Caster	\$
8" x 1.25" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$
8" x 2.00" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster <i>Urethane tire and silver anodized aluminum hub.</i>	\$
---	----

Frog Legs Suspension Forks

FrogLegs Suspension Forks <i>Only compatible with 16" and 24" rear wheels. Requires 6" soft-roll caster.</i>	\$
--	----

Rear Wheels

Composite Rear Wheel	No charge
12" 16" 20" 22" 24"	

Spoke Rear Wheel	\$
------------------	----

20" 22" 24"

Heavy Duty Spoke Wheel	\$
------------------------	----

24"

Rear Wheel Position: Rear Front (Default)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires | *Pneumatic tires may result in a finished seat height 1/2" greater than stated.*

Urethane Tires	No charge
Pneumatic Tires	\$
Pneumatic Tires with Flat-Free Inserts <i>Not available with spoke wheels.</i>	\$

Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$
Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Attendant Wheel Locks <i>Not compatible with wheel lock extensions.</i>	\$
Push to Lock Pull to Lock	
User and Attendant Wheel Locks <i>Not compatible with drum brakes on 20" rear wheels.</i>	\$
Push to Lock Pull to Lock	
Drum Brakes <i>Available with 12" composite or 20" spoke rear wheels only. Urethane tires only. Not compatible with plastic coated handrims or two trigger tilt system.</i>	\$
User Activated Unilateral Wheel Lock <i>Pull to lock. Available up to 20" seat widths only. Not compatible with hill holders.</i>	\$
Left Activation Right Activation	

Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	\$
Hill Holders <i>Requires pneumatic rear tires with flat-free inserts. Not compatible with unilateral wheel lock.</i>	\$

Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	\$
Front Anti-Tippers	\$
Heavy Duty Rear Anti-Tippers <i>Provides additional rearward stability for users with excessive movement. Not compatible with attendant wheel locks, rear floor jacks, IV pole or O2 holder.</i>	\$

Backrest

Back Post Style:	Stroller Handle Back Cane Handle Back	
	<i>Seat widths greater than 22" require stroller handle back. Cane handle back may reduce tilt range. Cane handle back not compatible with power tilt with manual override, dynamic rocker back or endomorph back.</i>	
Back Post Height:	Short (22") Tall (25")	
	<i>Measured from seat pan to handle.</i>	
Folding Style:	Fold-Down Non-Folding	
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>	
Seat to Back Angle:	Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°	
	<i>Back angles >98° in combination with seat depths >19" may limit range of tilt.</i>	

Backrest Options

Angle/Height Adjustable Stroller Handle: <i>Not compatible with cane handle back.</i>	\$
Button on Outside Button on Inside	
Endomorph Back <i>Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with flip-up armrests, dynamic rocker back, reclining backrest or cane handle back. Requires straight back post.</i>	\$
1" Offset Per Side 1.5" Offset Per Side 2" Offset Per Side	
Spring Back Suspension <i>Not compatible with reclining back or dynamic rocker back.</i>	No charge
Seating Dynamics Dynamic Rocker Back Interface <i>Not compatible with fold-down back, spring back suspension, endomorph back, 30° recline or cane handle. Not available with 20", 21" and 22" seat depths. Finished back post height is 4" taller than selected back post height.</i>	\$

Recline

Reclining Backrest (Manual) <i>30° maximum recline. Not compatible with cane handle back, spring back suspension, fold-down back, dynamic rocker back or endomorph back. Requires longer frame. Maximum seat depth is 19". Finished back post height is 3" taller than selected back post height.</i>	\$
--	----

Tilt Activation

Attendant Tilt Only | Chairs with 450lbs (204kg) weight capacity require two tilt triggers for tilt activation.

Attendant Tilt Only, Right Back Post	No charge
Attendant Tilt Only, Left Back Post	No charge
Two Trigger Setup A trigger mounted on each side of the stroller handle. Reduces grip force required to activate tilt. Not compatible with reclining backrest or drum brakes.	No charge

User Tilt Only | Not available for 450lbs (204kg) weight capacity. Full length arm pads only.

User Tilt Only, Right Armrest	\$
User Tilt Only, Left Armrest	\$

Attendant & User Tilt | Not available for 450lbs (204kg) weight capacity. Not compatible with tall (25") back posts or reclining backrest. Full length arm pads only.

Attendant & User Tilt Operation	\$
---------------------------------	----

Attendant Trigger Position:	Right Back Post	Left Back Post
User Trigger Position:	Right Armrest	Left Armrest

Power Tilt | 250lbs (113kg) weight capacity. 20° tilt range. Includes 2 batteries for off-chair charging. Power tilt is compatible with either O2 holder or IV pole, not both. Not available with 21" and 22" seat depths. Requires stroller handle back.

Standard Power Tilt	\$
Power Tilt with Manual Tilt Override Allows for instant manual release of tilt in case of power tilt failure. Not compatible with cane handle back.	\$
Toggle Switch Control	No charge
Installed on Right	No charge
Installed on Left	No charge
Button Switch Control	No charge
Installed on Right	No charge
Installed on Left	No charge
Stealth Egg Switch Control Not installed.	No charge
Stealth Swing-Away Mounting Hardware Only for Stealth egg switch. Not compatible with all headrest makes and models.	\$

Front Rigging

No Front Rigging	No charge
------------------	-----------

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 70°	\$
Universal Adjustable Angle (0°-100°) Select angle. Angle adjustment by loosening one bolt in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.	\$
0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°	
High Mount Option for Easy-Out Adjustable Angle For shorter lower leg lengths of 8"-12.5".	\$

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

Swing Away 60°	No charge
Swing Away 70°	No charge
Swing Away Elevating Legrests (ELR) Includes calf pads.	\$
1" Offset Hangers for Swing Away Front Rigging	\$

Plug In Front Rigging | Legacy design. Not compatible with Easy-Out.

Plug In 70°	No charge
-------------	-----------

Footplates

Composite Footplates	No charge
Aluminium Angle Adjustable Footplates	\$
One Piece Fixed Footplate Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates.	\$
One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates.	\$
Interlocking Angle Adjustable Footplates Available for seat widths up to 21"	\$

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short	No charge
Standard Extension Tubes, Long	\$
Heavy Duty Extension Tubes, Short	\$
Heavy Duty Extension Tubes, Long	\$

Footplate Height: Specify:

Ext Tube	Std Short	8" - 13"	11" - 16"	12" - 16"	13" - 16"
	Std Long	N/A	14" - 22"	14" - 19"	14" - 17"
	HD Short	N/A	N/A	N/A	13" - 17"
	HD Long	N/A	N/A	N/A	15" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
Legrest					

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging	Not compatible with gel padding. Includes reinforced premium calf pad mount.	\$
Calf Pads, Swing Away Front Rigging	Not compatible with gel padding.	\$
Gel Padding	Not compatible with calf pads or Easy-Out front rigging.	\$
Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging		\$
Gel Knee Pads, Swing Away ELR Front Rigging		\$
Length Adjustable Heel Loops		\$
Standard Calf Strap		\$
Miller's Dynamic Footrest Mounts, Plug In Front Rigging	Not compatible with Easy-Out front rigging.	\$
Miller's Dynamic Footrest Mounts, Swing Away Front Rigging	Not compatible with Easy-Out front rigging. Not available with 13" STF height.	\$

Positioning Belts

Q'Straint Lap Belt	Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.	\$
--------------------	--	----

Miscellaneous Options

One-Arm Drive	Available with 22" and 24" composite rear wheels. Minimum 17" STF and 18" STF respectively. Large hand setup includes 24" rear wheels with 22" driving handrim. The maximum achievable tilt could be as low as 5°-10° and be further reduced to 0° when transit tie down system is installed. Check with PDG prior to ordering. Not compatible with user and attendant wheel locks, drum brakes, quick release rear axles or pneumatic tires.	\$
Installed On:	Right Left	
Handrim Setup:	Standard Large Hand	
Large hand setup only available with 24" rear wheels.		
Spoke Guards	Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.	\$
Black Frame color		
Gel Hip Guides	Reduces available seat width by 2"	\$
IV Pole	Not compatible with heavy duty rear anti tippers, rear floor jacks or attendant wheel locks.	\$
Installed on Right Installed on Left		\$
O2 Holder	Not compatible with heavy duty rear anti tippers, rear floor jacks or attendant wheel locks.	\$
Installed on Right Installed on Left		\$
Medical Necessities Backpack	Black with Permobil logo.	\$
Adhesive Back Hook Strip	For cushion mounting, 2"x30".	\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes