



PDG Eclipse - Order Form



Dealer Name*:		Account #*:	
City, Post Code, Country:			
Phone:	Email:	P.O. Number	
Clinician:	Client Reference*:	Client Weight:	
Rep Name:			
*Required field Special considerations that need to be addressed here (e.g., diagnosis):			

Wheelchair & Frame

Eclipse Bariatric Manual Wheelchair

\$

Frame Type

600lbs (272kg) Weight Capacity

No charge

1000lbs (453kg) Weight Capacity *Not compatible with Transit Tie Down System*

\$

Frame Options

Transit Tie Down System (TTDS) TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. weight capacity of 600lbs (272kg). Not compatible with O2 holder or IV pole.

No charge

Seat Measurements

Seat Width | Overall chair width is +7.5"

Select Width:	20" Seat Width	22" Seat Width			\$
	24" Seat Width	26" Seat Width			\$
	28" Seat Width	30" Seat Width	32" Seat Width	36" Seat Width	\$
	40" Seat Width	44" Seat Width	48" Seat Width		\$

Seat Depth (Adjustable)

Select Depth:	16" Seat Depth	17" Seat Depth	18" Seat Depth	19" Seat Depth	
	20" Seat Depth	21" Seat Depth			\$
	22" Seat depth	23" Seat Depth	24" Seat Depth		\$

24" seat depth not available with TTDS.

Frame Finishes

Select Color:	High Gloss Black	Red Pepper	Neon Pink	Ocean Blue Metallic	Tangerine Metallic
	Acid Green Metallic	Midnight Blue Pearl	White Pearl	Black Lava	Titanium Grey
	Pomegranate Red	Blue Mist	Indigo Purple	Hyper Blue	Flame Orange
	OD Green	Silver Streak	Lavendar Dusk	Black Galaxy	Seafoam Green
	Daffodil Yellow				

Armrests

Armrest Type

No Armrests	No charge
Dual Post, Height Adjustable Armrests	\$
10-14" Range 12"-16" Range	
Dual Post, Fixed Height Armrests (12")	No charge

Armpads

Regular Armpads	No charge
Desk Length 2" x 10" Full Length 2" x 14"	
Gel Armpads <i>May interfere with wheels in some tilt configurations.</i>	\$
Narrow 2" x 12" Wide 3.5" x 12"	

Armrest Options

Forward Guards	\$
Armrest Extensions - 3" Forward <i>Not compatible with desk length armpads.</i>	\$

Seat Height

Seat Height | *Front and rear height can be maximum 1" apart.*

Front Seat Height:	14" STF Front	15" STF Front	16" STF Front	17" STF Front	18" STF Front	
	19" STF Front	20" STF Front				\$
Rear Seat Height:	14" STF Rear	15" STF Rear	16" STF Rear	17" STF Rear	18" STF Rear	
	19" STF Rear	20" STF Rear				\$

Seat Pan

Solid Aluminium Seat Pan	\$
--------------------------	----

Front Casters

5" Casters with Heavy Duty Steel Fork	No charge
---------------------------------------	-----------

Rear Wheels

24" Heavy Duty Spoke Wheels					No charge
Rear Wheel Position:	Position 1 (Rear)	Position 2	Position 3	Position 4 (Front)	
<i>Positions 3 and 4 not compatible with TTDS. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.</i>					

Axles

Quick Release Rear Axles	No charge
--------------------------	-----------

Rear Tires

Urethane Tires	No charge
----------------	-----------

Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims	\$

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Wheel Lock Extensions	\$

Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	\$
--	----

Backrest

Back Post Style:	Stroller Handle Back	Cane Handle Back		
	<i>Cane handle not compatible with headrest mounting plate.</i>			
Back Post Height:	Short (22")	Tall (25")		
	<i>Measured from seat pan to handle.</i>			
Seat to Back Angle:	Set at 85°	Set at 90°	Set at 95°	Set at 100°

Backrest Options

Endomorph Back Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not available with reclining back.

1" Offset Per Side 2" Offset Per Side

\$

Back Supports

Back Upholstery

No Back Upholstery No charge

Tension Adjustable Back Upholstery \$

Position of Back Upholstery

Uph Type	22" Back Post						25" Back Post								
	Regular														
Tall															
	17"	18"	19"	20"	21"	22"	17"	18"	19"	20"	21"	22"	23"	24"	25"
Upholstery Height															

Recline

Reclining Backrest (Manual) 30° maximum recline. Not compatible with endomorph back or cane handle back. Reclining backrest is 1" narrower than seat width. Back frame is not removable. Finished back post height is 3" taller than selected back post height. Not recommended for user weight above 450lbs (204kg). Maximum seat depth is 23".

\$

Front Rigging

No Front Rigging No charge

Center Mount Front Rigging Bracket \$

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 60° \$

Universal Adjustable Angle (0°-60°) Select Angle. Angle adjustment by loosening one bolt. Tool-free leg length adjustment. \$

0° 10° 20° 30° 40° 50° 60°

AA High Mount Option For shorter lower leg lengths of 8"-12.5". \$

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

Swing Away 60° No charge

Swing Away Elevating Legrests (ELR) Includes calf pads. \$

1" Offset Hangers for Swing Away Front Rigging \$

Footplates

Aluminum Angle Adjustable Footplates \$

Extension Tubes | No selection required for Easy-Out AA High-Mount option.

Standard Extension Tubes, Short No charge

Standard Extension Tubes, Long \$

Footplate Height: Specify:

Ext	Short	8" - 13"	11" - 16"	12" - 16"
	Long	N/A	14" - 22"	14" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away
Legrest				

Front Rigging Accessories

Oversized 9"x10" Calf Pad	<i>Not compatible with seat widths less than 25"</i>	\$
Calf Pads	<i>Not compatible with Easy-Out front rigging.</i>	\$
Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging		\$
Gel Knee Pads, Swing Away ELR Front Rigging		\$
Length Adjustable Heel Loops		\$
Standard Calf Strap		\$

Positioning Belts

Q'Straint Lap Belt	<i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	\$
---------------------------	---	----

Miscellaneous Options

IV Pole	<i>Not compatible with TTDS.</i>	\$
O2 Holder	<i>Not compatible with TTDS.</i>	\$
Medical Necessities Backpack	<i>Black with Permobil logo.</i>	\$
Adhesive Back Hook Strip	<i>For cushion mounting, 2"x30".</i>	\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes