



PDG Fuze M50 - Order Form



Dealer Name*:		Account #*:	
City, Post Code, Country:			
Phone:	Email:	P.O. Number	
Clinician:	Client Reference*:	Client Weight:	
Rep Name:			
*Required field Special considerations that need to be addressed here (e.g., diagnosis):			

Wheelchair & Frame

Fuze M50 Mid-Wheel Drive 50° Manual Tilt Wheelchair Articulating rear casters. Drive wheels are in the middle of the chair at the center of gravity. \$

Frame Type

Fixed Width Frame, 300lbs (136kg) Weight Capacity Select a fixed seat width below. No charge

Fixed Width Frame, 400lbs (181kg) Weight Capacity Not compatible with transit tie down system. Requires flat free rear tires or heavy duty spoke rear wheels. Select a fixed seat width below. No charge

Adjustable Width Frame. 250lbs (113kg) Weight Capacity Select an adjustable seat width below. \$

Frame Options

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Tested to maximum weight capacity of 300lb (136kg). No charge

Seat Measurements

Seat Width | Overall chair width is +8.5".

Fixed:	15"	16"	17"	18"		No charge
	19"	20"	21"	22"	23"	No charge
	24"					\$
Adjustable 15"-19":	15"	16"	17"	18"	19"	No charge

Seat Depth

Select Depth:	15"	16"	17"	18"	19"	No charge
	20"	21"				\$
	22"					\$

21"-22" seat depths are not adjustable and require stroller handle back - straight. Not compatible with reclining or fold-down back.

Frame Finishes

Select Color:	High Gloss Black	Red Pepper	Neon Pink	Ocean Blue Metallic	Tangerine Metallic
	Acid Green Metallic	Midnight Blue Pearl	White Pearl	Black Lava	Titanium Grey
	Pomegranate Red	Blue Mist	Indigo Purple	Hyper Blue	Flame Orange
	OD Green	Silver Streak	Lavendar Dusk	Black Galaxy	Seafoam Green
	Daffodil Yellow				

Armrests

Armrest Type

No Armrests																		No charge
Single Post, Height Adjustable Armrests																		\$
10-14" Range	12"-16" Range																	
Fixed Height Armrests	<i>12" height, full length armpads only.</i>																	No charge
Armrest Position:	Mid (Default)	Front																No charge
Locking Flip-Up Armrests	<i>Not recommended for users with excessive strength or tone. Not compatible with fold-down back. Only with regular armpads.</i>																	\$
Set at 10"	Set at 11"	Set at 12"	Set at 13"	Set at 14"														

Armpads

Regular Armpads																		No charge
Desk Length 2" x 10"	Full Length 2" x 14"																	
Waterfall Armpads																		\$
Desk Length 2" x 10"	Full Length 2" x 14"																	
Gel Armpads	<i>May interfere with wheels in some tilt configurations.</i>																	\$
Narrow 2" x 12"	Wide 3.5" x 12"																	

Seat Height

Seat to floor height is adjustable within the range of the chosen axle plate (Low/Standard). Chairs with swing-away front rigging are set up with 5° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out front rigging is not affected.

5" Casters					6" Casters & 6" Forks					6" Casters & 8" Forks					6" Frog Legs Casters & Suspension Forks		
		Low	Low	Low		Low	Low	Low	Std	Low	Low	Low	Std	Std	Low	Low	Std
	Low	Low	Low	Std	Low	Low	Low	Std	Std	Low	Low	Std	Std	Std	Low	Std	Std
Low	Low	Low	Std	Std	Low	Low	Std	Std	Std	Low	Std	Std	Std		Std	Std	Std
13"	14"	15"	16"	17"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	16"	17"	18"
Seat to Floor Height																	

Seat Pan

Seat Pan

No Seat Pan																		No charge
Solid Aluminium Seat Pan																		\$

Casters

Urethane Casters

5" x 1.25" Caster																		No charge
6" x 1.25" Caster																		No charge
6" x 1.50" Caster																		\$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster	<i>Urethane tire and silver anodized aluminum hub.</i>																	\$
-----------------------------------	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	----

Frog Legs Suspension Forks

Frog Legs Suspension Forks	<i>Not compatible with 12" rear wheels. Requires 6" soft-roll caster.</i>																	\$
-----------------------------------	---	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	--	----

Drive Wheels

Composite Rear Wheel	No charge
20" 22" 24"	
Spoke Rear Wheel	\$
20" 22" 24"	
Heavy Duty Spoke Wheel	\$
24"	

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Tires

Urethane Tires	No charge
Pneumatic Tires <i>May result in a finished seat height 1/2" greater than stated.</i>	\$
Pneumatic Tires with Flat-Free Inserts <i>May result in a finished seat height 1/2" greater than stated. Not available with spoke wheels.</i>	\$

Handrims

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$
Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Wheel Lock Extensions	\$

Backrest

Back Post Style:	Stroller Handle Back Cane Handle Back
	<i>Cane handle back not compatible with power tilt with manual override or headrest mounting plate.</i>
Back Post Type:	Straight 8° Bend
Back Post Height:	Short (17") Tall (23")
	<i>Measured from seat pan to handle.</i>
Folding Style:	Fold-Down Non-Folding
	<i>Fold-down back not compatible with reclining backrest.</i>
Seat to Back Angle:	Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle: <i>Not compatible with headrest mounting plate or cane handle back.</i>	\$
Button on Outside Button on Inside	
Endomorph Back <i>Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with flip-up armrests, dynamic rocker back or cane handle back.</i>	\$
1" Offset Per Side 1.5" Offset Per Side 2" Offset Per Side	

Recline

Reclining Backrest (Manual) <i>30° maximum recline. Not compatible with cane handle back or fold-down back. Not available with 21" and 22" seat depths. Finished back post height is 3" taller than selected back post height.</i>	\$
---	----

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post

No charge

Attendant Tilt Only, Left Back Post

No charge

Tilt Limit | Optional. If no tilt limit is selected the full tilt range will be available. See below chart for maximum tilt limits for each rear wheel size and seat to floor height.

Tilt Limit:

Set at 20°

Set at 30°

Set at 40°

Maximum Tilt Limit

RW	20"	30°	40°	50°	50°	50°	50°	
	22"		30°	40°	50°	50°	50°	50°
	24"			30°	40°	50°	50°	50°
		13"	14"	15"	16"	17"	18"	19"
Seat to Floor Height								

Tip Configuration

Front caster lift.

Attendant Tip Only Select trigger mounting location. Omit trigger will allow no front caster lift.

Right Back Post Left Back Post Omit Trigger

Tip Clearance Limit Actual tip may vary by up to 0.5". No selection required if trigger is omitted.

4" Caster Lift (Default) 2" Caster Lift

Front Rigging

No Front Rigging

No charge

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 70°

\$

Universal Adjustable Angle (0°-100°) Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.

\$

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle For shorter lower leg lengths of 8"-12.5".

\$

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out or push to lock wheel locks. Requires pre-tilt of up to 5°.

Swing Away 60°

No charge

Swing Away 70°

No charge

Swing Away Elevating Legrests (ELR) Includes calf pads.

\$

Footplates

Composite Footplates

No charge

Aluminium Angle Adjustable Footplates

\$

One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates.

\$

Interlocking Angle Adjustable Footplates Available for seat widths up to 21"

\$

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short

No charge

Standard Extension Tubes, Long

\$

Ext	Short	8" - 13"	11" - 16"	12" - 16"
	Long	N/A	14" - 22"	14" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away
Legrest				

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging	<i>Not compatible with gel padding. Includes reinforced premium calf pad mount.</i>	\$
Calf Pads, Swing Away Front Rigging	<i>Not compatible with gel padding.</i>	\$
Gel Padding	<i>Not compatible with calf pads or Easy-Out front rigging.</i>	\$
Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging		\$
Gel Knee Pads, Swing Away ELR Front Rigging		\$
Length Adjustable Heel Loops		\$
Standard Calf Strap		\$

Positioning Belts

Q'Straint Lap Belt	<i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	\$
---------------------------	---	----

Miscellaneous Options

Spoke Guards	<i>Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.</i>	\$
Black	Frame Color	
Gel Hip Guides	<i>Reduces available seat width by 2"</i>	\$
Hip Belt Mount	<i>Provides improved mounting options for Bodypoint hip belt or other third-party belts.</i>	\$
Clear Lap Tray (Slide On)	<i>Available for seat widths 17"-22" only. Will interfere with drive wheels when in tilt.</i>	\$
Medical Necessities Backpack	<i>Black with Permobil logo.</i>	\$
Adhesive Back Hook Strip	<i>For cushion mounting, 2"x30".</i>	\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes