



PDG Fuze T50 - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

Wheelchair & Frame

Fuze T50 50° Manual Tilt Wheelchair Also available: the narrow Fuze T50n with an overall chair width of only seat width + 4". See separate order form.

\$

Frame Type

Fixed Width Frame, 300lbs (136kg) Weight Capacity Select a fixed seat width below.

No charge

Fixed Width Frame, 400lbs (181kg) Weight Capacity Not compatible with transit tie down system. Requires flat free rear tires or heavy duty spoke rear wheels. Select a fixed seat width below.

No charge

Adjustable Width Frame. 250lbs (113kg) Weight Capacity Select an adjustable seat width below.

\$

Frame Options

Take-Apart Frame The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. This feature reduces the weight of the heaviest piece of the wheelchair to 29lbs (13kg) or less, perfect for transporting the chair in the trunk of a car or SUV. Not compatible with power tilt, drum brakes, IV pole, O2 holder, ventilator tray, ventilator battery tray, dynamic rocker back or 13" seat to floor height.

\$

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Tested to maximum weight capacity of 300lbs (136kg).

No charge

Seat Measurements

Seat Width | Overall chair width is +7.5" with 12/16" rear wheels and +8.5" with 20/22/24" rear wheels.

Fixed:	15"	16"	17"	18"	19"	No charge
	20"	21"	22"	23"		\$
	24"					\$
Adjustable 12"-16":	12"	13"	14"	15"	16"	No charge
Adjustable 15"-20":	15"	16"	17"	18"	19"	No charge
	20"					\$

Seat Depth

Select Depth:	15"	16"	17"	18"	19"	No charge
	20"	21"				\$
	22"					\$

21"-22" seat depths are not adjustable and require straight stroller handle back. Not compatible with reclining or fold-down back.

Frame Finishes

Select Color:	High Gloss Black	Red Pepper	Neon Pink	Ocean Blue Metallic	Tangerine Metallic
	Acid Green Metallic	Midnight Blue Pearl	White Pearl	Black Lava	Titanium Grey
	Pomegranate Red	Blue Mist	Indigo Purple	Hyper Blue	Flame Orange
	OD Green	Silver Streak	Lavendar Dusk	Black Galaxy	Seafoam Green
	Daffodil Yellow				

Armrests

Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests \$

10-14" Range 12"-16" Range

Fixed Height Armrests 12" height, full length armpads only. No charge

Armrest Position: Mid (Default) Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Locking Flip-Up Armrests Not recommended for users with excessive strength or tone. Not compatible with fold-down back. Only with regular armpads. \$

Set at 10" Set at 11" Set at 12" Set at 13" Set at 14"

Armpads

Regular Armpads No charge

Desk Length 2" x 10" Full Length 2" x 14"

Waterfall Armpads \$

Desk Length 2" x 10" Full Length 2" x 14"

Gel Armpads May interfere with wheels in some tilt configurations. \$

Narrow 2" x 12" Wide 3.5" x 12"

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change.

Armrest Wideners, 1" Total Additional Width \$

Armrest Wideners, 2" Total Additional Width \$

Seat Height

Seat to floor height is adjustable within the range of the chosen axle plate (Low/Standard/High). Chairs with swing-away front rigging are set with 3° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out and standard plug-in front rigging are not affected.

		5" Casters					6" Casters & 6" Forks					6" Casters & 8" Forks					8" Casters & 8" Forks					6" Frog Legs Casters & Suspension Forks						
				Low	Low	Low		Low	Low	Low	Std		Low	Low	Low	Std	Std		Low	Low	Std	Std	Std		Low	Low	Std	
RW	24"			Low	Low	Low		Low	Low	Low	Std		Low	Low	Low	Std	Std		Low	Low	Std	Std			Low	Low	Std	
	22"		Low	Low	Low	Std		Low	Low	Low	Std	Std		Low	Low	Std	Std	Std		Low	Std	Std	Std		Low	Std	Std	
	20"	Low	Low	Low	Std	Std		Low	Low	Std	Std	Std		Low	Std	Std	Std		Std	Std	Std		High		Std	Std	Std	
	16"		Std	Std	Std			Std	Std	Std		High		Std	Std		High	High		Std		High	High	High		Std		High
	12"	Std	Std		High	High		Std		High	High	High			High	High	High		High	High	High	High	High		High	High	High	
		13"	14"	15"	16"	17"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	16"	17"	18"	19"	20"	16"	17"	18"				
Seat to Floor Height																												

Seat Pan

Seat Pan

No Seat Pan	No charge
Solid Aluminium Seat Pan	\$

Front Casters

Urethane Casters

5" x 1.25" Caster	No charge
6" x 1.25" Caster	No charge
8" x 1.25" Caster	No charge
6" x 1.50" Caster	\$
8" x 2.00" Caster	\$

Pneumatic Casters

6" x 1.25" Pneumatic Caster	\$
8" x 1.25" Pneumatic Caster	\$
8" x 2.00" Pneumatic Caster	\$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster	Urethane tire and silver anodized aluminum hub.	\$
----------------------------	---	----

Frog Legs Suspension Forks

FrogLegs Suspension Forks	Requires 6" soft-roll caster.	\$
---------------------------	-------------------------------	----

Rear Wheels

Composite Rear Wheel	No charge
12" 16" 20" 22" 24"	
Spoke Rear Wheel	\$
20" 22" 24"	
Heavy Duty Spoke Wheel	\$
24"	

Rear Wheel Position: Rear (More Stable) Mid (Default) Front (Less Stable)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires

Urethane Tires		No charge
Pneumatic Tires	<i>May result in a finished seat height 1/2" greater than stated.</i>	\$
Pneumatic Tires with Flat-Free Inserts	<i>May result in a finished seat height 1/2" greater than stated. Not available with spoke wheels.</i>	\$

Handrims | Handrims are only available with 20", 22" and 24" rear wheels.

No Handrims	No charge	
Standard Aluminium Handrims	No charge	
Plastic Coated Handrims	<i>Not compatible with drum brakes.</i>	\$
Projection Handrims	<i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$

Wheel Locks

User Wheel Locks

No charge

Push to Lock Pull to Lock

Attendant Wheel Locks *Not compatible with drum brakes or extensions.*

\$

Push to Lock Pull to Lock

User and Attendant Wheel Locks *Not compatible with drum brakes.*

\$

Push to Lock Pull to Lock

Foot Operated Wheel Lock *Rear-mounted unilateral single-activation lock. Not available with power tilt in combination with 12" or 16" rear wheels. Not compatible with O2 holder, IV pole or ventilator tray.*

\$

Drum Brakes *Available with 12" composite or 20" spoke rear wheels only. Urethane tires only. Not with plastic coated handrims. 12" wheels add 1" to overall chair width. 12" compatible with 13", 14", 17.5", 18.5", 19.5" STF heights only. 20" wheels compatible with 16", 17", 18" STF heights only.*

\$

User Activated Unilateral Wheel Lock *Pull to lock, Available up to 20" seat widths only.*

\$

Left Activation Right Activation

Wheel Lock Options

Wheel Lock Extensions *Not compatible with attendant wheel lock.*

\$

Anti-Tippers

Rear Anti-Tippers *Required safety option.*

\$

Backrest

Back Post Style: Stroller Handle Back Cane Handle Back

Cane handle back not compatible with power tilt with manual override.

Back Post Type: Straight 8° Bend

Back Post Height: Short (17") Tall (23")

Measured from seat pan to handle.

Folding Style: Fold-Down Non-Folding

Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.

Seat to Back Angle: Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle: *Not compatible with cane handle back.*

\$

Button on Outside Button on Inside

Endomorph Back *Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with flip-up armrests, dynamic rocker back or cane handle back.*

\$

1" Offset Per Side 1.5" Offset Per Side 2" Offset Per Side

Seating Dynamics Dynamic Rocker Back Interface *Not compatible with Take-Apart frame, fold-down back, endomorph back, 30° recline or cane handle back. Not available with 21" and 22" seat depths. Only available with stroller handle - 8° bend. Finished back post height is 7" taller than selected back post height. Maximum compatible seat depth is 20". 90° and 98° back angle only.*

\$

Recline

Reclining Backrest (Manual) *30° maximum recline. Not compatible with cane handle back or fold-down back. Not available with 21" and 22" seat depths. Finished back post height is 3" taller than selected back post height.*

\$

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post

No charge

Attendant Tilt Only, Left Back Post

No charge

Power Tilt | 40° tilt range. Includes 2 batteries for off-chair charging. Requires stroller handle back.

Standard Power Tilt, 250lb (113kg) Weight Capacity

\$

Standard Power Tilt, 350lb (158kg) Weight Capacity *Requires fixed frame.*

\$

Power Tilt with Manual Tilt Override, 250lb (113kg) Weight Capacity *Allows for instant manual release of tilt in case of power tilt failure.*

\$

Power Tilt with Manual Tilt Override, 350lb (158kg) Weight Capacity *Allows for instant manual release of tilt in case of power tilt failure. Requires fixed frame.*

\$

Toggle Switch Control		No charge
Installed on Right	Installed on Left	No charge
Button Switch Control		No charge
Installed in Right	Installed in Left	No charge
Stealth Egg Switch Control	<i>Not installed.</i>	No charge
Swing-Away Mounting Hardware	<i>Only for Stealth egg switch. Not compatible with all headrest makes and models.</i>	\$

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit:	Set at 20°	Set at 30°	Set at 40°
--------------------	------------	------------	------------

Front Rigging

No Front Rigging	No charge
-------------------------	-----------

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70°		\$									
Universal Adjustable Angle (0°-100°)	<i>Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.</i>	\$									
0°	10°	20°	30°	40°	50°	60°	70°	80°	90°	100°	
High Mount Option for Easy-Out Adjustable Angle	<i>For shorter lower leg lengths of 8"-12.5".</i>	\$									

Swing Away Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Swing Away 60°	No charge
Swing Away 70°	No charge
Swing Away 90° <i>Includes 11"-16" extension tubes. Not compatible with 8" casters.</i>	No charge
Swing Away Elevating Legrests (ELR) <i>Includes calf pads.</i>	\$
1" Offset Hangers for Swing Away Front Rigging	\$

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Plug In 70°	No charge
--------------------	-----------

Footplates

Composite Footplates	No charge
Alumnium Angle Adjustable Footplates	\$
One Piece Fixed Footplate Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates.	\$
One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates.	\$
Interlocking Angle Adjustable Footplates Available for seat widths up to 21"	\$

Extension Tubes | *Consider the height of your seat cushion.*

Standard Extension Tubes, Short	No charge
Standard Extension Tubes, Medium <i>Only with 90° swing-away front rigging.</i>	No charge
Standard Extension Tubes, Long	\$

Footplate Height:		Specify:				
Ext Tube	Short	8" - 13"	11" - 16"	12" - 16"	4" - 9"	13" - 16"
	Medium	N/A	N/A	N/A	10" - 12"	N/A
	Long	N/A	14" - 22"	14" - 19"	12" - 16"	14" - 17"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	90° Swing-Away	Standard Plug-In
Legrest						

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging	<i>Not compatible with gel padding. Includes reinforced premium calf pad mount.</i>	\$
Calf Pads, Swing Away Front Rigging	<i>Not compatible with gel padding.</i>	\$
Gel Padding	<i>Not compatible with calf pads or Easy-Out front rigging.</i>	\$
Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging		\$
Gel Knee Pads, Swing Away ELR Front Rigging		\$
Length Adjustable Heel Loops		\$
Standard Calf Strap		\$

Positioning Belts

Q'Straint Lap Belt	<i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	\$
---------------------------	---	----

Miscellaneous Options

One-Arm Drive	<i>Available with 22" and 24" composite rear wheels. Minimum 17" and 18" seat to floor height respectively. Large hand setup includes 24" rear wheels with 22" driving handrim. Not compatible with drum brakes, quick release axles or pneumatic tires.</i>	\$
Installed On:	Right Left	\$
Handrim Setup:	Standard Large Hand	
	<i>Large hand setup only available with 24" rear wheels.</i>	
Spoke Guards	<i>Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.</i>	\$
	Black Frame color	
Hip Guide, Bracket Only	<i>Installed with receivers. Accessory mount.</i>	\$
Gel Hip Guides	<i>Reduces available seat width by 2"</i>	\$
Hip Belt Mount	<i>Provides improved mounting options for Bodypoint hip belt or other third-party belts.</i>	\$
IV Pole		\$
	Installed on Right Installed on Left	\$
O2 Holder	<i>May add to overall width of chair.</i>	\$
	Installed on Right Installed on Left	\$
Ventilator Tray	<i>Reduced ground clearance. Chair may be unstable when unoccupied. Minimum seat width is 16". Not compatible with O2 holder, power tilt or foot operated wheel locks.</i>	\$
Clear Lap Tray (Slide On)	<i>Available for seat widths 17"-22" only. Not compatible with contour arm trough.</i>	\$
Medical Necessities Backpack	<i>Black with Permobil logo.</i>	\$
Adhesive Back Hook Strip	<i>For cushion mounting, 2"x30".</i>	\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes