



PDG Fuze T50n - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

Wheelchair & Frame

Fuze T50n 50° Manual Tilt Wheelchair Narrow frame.

\$

Frame Type

Fixed Width Frame, 300lbs (136kg) Weight Capacity Select a fixed seat width below.

No charge

Fixed Width Frame, 400lbs (181kg) Weight Capacity Not compatible with transit tie down system. Requires flat free rear tires. Select a fixed seat width below.

No charge

Adjustable Width Frame. 250lbs (113kg) Weight Capacity Select an adjustable seat width below.

\$

Frame Options

Take-Apart Frame The Take-Apart frame allows attendants to take the Fuze apart in seconds without tools. This feature reduces the weight of the heaviest piece of the wheelchair to 29lbs (13kg) or less, perfect for transporting the chair in the trunk of a car or SUV. Not compatible with power tilt, drum brakes, ventilator tray, ventilator battery tray, dynamic rocker back or user wheel locks.

\$

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Tested to maximum weight capacity of 300lbs (136kg).

No charge

Seat Measurements

Seat Width | Overall chair width is +4"

Fixed:	16"	17"	18"	19"	No charge
	20"	21"	22"	23"	\$
	24"				\$
Adjustable 16"-20":	16"	17"	18"	19"	No charge
	20"				\$

Seat Depth

Select Depth:	15"	16"	17"	18"	19"	No charge
	20"	21"				\$
	22"					\$

21"-22" seat depths are not adjustable and require straight stroller handle back. Not compatible with reclining or fold-down back.

Frame Finishes

Select Color:	High Gloss Black	Red Pepper	Neon Pink	Ocean Blue Metallic	Tangerine Metallic
	Acid Green Metallic	Midnight Blue Pearl	White Pearl	Black Lava	Titanium Grey
	Pomegranate Red	Blue Mist	Indigo Purple	Hyper Blue	Flame Orange
	OD Green	Silver Streak	Lavendar Dusk	Black Galaxy	Seafoam Green
	Daffodil Yellow				

Armrests

Armrest Type

No Armrests No charge

Single Post, Height Adjustable Armrests \$

10-14" Range 12"-16" Range

Fixed Height Armrests 12" height, full length armpads only. No charge

Armrest Position: Mid (Default) Front

If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.

Locking Flip-Up Armrests Not recommended for users with excessive strength or tone. Not compatible with fold-down back. Only with regular armpads. \$

Set at 10" Set at 11" Set at 12" Set at 13" Set at 14"

Armpads

Regular Armpads No charge

Desk Length 2" x 10" Full Length 2" x 14"

Waterfall Armpads \$

Desk Length 2" x 10" Full Length 2" x 14"

Gel Armpads May interfere with wheels in some tilt configurations. \$

Narrow 2" x 12" Wide 3.5" x 12"

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change.

Armrest Wideners, 1" Total Additional Width \$

Armrest Wideners, 2" Total Additional Width \$

Seat Height

Chairs with swing-away front rigging are set up with 3° of pre-tilt to avoid interference. This increases the actual seat to floor height at the front of the seat by up to 1" compared to the selected value. Easy-Out and standard plug-in front rigging are not affected.

Caster Size	12" Rear Wheels				
	8"				
	6"				
	6" Soft Roll				
	5"				
	16"	17"	18"	19"	20"
Seat to Floor Height					

Seat Pan

No Seat Pan No charge

Solid Aluminium Seat Pan \$

Front Casters

Urethane Casters

5" x 1.25" Caster No charge

6" x 1.25" Caster No charge

8" x 1.25" Caster No charge

6" x 1.50" Caster \$

8" x 2.00" Caster \$

Pneumatic Casters

6" x 1.25" Pneumatic Caster \$

8" x 1.25" Pneumatic Caster \$

8" x 2.00" Pneumatic Caster \$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster Urethane tire and silver anodized aluminum hub. \$

Rear Wheels

Composite Rear Wheel, 12"			No charge
Rear Wheel Position:	Rear (More Stable)	Mid (Default)	Front (Less Stable)
<i>If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.</i>			

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires

Urethane Tires	No charge
Pneumatic Tires <i>May result in a finished seat height 1/2" greater than stated.</i>	\$
Pneumatic Tires with Flat-Free Inserts <i>May result in a finished seat height 1/2" greater than stated.</i>	\$

Handrims

No Handrims	No charge
--------------------	-----------

Wheel Locks

User Wheel Locks		No charge
Push to Lock	Pull to Lock	
Attendant Wheel Locks		\$
Push to Lock	Pull to Lock	
Foot Operated Wheel Lock	<i>Rear-mounted unilateral single-activation lock. Not available with power tilt. Not compatible with ventilator tray or user wheel locks.</i>	\$
Drum Brakes	<i>Available with urethane tires only. Compatible with 17.5", 18.5", 19.5" STF heights only. Adds 1" to the overall chair width.</i>	\$
User Activated Unilateral Wheel Lock	<i>Pull to lock, Available up to 20" seat widths only.</i>	\$
Left Activation	Right Activation	

Wheel Lock Options

Wheel Lock Extensions	\$
------------------------------	----

Anti-Tippers

Rear Anti-Tippers	<i>Required safety option.</i>	\$
--------------------------	--------------------------------	----

Backrest

Back Post Style:	Stroller Handle Back	Cane Handle Back			
	<i>Cane handle back not compatible with power tilt with manual override.</i>				
Back Post Type:	Straight	8° Bend			
Back Post Height:	Short (17")	Tall (23")			
	<i>Measured from seat pan to handle.</i>				
Folding Style:	Fold-Down	Non-Folding			
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>				
Seat to Back Angle:	Set at 90°	Set at 98°	Set at 106°	Set at 114°	Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle:		<i>Not compatible with cane handle back.</i>	\$
Button on Outside	Button on Inside		
Endomorph Back		<i>Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with flip-up armrests, dynamic rocker back or cane handle back.</i>	\$
1" Offset Per Side	1.5" Offset Per Side	2" Offset Per Side	
Seating Dynamics Dynamic Rocker Back Interface		<i>Not compatible with Take-Apart frame, fold-down back, endomorph back, 30° recline or cane handle back. Not available with 21" and 22" seat depths. Only available with stroller handle - 8° bend. Finished back post height is 7" taller than selected back post height. Maximum compatible seat depth is 20". 90° and 98° back angle only.</i>	\$

Recline

Reclining Backrest (Manual) 30° maximum recline. Not compatible with cane handle back or fold-down back. Not available with 21" and 22" seat depths. \$
Finished back post height is 3" taller than selected back post height.

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post No charge
Attendant Tilt Only, Left Back Post No charge

Power Tilt | 40° tilt range. Includes 2 batteries for off-chair charging. Requires stroller handle back.

Standard Power Tilt, 250lb (113kg) Weight Capacity \$

Standard Power Tilt, 350lb (158kg) Weight Capacity Requires fixed frame. \$

Power Tilt with Manual Tilt Override, 250lb (113kg) Weight Capacity Allows for instant manual release of tilt in case of power tilt failure. Not compatible with cane handle back. \$

Power Tilt with Manual Tilt Override, 350lb (158kg) Weight Capacity Allows for instant manual release of tilt in case of power tilt failure. Requires fixed frame. Not compatible with cane handle back. \$

Toggle Switch Control No charge
Installed on Right Installed on Left

Button Switch Control No charge
Installed in Right Installed in Left

Stealth Egg Switch Control Not installed. No charge

Swing-Away Mounting Hardware Only for Stealth egg switch. Not compatible with all headrest makes and models. \$

Tilt Limit | Optional. If no tilt limit is selected the full tilt range will be available.

Tilt Limit: Set at 20° Set at 30° Set at 40°

Front Rigging

No Front Rigging No charge

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 70° \$

Universal Adjustable Angle (0°-100°) Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings. \$

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle For shorter lower leg lengths of 8"-12.5". \$

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

Swing Away 60° No charge

Swing Away 70° No charge

Swing Away 90° Includes 11"-16" extension tubes. Not compatible with 8" casters. No charge

Swing Away Elevating Legrests (ELR) Includes calf pads. \$

1" Offset Hangers for Swing Away Front Rigging \$

Plug In Front Rigging | Legacy design. Not compatible with Easy-Out.

Plug In 70° No charge

Footplates

Composite Footplates No charge

Aluminium Angle Adjustable Footplates \$

One Piece Fixed Footplate Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates. \$

One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates. \$

Interlocking Angle Adjustable Footplates Available for seat widths up to 21" \$

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short

No charge

Standard Extension Tubes, Medium Only with 90° swing-away front rigging.

No charge

Standard Extension Tubes, Long

\$

Footplate Height:

Specify:

Ext Tube	Short	8" - 13"	11" - 16"	12" - 16"	4" - 9"	13" - 16"
	Medium	N/A	N/A	N/A	10" - 12"	N/A
	Long	N/A	14" - 22"	14" - 19"	12" - 16"	14" - 17"

High Mount Easy-Out

Standard Easy-Out

60°, 70°, ELR Swing-Away

90° Swing-Away

Standard Plug-In

Legrest

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging Not compatible with gel padding. Includes reinforced premium calf pad mount.

\$

Calf Pads, Swing Away Front Rigging Not compatible with gel padding.

\$

Gel Padding Not compatible with calf pads or Easy-Out front rigging.

\$

Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging

\$

Gel Knee Pads, Swing Away ELR Front Rigging

\$

Length Adjustable Heel Loops

\$

Standard Calf Strap

\$

Positioning Belts

Q'Straint Lap Belt Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.

\$

Miscellaneous Options

Hip Guide, Bracket Only Installed with receivers. Accessory mount.

\$

Gel Hip Guides Reduces available seat width by 2"

\$

Hip Belt Mount Provides improved mounting options for Bodypoint hip belt or other third-party belts.

\$

Ventilator Tray Reduced ground clearance. Chair may be unstable when unoccupied. Minimum seat width is 20". Not compatible with power tilt or foot operated wheel locks.

\$

Clear Lap Tray (Slide On) Available for seat widths 17"-22" only. Not compatible with contour arm trough.

\$

Medical Necessities Backpack Black with Permobil logo.

\$

Adhesive Back Hook Strip For cushion mounting, 2"x30".

\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes