



PDG Stellar GL - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

Wheelchair & Frame

Stellar GL 45° Manual Tilt Wheelchair Also available: the lightweight Stellar GLT with aluminium frame. See separate order form.

\$

Frame Type

250lbs (113kg) Weight Capacity Available for 14" - 22" seat widths.

No charge

300lbs (136kg) Weight Capacity Available for 14" - 22" seat widths.

\$

450lbs (204kg) Weight Capacity Available for 20" - 32" seat widths.

\$

User Weight | Required for setup of tilt uplift force.

User Weight Range:

Up to 75lbs (34kg)

75lbs - 100lbs (34kg - 45kg)

100lbs - 115lbs (45kg - 52kg)

115lbs - 150lbs (52kg - 68kg)

150lbs - 200lbs (68kg - 90kg)

200lbs - 250lbs (90kg - 113kg)

250lbs - 300lbs (113kg - 136kg)

300lbs - 375lbs (136kg - 170kg)

375lbs - 450lbs (170kg - 204kg)

Frame Options

Extended Base with Rear Floor Jacks Lifts rear wheels off the ground for increased stability. Limited ground clearance. Not compatible with attendant wheel lock, HD anti-tipper, IV pole or O2 holder.

\$

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Not compatible with seat widths > 22". Tested to maximum weight capacity of 300lbs (136kg).

No charge

Seat Measurements

Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

Select Width:

14"

15"

16"

17"

18"

19"

No charge

20"

21"

22"

23"

\$

24"

25"

26"

27"

\$

28"

29"

30"

31"

32"

\$

Seat widths 23"-32" require 450lb (204kg) frame. Extended delivery lead times may apply.

Seat Depth

Select Depth:

16"

17"

18"

19"

No charge

20"

21"

\$

22"

\$

16" - 20" seat depths are adjustable. 21" and 22" seat depths require solid seat pan.

Frame Finishes

Select Color:

High Gloss Black

Red Pepper

Neon Pink

Ocean Blue Metallic

Tangerine Metallic

Acid Green Metallic

Midnight Blue Pearl

White Pearl

Black Lava

Titanium Grey

Pomegranate Red

Blue Mist

Indigo Purple

Hyper Blue

Flame Orange

OD Green

Silver Streak

Lavendar Dusk

Black Galaxy

Seafoam Green

Daffodil Yellow

Armrests

Armrest Type

No Armrests	No charge
Single Post, Height Adjustable Armrests	\$
8"-11" Range 10-14" Range 12"-16" Range	
Dual Post, Height Adjustable Armrests	\$
8"-11" Range 10-14" Range 12"-16" Range	
Fixed Height Armrests 12" height, full length armpads only.	No charge
Single Post Dual Post	
Armrest Position: Rear Front	
<i>If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.</i>	
Locking Flip-Up Armrests Not recommended for users with excessive strength or tone. Not compatible with fold-down back or clothing guards.	\$
Set at 10" Set at 11" Set at 12" Set at 13" Set at 14"	

Armpads

Regular Armpads	No charge
Desk Length 2" x 10" Full Length 2" x 14"	
Waterfall Armpads	\$
Desk Length 2" x 10" Full Length 2" x 14"	
Gel Armpads May interfere with wheels in some tilt configurations.	\$
Narrow 2" x 12" Wide 3.5" x 12"	

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change.

Armrest Wideners, 1" Total Additional Width	\$
Armrest Wideners, 2" Total Additional Width	\$

Clothing Guards | Not compatible with flip-up armrests.

Clothing Guards, Single Post Armrests Not compatible with 8"-11" height range.	\$
Clothing Guards, Dual Post Armrests Not compatible with 8"-11" height range.	\$

Armrest Extensions | Not compatible with flip-up armrests or desk length armpads.

Armrest Extensions - 3" Forward	\$
Armrest Extensions - 3" Rearward	\$

Seat Height

Seat Height - 6" Casters | Seat height is adjustable within the range of the chosen lower frame type (Low/High)

RW	Low						High					
	24"											
	22"											
	20"											
	16"											
	12"											
	13"	14"	15"	16"	17"	18"	15"	16"	17"	18"	19"	20"
Seat to Floor Height												

Seat Height - 8" Casters | Seat height is adjustable within the range of the chosen lower frame type (Low/High)

RW	Low				High			
	24"							
	16"							
	15"	16"	17"	18"	17"	18"	19"	20"
Seat to Floor Height								

Seat Pan

No Seat Pan	No charge
Solid Aluminium Seat Pan <i>Required for 21" and 22" seat depths.</i>	\$

Front Casters

Urethane Casters

6" x 1.25" Caster	No charge
8" x 1.25" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	No charge
6" x 1.50" Caster	\$
8" x 2.00" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Pneumatic Casters

6" x 1.25" Pneumatic Caster	\$
8" x 1.25" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$
8" x 2.00" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster <i>Urethane tire and silver anodized aluminum hub.</i>	\$
---	----

Frog Legs Suspension Forks

FrogLegs Suspension Forks <i>Only compatible with 16" and 24" rear wheels. Requires 6" soft-roll caster.</i>	\$
--	----

Rear Wheels

Composite Rear Wheel	No charge
12" 16" 20" 22" 24"	

Spoke Rear Wheel	\$
------------------	----

20" 22" 24"

Heavy Duty Spoke Wheel	\$
------------------------	----

24"

Rear Wheel Position: Rear Front (Default)

If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires

Urethane Tires	No charge
Pneumatic Tires <i>May result in a finished seat height 1/2" greater than stated.</i>	\$
Pneumatic Tires with Flat-Free Inserts <i>May result in a finished seat height 1/2" greater than stated. Not available with spoke wheels.</i>	\$

Handrims | Handrims are only available with 20", 22" and 24" rear wheels.

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$
Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Attendant Wheel Locks <i>Not compatible with wheel lock extensions.</i>	\$
Push to Lock Pull to Lock	
User and Attendant Wheel Locks <i>Not compatible with drum brakes on 20" rear wheels.</i>	\$
Push to Lock Pull to Lock	
Drum Brakes <i>Available with 12" composite or 20" spoke rear wheels only. Urethane tires only. Not compatible with plastic coated handrims or two trigger tilt setup.</i>	\$
User Activated Unilateral Wheel Lock <i>Pull to lock. Available up to 20" seat widths only. Not compatible with hill holders.</i>	\$
Left Activation Right Activation	

Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	\$
Hill Holders <i>Requires pneumatic rear tires with flat-free inserts.</i>	\$

Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	\$
Front Anti-Tippers	\$
Heavy Duty Rear Anti-Tippers <i>Provides additional rearward stability for users with excessive movement. Not compatible with attendant wheel locks, rear floor jacks, IV pole or O2 holder.</i>	\$

Backrest

Back Post Style:	Stroller Handle Back Cane Handle Back
	<i>Seat widths greater than 22" require stroller handle back. Cane handle back may reduce tilt range. Cane handle back not compatible with power tilt with manual override, dynamic rocker back or endomorph back.</i>
Back Post Height:	Short (22") Tall (25")
	<i>Measured from seat pan to handle.</i>
Folding Style:	Fold-Down Non-Folding
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>
Seat to Back Angle:	Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle:	
Button on Inside	
Button on Outside	
Endomorph Back <i>Backrest will be narrower than seat by total of 2x offset selected. Minimum backrest width is 12". Not compatible with flip-up armrests, dynamic rocker back, reclining backrest or cane handle back.</i>	\$
1" Offset Per Side 1.5" Offset Per Side 2" Offset Per Side	
Spring Back Suspension <i>Not compatible with reclining back or dynamic rocker back.</i>	No charge
Seating Dynamics Dynamic Rocker Back Interface <i>Not compatible with fold-down back, spring back suspension, endomorph back, 30° recline or cane handle back. Not available with 20", 21" and 22" seat depths. Finished back post height is 4" taller than selected back post height.</i>	\$

Recline

Reclining Backrest (Manual) <i>30° maximum recline. Not compatible with cane handle back, spring back suspension, fold-down back, dynamic rocker back or endomorph back. Finished back post height is 3" taller than selected back post height. Limits tilt by up to 20°.</i>	\$
--	----

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post	No charge
Attendant Tilt Only, Left Back Post	No charge

User Tilt Only | *Not available for 450lbs (204kg) weight capacity. Full length armpads only.*

User Tilt Only, Right Armrest	\$
User Tilt Only, Left Armrest	\$

Attendant & User Tilt | *Not available for 450lbs (204kg) weight capacity. Not compatible with tall (25") back posts or reclining backrest. Full length arm pads only.*

Attendant & User Tilt Operation \$

Attendant Trigger Position: Right Back Post Left Back Post

User Trigger Position: Right Armrest Left Armrest

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit: Set at 20° Set at 30° Set at 40°

Anterior Pre-Tilt | *Add up to 5° of Anterior Tilt to the tilt range.*

Anterior Pre-Tilt Set at 0° (Default) Set at 5°

Power Tilt | *250lbs (113kg) weight capacity. 45° tilt range. Includes 2 batteries for off-chair charging. Minimum seat width is 16". Power tilt is compatible with either O2 holder or IV pole, not both. Not available with 21" and 22" seat depths. Requires stroller handle.*

Standard Power Tilt \$

Power Tilt with Manual Tilt Override *Allows for instant manual release of tilt in case of power tilt failure.* \$

Toggle Switch Control No charge

Installed on Right Installed on Left

Button Switch Control No charge

Installed in Right Installed in Left

Stealth Egg Switch Control *Not installed.* No charge

Swing-Away Mounting Hardware *Only for Stealth egg switch. Not compatible with all headrest makes and models.* \$

Front Rigging

No Front Rigging No charge

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70° \$

Universal Adjustable Angle (0°-100°) *Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.* \$

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle *For shorter lower leg lengths of 8"-12.5".* \$

Swing Away Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Swing Away 60° No charge

Swing Away 70° No charge

Swing Away Elevating Legrests (ELR) *Includes calf pads.* \$

1" Offset Hangers for Swing Away Front Rigging \$

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Plug In 70° No charge

Footplates

Composite Footplates No charge

Alumnium Angle Adjustable Footplates \$

One Piece Fixed Footplate *Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates.* \$

One Piece Flip-Up Footplate *Requires selection of composite or aluminum angle adjustable footplates.* \$

Interlocking Angle Adjustable Footplates *Available for seat widths up to 21"* \$

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short	No charge
Standard Extension Tubes, Long	\$
Heavy Duty Extension Tubes, Short	\$
Heavy Duty Extension Tubes, Long	\$

Footplate Height: Specify:

Ext Tube	Std Short	8" - 13"	11" - 16"	12" - 16"	13" - 16"
	Std Long	N/A	14" - 22"	14" - 19"	14" - 17"
	HD Short	N/A	N/A	N/A	13" - 17"
	HD Long	N/A	N/A	N/A	15" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
Legrest					

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging	Not compatible with gel padding. Includes reinforced premium calf pad mount.	\$
Calf Pads, Swing Away Front Rigging	Not compatible with gel padding.	\$
Gel Padding	Not compatible with calf pads or Easy-Out front rigging.	\$
Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging		\$
Gel Knee Pads, Swing Away ELR Front Rigging		\$
Length Adjustable Heel Loops		\$
Standard Calf Strap		\$
Miller's Dynamic Footrest Mounts, Plug In Front Rigging	Not compatible with Easy-Out front rigging.	\$
Miller's Dynamic Footrest Mounts, Swing Away Front Rigging	Not compatible with Easy-Out front rigging. Not available on 13" STF height.	\$

Positioning Belts

Q'Straint Lap Belt	Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.	\$
--------------------	--	----

Miscellaneous Options

Spoke Guards	Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.	\$
Black	Frame color	
Gel Hip Guides	Reduces available seat width by 2"	\$
IV Pole	Not compatible with attendant wheel locks.	\$
Installed on Right	Installed on Left	
O2 Holder	Not compatible with attendant wheel locks.	\$
Installed on Right	Installed on Left	
Medical Necessities Backpack	Black with Permobil logo.	\$
Adhesive Back Hook Strip	For cushion mounting, 2"x30".	\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes