



PDG Stellar HD - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

**Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

Wheelchair & Frame

Stellar HD 45° Manual Tilt Wheelchair *Heavy gauge steel reinforcement of frame for bariatric users.*

\$

Frame Type

450lbs (204kg) Weight Capacity

No charge

600lbs (272kg) Weight Capacity *Tilt is limited to maximum 25°.*

\$

User Weight | *Required for setup of tilt uplift force.*

User Weight Range: Up to 250lbs (113kg)
375lbs - 450lbs (170kg - 204kg)

250lb - 300lbs (113kg - 136kg)
450lbs - 600lbs (204kg - 272kg)

300lbs - 375lbs (136kg - 170kg)

Seat Measurements

Seat Width | *Overall chair width is +11".*

Select Width:

22"

\$

24" 26"

\$

28" 30" 32"

\$

Seat Depth

Select Depth:

18" 19"

No charge

20" 21"

\$

22"

\$

Frame Finishes

Select Color:

High Gloss Black
Acid Green Metallic
Pomegranate Red
OD Green
Daffodil Yellow

Red Pepper
Midnight Blue Pearl
Blue Mist
Silver Streak

Neon Pink
White Pearl
Indigo Purple
Lavendar Dusk

Ocean Blue Metallic
Black Lava
Hyper Blue
Black Galaxy

Tangerine Metallic
Titanium Grey
Flame Orange
Seafoam Green

Armrests

Armrest Type

No Armrests	No charge
Single Post, Height Adjustable Armrests	\$
8"-11" Range 10-14" Range 12"-16" Range	
Dual Post, Height Adjustable Armrests	\$
8"-11" Range 10-14" Range 12"-16" Range	
Fixed Height Armrests 12" height, full length armpads only.	No charge
Single Post Dual Post	
Armrest Position: Rear (Default) Front	
If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.	

Armpads

Regular Armpads	No charge
Desk Length 2" x 10" Full Length 2" x 14"	
Waterfall Armpads	\$
Desk Length 2" x 10" Full Length 2" x 14"	
Gel Armpads May interfere with wheels in some tilt configurations.	\$
Narrow 2" x 12" Wide 3.5" x 12"	

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change.

Armrest Wideners, 1" Total Additional Width	\$
Armrest Wideners, 2" Total Additional Width	\$

Clothing Guards | Not compatible with flip-up armrests.

Clothing Guards, Single Post Armrests Not compatible with 8"-11" height range.	\$
Clothing Guards, Dual Post Armrests Not compatible with 8"-11" height range.	\$

Armrest Extensions | Not compatible with flip-up armrests or desk length armpads.

Armrest Extensions - 3" Forward	\$
Armrest Extensions - 3" Rearward	\$

Seat Height

Seat to Floor Height:	Set at 16"	Set at 17"	Set at 18"
-----------------------	------------	------------	------------

Seat Pan

Seat Pan

Solid Aluminium Seat Pan	\$
--------------------------	----

Front Casters

5" Casters with Heavy Duty Steel Fork	No charge
---------------------------------------	-----------

Rear Wheels

24" Heavy Duty Spoke Wheel	No charge
Rear Wheel Position: Rear Front (Default)	
If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.	

Axles

Fixed Rear Axles	No charge
------------------	-----------

Rear Tires

Urethane Tires	No charge
----------------	-----------

Handrims | Handrims are only available with 20", 22" and 24" rear wheels.

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims Not compatible with drum brakes.	\$

Wheel Locks

User Wheel Locks

No charge

Push to Lock Pull to Lock

Wheel Lock Options

Wheel Lock Extensions *Not compatible with attendant wheel lock.*

\$

Anti-Tippers

Rear Anti-Tippers *Required safety option.*

\$

Backrest

Stroller Handle Back

Back Post Height *Measured from seat pan to handle.*

Short (22") Tall (25")

Seat to Back Angle 90°

Back Options

Endomorph Back *Backrest will be narrower than seat by total of 2x offset selected. Not compatible with dynamic rocker back or reclining backrest.*

\$

1" Offset Per Side 1.5" Offset Per Side 2" Offset Per Side

Seating Dynamics Dynamic Rocker Back Interface *Not compatible with endomorph back or 30° recline. Not available with 21" and 22" seat depths. Finished back post height is 4" taller than selected back post height.*

\$

Recline

Reclining Backrest (Manual) *30° maximum recline. Not compatible with 600lb (272kg) weight capacity, dynamic rocker back or endomorph back. Finished back post height is 3" taller than selected back post height.*

\$

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only *Two tilt triggers are required for tilt activation.*

No charge

Tilt Limit | *Optional. If no tilt limit is selected the full tilt range will be available.*

Tilt Limit: Set at 20° Set at 30° Set at 40°

Maximum tilt limit for 600lb (272kg) weight capacity is 25°.

Front Rigging

No Front Rigging

No charge

Easy-Out Front Rigging | *Effortless insertion and removal. No hanger brackets remain on the frame.*

Fixed 70°

\$

Universal Adjustable Angle (0°-100°) *Select angle. Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.*

\$

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle *For shorter lower leg lengths of 8"-12.5".*

\$

Swing Away Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Swing Away 60°

No charge

Swing Away 70°

No charge

Swing Away Elevating Legrests (ELR) *Includes calf pads.*

\$

1" Offset Hangers for Swing Away Front Rigging

\$

Plug In Front Rigging | *Legacy design. Not compatible with Easy-Out.*

Plug In 70°

No charge

Footplates

Aluminium Angle Adjustable Footplates

\$

Extension Tubes | Consider the height of your seat cushion.

Standard Extension Tubes, Short

No charge

Standard Extension Tubes, Long

\$

Heavy Duty Extension Tubes, Short

\$

Heavy Duty Extension Tubes, Long

\$

Footplate Height:

Specify:

Ext Tube	Std Short	8" - 13"	11" - 16"	12" - 16"	13" - 16"
	Std Long	N/A	14" - 22"	14" - 19"	14" - 17"
	HD Short	N/A	N/A	N/A	13" - 17"
	HD Long	N/A	N/A	N/A	15" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In

Legrest

Front Rigging Accessories

Calf Pads, Swing Away Front Rigging *Not compatible with gel padding or Easy-Out front rigging.*

\$

Gel Padding *Not compatible with calf pads or Easy-Out front rigging.*

\$

Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging

\$

Gel Knee Pads, Swing Away ELR Front Rigging

\$

Length Adjustable Heel Loops

\$

Standard Calf Strap

\$

Miscellaneous Options

Gel Hip Guides *Reduces available seat width by 2"*

\$

IV Pole *Not compatible with attendant wheel locks.*

\$

O2 Holder *Not compatible with attendant wheel locks.*

\$

Medical Necessities Backpack *Black with Permobil logo.*

\$

Adhesive Back Hook Strip *For cushion mounting, 2"x30".*

\$

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes