



PDG Stellar Impact - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

***Required field**

Special considerations that need to be addressed here (e.g., diagnosis):

Wheelchair & Frame

Stellar Impact 45° Manual Tilt Wheelchair Ultra-durable steel reinforcement of frame for highly agitated users. Extended lead times may apply.

\$

Frame Type

250lbs (113kg) Weight Capacity

No charge

Frame Options

Extended Base with Rear Floor Jacks Lifts rear wheels off the ground for increased stability. Limited ground clearance. Not compatible with attendant wheel lock, HD anti-tipper, IV pole or O2 holder.

\$

Transit Tie Down System (TTDS) Includes 4 factory installed transport brackets. TTDS has been tested according to ANSI/RESNA WC19. See Owners Manual for complete details. Not compatible with seat widths > 22".

No charge

Seat Measurements

Seat Width | Overall chair width is +9.5" with 12"/16" rear wheels, +10.5" with 20"/22"/24" rear wheels.

Select Width:	14"	15"	16"	17"	18"	19"	No charge
	20"	21"	22"				\$

Seat Depth

Select Depth:	17"	18"	19"	No charge
	20"	21"		\$
	22"			\$

16"-20" seat depths are adjustable.

Frame Finishes

Select Color:	High Gloss Black	Red Pepper	Neon Pink	Ocean Blue Metallic	Tangerine Metallic
	Acid Green Metallic	Midnight Blue Pearl	White Pearl	Black Lava	Titanium Grey
	Pomegranate Red	Blue Mist	Indigo Purple	Hyper Blue	Flame Orange
	OD Green	Silver Streak	Lavendar Dusk	Black Galaxy	Seafoam Green
	Daffodil Yellow				

Armrests

Armrest Type

No Armrests			No charge
Single Post, Height Adjustable Armrests			\$
8"-11" Range	10-14" Range	12"-16" Range	
Dual Post, Height Adjustable Armrests			\$
8"-11" Range	10-14" Range	12"-16" Range	
Fixed Height Armrests	12" height, full length armpads only.		No charge
Single Post	Dual Post		
Armrest Position:	Rear (Default)	Front	
If left blank default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions PDG may not be able to accommodate stated preference.			

Armpads

Regular Armpads		No charge
Desk Length 2" x 10"	Full Length 2" x 14"	
Waterfall Armpads		\$
Desk Length 2" x 10"	Full Length 2" x 14"	
Gel Armpads	May interfere with wheels in some tilt configurations.	\$
Narrow 2" x 12"	Wide 3.5" x 12"	

Armrest Wideners | Moves armrests outwards by up to 1" each side. Seat size does not change.

Armrest Wideners, 1" Total Additional Width	\$
Armrest Wideners, 2" Total Additional Width	\$

Clothing Guards | Not compatible with flip-up armrests.

Clothing Guards, Single Post Armrests	Not compatible with 8"-11" height range.	\$
Clothing Guards, Dual Post Armrests	Not compatible with 8"-11" height range.	\$

Armrest Extensions | Not compatible with flip-up armrests or desk length armpads.

Armrest Extensions - 3" Forward	\$
Armrest Extensions - 3" Rearward	\$

Seat Height

Seat Height - 6" Casters | Seat height is adjustable within the range of the chosen lower frame type (Low/High)

RW	Low					High				
	24"									
	22"									
	20"									
	16"									
	12"									
	14"	15"	16"	17"	18"	17"	18"	19"	20"	21"
Seat to Floor Height										

Seat Height - 8" Casters | Seat height is adjustable within the range of the chosen lower frame type (Low/High)

RW	Low			High		
	24"					
	16"					
	16"	17"	18"	19"	20"	21"
Seat to Floor Height						

Seat Pan

No Seat Pan	No charge
Solid Aluminium Seat Pan	\$

Front Casters

Urethane Casters

6" x 1.25" Caster	No charge
8" x 1.25" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	No charge
6" x 1.50" Caster	\$
8" x 2.00" Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Pneumatic Casters | 8" casters are only compatible with 16" and 24" rear wheels.

6" x 1.25" Pneumatic Caster	\$
8" x 1.25" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$
8" x 2.00" Pneumatic Caster <i>Only compatible with 16" and 24" rear wheels.</i>	\$

Soft-Roll Casters

6" x 1.5" Soft-Roll Caster <i>Urethane tire and silver anodized aluminum hub.</i>	\$
---	----

Frog Legs Suspension Forks

FrogLegs Suspension Forks <i>Only compatible with 16" and 24" rear wheels. Requires 6" soft-roll caster.</i>	\$
--	----

Rear Wheels

Composite Rear Wheel	No charge
12" 16" 20" 22" 24"	
Spoke Rear Wheel	\$
20" 22" 24"	
Heavy Duty Spoke Wheel	\$
24"	
Rear Wheel Position: Rear Front (Default)	No charge
<i>If left blank, default setting will be used. In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.</i>	

Axles

Quick Release Rear Axles	No charge
Fixed Rear Axles	No charge

Rear Tires

Urethane Tires	No charge
Pneumatic Tires <i>May result in a finished seat height 1/2" greater than stated.</i>	\$
Pneumatic Tires with Flat-Free Inserts <i>May result in a finished seat height 1/2" greater than stated. Not available with spoke wheels.</i>	\$

Handrims | Handrims are only available with 20", 22" and 24" rear wheels.

No Handrims	No charge
Standard Aluminium Handrims	No charge
Plastic Coated Handrims <i>Not compatible with drum brakes.</i>	\$
Projection Handrims <i>Composite wheels only. 10 projections, 30°, 1-1.8" length.</i>	\$

Wheel Locks

User Wheel Locks	No charge
Push to Lock Pull to Lock	
Attendant Wheel Locks, Push to Lock <i>Not compatible with wheel lock extensions.</i>	\$
Push to Lock Pull to Lock	
User and Attendant Wheel Locks <i>Not compatible with drum brakes on 20" rear wheels.</i>	\$
Push to Lock Pull to Lock	
Drum Brakes <i>Available with 12" composite or 20" spoke rear wheels only. Urethane tires only.</i>	\$
User Activated Unilateral Wheel Lock <i>Pull to lock. Available up to 20" seat widths only. Not compatible with hill holders.</i>	\$
Left Activation Right Activation	

Wheel Lock Options

Wheel Lock Extensions <i>Not compatible with attendant wheel lock.</i>	\$
Hill Holders <i>Requires pneumatic rear tires with flat-free inserts.</i>	\$

Anti-Tippers

Rear Anti-Tippers <i>Required safety option.</i>	\$
Front Anti-Tippers	\$
Heavy Duty Rear Anti-Tippers <i>Provides additional rearward stability for users with excessive movement. Not compatible with attendant wheel locks, rear floor jacks, IV pole or O2 holder.</i>	\$

Backrest

Back Post Style:	Stroller Handle Back Cane Handle Back
	<i>Cane handle back may reduce tilt range. Cane handle back not compatible with dynamic rocker back or reclining back.</i>
Back Post Height:	Short (22") Tall (25")
	<i>Measured from seat pan to handle.</i>
Folding Style:	Fold-Down Non-Folding
	<i>Fold-down back not compatible with reclining backrest or dynamic rocker back and will be substituted with a non-folding back.</i>
Seat to Back Angle:	Set at 90° Set at 98° Set at 106° Set at 114° Set at 122°

Backrest Options

Angle/Height Adjustable Stroller Handle: <i>Not compatible with cane handle back.</i>	\$
Button on Outside Button on Inside	
Spring Back Suspension <i>Not compatible with reclining back.</i>	No charge
Seating Dynamics Dynamic Rocker Back Interface <i>Not compatible with fold-down back, spring back suspension, 30° recline or cane handle back. Not available with 21" and 22" seat depths. Finished back post height is 4" taller than selected back post height.</i>	\$

Recline

Reclining Backrest (Manual) <i>30° maximum recline. Not compatible with cane handle back, spring back suspension, fold-down back or dynamic rocker back. Finished back post height is 3" taller than selected back post height. Limits tilt by up to 20°.</i>	\$
--	----

Tilt Activation

Attendant Tilt Only

Attendant Tilt Only, Right Back Post	No charge
Attendant Tilt Only, Left Back Post	No charge

User Tilt Only | Full length armpads only.

User Tilt Only, Right Armrest	\$
User Tilt Only, Left Armrest	\$

Attendant & User Tilt | Not compatible with tall (25") back posts or reclining backrest. Full length armpads only.

Attendant & User Tilt Operation	\$
--	----

Attendant Trigger Position:	Right Back Post Left Back Post
User Trigger Position:	Right Armrest Left Armrest

Tilt Limit | Optional. If no tilt limit is selected the full tilt range will be available.

Tilt Limit:	Set at 20° Set at 30° Set at 40°
--------------------	----------------------------------

Front Rigging

No Front Rigging

No charge

Easy-Out Front Rigging | Effortless insertion and removal. No hanger brackets remain on the frame.

Fixed 70°

\$

Universal Adjustable Angle (0°-100°) Angle adjustment by loosening one bolt, in 10° increments. Tool-free leg length adjustment. Front rigging may interfere with castors in the 90°-100° settings.

\$

0° 10° 20° 30° 40° 50° 60° 70° 80° 90° 100°

High Mount Option for Easy-Out Adjustable Angle For shorter lower leg lengths of 8"-12.5".

\$

Swing Away Front Rigging | Legacy design. Not compatible with Easy-Out.

Swing Away 60°

No charge

Swing Away 70°

No charge

Swing Away Elevating Legrests (ELR) Includes calf pads.

\$

1" Offset Hangers for Swing Away Front Rigging

\$

Plug In Front Rigging | Legacy design. Not compatible with Easy-Out.

Plug In 70°

No charge

Footplates

Composite Footplates

No charge

Aluminum Angle Adjustable Footplates

\$

One Piece Fixed Footplate Available on plug in front rigging only. Requires selection of composite or aluminum angle adjustable footplates.

\$

One Piece Flip-Up Footplate Requires selection of composite or aluminum angle adjustable footplates.

\$

Interlocking Angle Adjustable Footplates Available for seat widths up to 21"

\$

Extension Tubes | Consider the height of your seat cushion.

Regular Extension Tubes, Short 8"-12.5" for Easy-Out high mount; 14"-18" for Easy-Out; 13"-17" for swing away and plug in.

No charge

Regular Extension Tubes, Long 11"-16" for Easy-Out high mount; 18"-24" for Easy-Out; 15"-19" for swing away and plug in.

\$

Heavy Duty Extension Tubes, Short 13"-17". Available on plug in front rigging only.

\$

Heavy Duty Extension Tubes, Long 15"-19". Available on plug in front rigging only.

\$

Footplate Height:

Specify:

Ext Tube	Std Short	8" - 13"	11" - 16"	12" - 16"	13" - 16"
	Std Long	N/A	14" - 22"	14" - 19"	14" - 17"
	HD Short	N/A	N/A	N/A	13" - 17"
	HD Long	N/A	N/A	N/A	15" - 19"
		High Mount Easy-Out	Standard Easy-Out	60°, 70°, ELR Swing-Away	Standard Plug-In
Legrest					

Front Rigging Accessories

Calf Pads, Easy-Out Front Rigging Not compatible with gel padding. Includes reinforced premium calf pad mount.

\$

Calf Pads, Swing Away Front Rigging Not compatible with gel padding.

\$

Gel Padding Not compatible with calf pads or Easy-Out front rigging.

\$

Gel Knee Pads, Easy-Out Angle Adjustable Front Rigging

\$

Gel Knee Pads, Swing Away ELR Front Rigging

\$

Length Adjustable Heel Loops

\$

Standard Calf Strap

\$

Miller's Dynamic Footrest Mounts, Plug In Front Rigging Not compatible with Easy-Out front rigging.

\$

Miller's Dynamic Footrest Mounts, Swing Away Front Rigging Not compatible with Easy-Out front rigging. Not available on 13" STF height.

\$

Positioning Belts

Q'Straint Lap Belt	<i>Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" belt length.</i>	\$
---------------------------	---	----

Miscellaneous Options

Strut Extension Damper	<i>Required seat damper for improved rocking dynamics.</i>	\$
-------------------------------	--	----

Spoke Guards	<i>Available for 20"-24" composite rear wheels. Made from durable and lightweight aluminium.</i>	\$
---------------------	--	----

Black	Frame color	\$
-------	-------------	----

Gel Hip Guides	<i>Reduces available seat width by 2"</i>	\$
-----------------------	---	----

IV Pole	<i>Not compatible with attendant wheel locks.</i>	\$
----------------	---	----

O2 Holder	<i>Not compatible with attendant wheel locks.</i>	\$
------------------	---	----

Medical Necessities Backpack	<i>Black with Permobil logo.</i>	\$
-------------------------------------	----------------------------------	----

Adhesive Back Hook Strip	<i>For cushion mounting, 2"x30".</i>	\$
---------------------------------	--------------------------------------	----

General

Permobil manufactures PDG wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes