

Panthera Bambino 3 - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

**Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

Chassis

Standard configurations include front frame angle 105°, seat angle 4°, standard colour option, fold-down backrest, carbon rear axle with 6° camber, front castors, flip-up footplate, push to lock brakes, standard rear wheels with aluminium handrims and anti-tips. Seat height with 24" wheels: front 41cm, rear 39cm. Seat height can be adjusted with different size rear wheels or rear axle raising kit. 60kg weight limit.

G517	Bambino 3	Seat width 24-33cm, seat depth 25-30cm.	\$
G518	Bambino 3, Short	Chassis is 5cm shorter. Seat width 24-27cm, seat depth 20-25cm	\$

Chassis Colour

Standard Colours:	Racing Blue	Anodic Black						Standard
Optional Colours:	Blue Classic	Mean Green	Panorange	Brave Bronze	White Pearl	Telemagenta		No Charge
	Grey Metallic	Ruby Red	Sting Yellow	Horizon Grey	Pine Green	Pearl Violet		

Non-standard colours require additional delivery time of 20 days. Additional information [here](#).

Seat Measurements

Seat Width:	24cm	27cm	30cm	33cm
	24cm Short	27cm Short		

Backrest

523BBHH-C	Height Adjustable Backrest, Blue Trim		Standard
	20-28cm 27-35cm		
523BBHHBL-C	Height Adjustable Backrest, Black Trim		\$ 0.00
	20-28cm 27-35cm		
523BBHH VER2-C	S3 Style Backrest Select height.		\$ 0.00
	20cm 25cm 30cm 35cm		
	40cm (for seat widths ≥27cm)		

Rear Wheels

Rear Axle

33030BB-C	Carbon Axle, 6°		Standard
33031BB-C	Carbon Axle, 6° + 10mm	Only for 120mm castors and 20" rear wheels.	\$ 0.00
33033BB-C	Carbon Axle, 8°		\$

Increased Seat Height

10mm	20mm	30mm	40mm	50mm
------	------	------	------	------

\$

Rear Wheels

Standard Wheels:	20"	22"	24"	Standard
Spinergy X Wheels:	24"			\$
	Black Spokes	White Spokes		
One Arm Drive:	22"	24"		\$
	Aluminium Pushrim			Included
	Friction Pushrim			\$

Not with 24cm seat width. Includes special rear axle. Additional delivery time of 5 days.

Tyres

Schwalbe Right Run 20"(451), 22"(501), 24"(540)	Standard
Rugged 24" (540)	\$
Schwalbe Marathon Plus 22"(501), 24"(540)	\$
Schwalbe One 24"(540)	\$
PU 20"(451), 22"(501), 24"(540) <i>Not compatible with reTyre skins.</i>	\$

Handrims

190DD30-C	Aluminium	Standard
190DD11-C	Titanium	\$
19051DDBT-C	Paragrip	\$
19052DDBT-C	Tetragrip	\$
19030DD-C	Maxgrip	\$

Front Wheels

Castor Selection:	90mm (3.5")	120mm (4.7")
--------------------------	-------------	--------------

90mm standard for 20" rear wheels.

Footrest

42341BB-C	Footrest, Flip-Up Adjustable	Standard
65040BB-C	Carbon Footplate	\$
65140BB-C	Angle-Adjustable Footplate	\$

Wheel Locks

4630200-C	Push to Lock	Standard
4630203-C	Push to Lock, Extended Lever	\$
4630234-C	Wheel Lock Blades <i>For solid PU tyres.</i>	\$

Side Guards / Armrests

	No Side Guards	Standard
4235000-C	Side Guards, Standard	\$
4235010-C	Side Guards, Short Frame	\$
4401400-C	Arm Supports with Side Guards, Standard	\$
4401410-C	Arm Supports with Side Guards, Short Frame	\$
4235020-C	Seat Width Adaptors, Pair <i>Narrows the seat width 7.5mm per side.</i>	\$

Push Handles

4233000-C	Height/Angle Adjustable Push Bar <i>Not for S3 backrest.</i>	\$
4252060 VER2-C	Height Adjustable Push Handles <i>Only for S3 backrest.</i>	\$
4252063 VER2-C	Height Adjustable Push Handles +10cm <i>Only for S3 backrest. Additional 10cm height.</i>	\$
42400BB VER2-C	Push Bar, S3 Style <i>Only for S3 backrest.</i>	\$

Accessories

Cushions

572BB34-C	Seat Cushion, Standard,2.5cm	\$
575BB34-C	Seat Cushion, Standard, 5cm	\$
572BB25-C	Seat Cushion, Short, 2.5cm <i>For short frame.</i>	\$
575BB25-C	Seat Cushion, Short, 5cm <i>For short frame.</i>	\$
57710BB-C	Wedge Cushion <i>Used to alter seat angle.</i>	\$

Safety

4353300-C	Anti-Tip, Pair	Standard
-----------	----------------	----------

Spoke Protectors

4500D11-C	Spoke Protectors, Clear <i>22" and 24" wheels only.</i>	\$
4500DDD-C	Spoke Protectors, with Logo	\$
45003DD-C	Spoke Protectors, Splash	\$
45001DD-C	Spoke Protectors, Car Rim	\$
45004DD-C	Spoke Protectors, Wild Animals	\$
45002DD-C	Spoke Protectors, Smiley	\$

Torso Supports

4236000-C	Torso Support, Left	\$
	Standard Narrow (9mm)	
4236000-C	Torso Support, Right	\$
	Standard Narrow (9mm)	
4380238-C	Torso Strap	\$
4710012-C	Backrest Wedge <i>To fit behind back upholstery for additional lateral support.</i>	\$

Head Support

4710100-C	Head Support, Blue Trim	\$
4710100BL-C	Head Support, Black Trim	\$

Harness Attachments

9900247-C	Attachment Plates, Two Point Seatbelt	\$
9900250-C	Attachment Plates, Four Point Seatbelt	\$
9900248-C	Clamps for Attachment Plates, Two Point	\$
9900248-C	Clamps for Attachment Plates, Four Point	\$

Tray Table

6810000-C	Tray Table	\$
6810010-C	Magnet Holder for Tray Table	\$

Other Accessories

9900251-C	Tube Protection, Front Frame	\$
9900251-C	Tube Protection, Rear Frame	\$
9900301-C	Crutch Holder	\$
9900261-C	Backpack	\$
1100300-C	reTyre Skins, 24" Pair	\$

FreeWheel

260010X-C	FreeWheel Attachment	\$
42342BB-C	FreeWheel Reinforcement Kit <i>Required for FreeWheel</i>	\$

Set-Up Adjustments

Backrest Angle:	Specify (mm) A = 14-28mm, 18mm standard.		
Rear Axle Adjustment:	Specify (mm) B = 77mm standard.		
Footrest Heights:	Specify (mm) A =	Specify (mm) B =	Specify (mm) K =
	'B' only required for angle adjustable footrests.		
Footrest Depth Position:	Standard	Forwards 30mm	
Footrest Height Position:	High	Low	

General

Permobil manufactures Panthera wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes


Reference Diagrams



Backrest Angle



Rear Wheel Position



Footrest Height