

Panthera Micro 3 - Order Form



Dealer Name*:		Account #*:
City, Post Code, Country:		
Phone:	Email:	P.O. Number
Clinician:	Client Reference*:	Client Weight:
Rep Name:		
*Required field Special considerations that need to be addressed here (e.g., diagnosis):		

Chassis

Standard configuration includes chassis, seat angle 0°, standard colour option, height adjustable backrest, 10° camber, 90mm front caster, footrest, anti-tip and rear wheels. Seat height with 18" wheels 30.5cm, 20" wheels 33cm.

G515	Micro 3	Seat width 24-27cm. Seat depth 15-20cm.	\$
G515	Micro 3 Long	Chassis is 5cm longer. Seat width 24-27cm. Seat depth 20-25cm.	\$

Chassis Colour

Standard Colour:	Racing Blue						Standard
Optional Colours:	Anodic Black	Blue Classic	Mean Green	Panorange	Brave Bronze	White Pearl	No Charge
	Telemagenta	Grey Metallic	Ruby Red	Sting Yellow	Horizon Grey	Pine Green	
	Pearl Violet						

Non-standard colours require additional delivery time of 20 days. Additional information [here](#).

Seat Measurements

Seat Width:	24cm	27cm
	24cm Long	27cm Long

Backrest Upholstery

536BB17-C	Black Upholstery, Blue Trim	Standard
536BB17BL-C	Black Upholstery, Black Trim	\$ 0.00

Rear Wheels

1000018-C	Micro 18" Rear Wheel	Seat height 30cm without cushion.	Standard
1000120-C	Micro 20" Rear Wheel	Seat height 32.5cm without cushion.	\$ 0.00

Footrest

43800BB-C	Footrest, Micro 3	Standard
-----------	--------------------------	----------

Side Guards

	No Side Guards	Standard
4400050-C	Side Guards, Standard	\$

Push Handles

Push Handles

4252167-C	Push Handle, Micro 3, Right	\$
4252168-C	Push Handle, Micro 3, Left	\$

Accessories

Seat Cushion

572BB20-C	Seat Cushion, Standard, 2.5cm	\$
575BB20-C	Seat Cushion, Standard, 5cm	\$
572BB25-C	Seat Cushion, Long, 2.5cm <i>For long frame.</i>	\$
575BB25-C	Seat Cushion, Long, 5cm <i>For long frame.</i>	\$

Safety

4640000-C	Anti-Roll Device, Micro 3	\$
-----------	---------------------------	----

Spoke Protectors

45000DD-C	Spoke Protectors, with Logo	\$
45005DD-C	Spoke Protectors, Ocean	\$
45006DD-C	Spoke Protectors, Cars	\$
45007DD-C	Spoke Protectors, Animals	\$
45008DD-C	Spoke Protectors, Splash	\$

Set-Up Adjustments

Footrest Height: Specify (mm) | A = Specify (mm) | K =

General

Permobil manufactures Panthera wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

Reference Diagrams

