

## Panthera S3 Swing - Order Form



Dealer Name\*:

Account #\*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference\*:

Client Weight:

Rep Name:

*\*Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

### Chassis

Standard configurations include front frame angle 105°, seat angle 4°, calf strap, standard colour option, backrest, carbon rear axle, front castors, swing-away footrests with foldable footplates, push to lock brakes and standard rear wheels with aluminium handrims. Seat height with 24" wheels: front 45cm, rear 43cm. Seat height can be adjusted with different size rear wheel or rear axle raising kit. 100kg weight limit for seat widths 33-42cm. 150kg weight limit for seat widths 45-50cm.

G547	<b>Panthera S3 Swing</b>	Seat width 36-50cm, seat depth 40cm.	\$
G553	<b>Panthera S3 Swing Short</b>	Chassis is 5cm shorter compared to S3 Swing. Seat width 33-39cm, seat depth 35cm.	\$
G5837	<b>Panthera S3 Swing Low 0°</b>	Chassis has flat 0° seat angle and lower seat height compared to S3 Swing. For seat width 36-45cm, seat depth 35, 37.5, 40. For seat width 50cm, seat depth 37.5, 40cm. 39cm front seat height with 24" rear wheels. <i>Extended model not kept in stock, longer leadtime.</i>	\$
G5818	<b>Panthera S3 Swing Long</b>	Chassis is 5cm longer compared to S3 Swing. Seat width 36-50cm, seat depth 37.5, 40, 42.5, 45cm. <i>Extended model not kept in stock, longer leadtime.</i>	\$

### Chassis Colour

<b>Standard Colours:</b>	Anodic Black							Standard
<b>Optional Colours:</b>	Blue Classic	Mean Green	Panorange	Brave Bronze	White Pearl	Telemagenta		No Charge
	Grey Metallic	Racing Blue	Ruby Red	Sting Yellow	Horizon Grey	Pine Green		
	Pearl Violet							

*Non-standard colours require additional delivery time of 20 days. Additional information [here](#).*

### Seat Measurements

**Seat Width:** 33cm 36cm 39cm 42cm 45cm 50cm

*Check allowed seat width based on chassis selection.*

**Seat Depth:** 35cm 37.5cm 40cm 42.5cm 45cm

*Check allowed seat depth based on chassis selection.*

### Backrest

<b>Standard Backrest</b>	20-45cm range. 20cm back height not available for 50cm seat width.	Standard
<b>Tapered Backrest 3cm</b>	30-40cm range.	\$
<b>Tapered Backrest 6cm</b>	30-40cm range. <i>Additional delivery time of 20 days.</i>	\$
<b>Flared Backrest 3cm</b>	25-40cm range. <i>Additional delivery time of 20 days.</i>	\$
<b>Curved Backrest</b>	25-40cm range.	\$
<b>Curved Tapered Backrest</b>	35-40cm range. <i>Additional delivery time of 20 days.</i>	\$
<b>Frame for Solid Backrests</b>	30cm. Backrest frame still folds.	\$
<b>Backrest Height:</b>	20cm 25cm 30cm 35cm 40cm 45cm	
	<i>Check allowed height based on backrest selection.</i>	

## Rear Wheels

### Rear Axle

33023BB-C	Carbon Standard	Standard
33023BB-C	Carbon +15mm <i>Use this rear axle when choosing arm supports and/or increased seat height.</i>	No charge
33035BB-C	Carbon 4° Camber <i>Only for seat widths 33-45cm.</i>	\$
33033BB-C	Carbon 8° Camber <i>Only for seat widths 33-39cm.</i>	\$
33025HH-C	Increased Seat Height	\$
	10mm    20mm    30mm    40mm    50mm	

### Rear Wheels

Standard Wheels:	22"	24"		Standard
Spinergy X Wheels:	24"	25"	26"	\$
	Black Spokes	White Spokes		
One Arm Drive:	22"	24"	26"	\$
	Aluminium Pushrim			Included
	Friction Pushrim			\$
<i>Includes special rear axle. Additional delivery time of 5 days.</i>				

### Tyres

Schwalbe Right Run 22"(501), 24"(540), 25" (559), 26"(590)	Standard
Rugged 24" (540)	\$
Schwalbe Marathon Plus 24"(540), 25" (559), 26"(590)	\$
Schwalbe One 24"(540), 25" (559)	\$
PU 24" <i>Not compatible with reTyre skins.</i>	\$

### Handrims

190DD30-C	Aluminium	Standard
190DD11-C	Titanium	\$
19051DDBT-C	Paragrip	\$
19052DDBT-C	Tetragrip	\$
19030DD-C	Maxgrip <i>22" and 24" only.</i>	\$

## Front Wheels

Castor Selection:	90mm (3.5")	120mm (4.7")	150mm (6")
-------------------	-------------	--------------	------------

### Footrest

670BB00-C	Foldable Footplates	Standard
672BB00-C	Foldable Footplates, Extended	\$
673BB00-C	Foldable Footplates, Large	\$
674BB00-C	Foldable Footplates, Large Extended	\$
	Heel Strap for Footplates	No charge

## Wheel Locks

4630000-C	Push to Lock	Standard
4630020-C	Push to Lock, Extended Lever	\$
4630150-C	Wheel Lock Blades <i>For solid PU tyres. Not with one hand brakes.</i>	\$
367BBBB-C	One Hand Brake, Right Activation	\$
368BBBB-C	One Hand Brake, Left Activation	\$

## Side Guards / Armrests

No Side Guards	Standard
Standard Side Guards	\$
Carbon Side Guards, Without Fender	\$
Carbon Side Guards, Small Fender	\$
Carbon Side Guards, Large Fender	\$
Side Guards With Arm Supports	\$

## Push Handles

No Push Handles	Standard
4252060 VER2-C Height Adjustable Push Handles	\$
4252063 VER2-C Height Adjustable Push Handles, +10cm Additional 10cm height.	\$
4245000-C Foldable Push Handles	\$
42400BB VER2-C Push Bar	\$

## Accessories

### Anti-Tips

4353000-C Anti-Tips, Standard Seat Heights	\$
4355000-C Anti-Tips, Short/Low For short and low chassis.	\$

### Cushions

572BB40-C Seat Cushion, 2.5cm	\$
575BB40-C Seat Cushion, 5cm	\$
57710BB-C Wedge Cushion Used to alter seat angle.	\$
990021X-C Seat Bag	\$
Mounted On Seat Upholstery Mounted Under Cushion	

### Spoke Protectors | Only for 22" and 24" rear wheels.

4500D11-C Spoke Protectors, Clear	\$
4500D00-C Spoke Protectors, with Logo	\$
45003DD-C Spoke Protectors, Splash	\$
45001DD-C Spoke Protectors, Car Rim	\$
45004DD-C Spoke Protectors, Wild Animals	\$
45002DD-C Spoke Protectors, Smiley	\$

### Torso Supports

6770305-C Torso Support Adjustable, Left	\$
6770306-C Torso Support Adjustable, Right	\$
4710012-C Backrest Wedge To fit behind back upholstery for additional lateral support.	\$

### Other Accessories

9900251-C Tube Protection, Front	\$
9900251-C Tube Protection, Rear	\$
6810020-C Tray Table Only together with arm supports.	\$
9900301-C Crutch Holder S3	\$
9900261-C Backpack	\$
reTyre Skins	\$
24" 25" 26"	

## Set-Up Adjustments

**Backrest Angle:** Specify (mm) | A =  
14-28mm, 20mm standard.

**Rear Axle Adjustment:** Specify (mm) | B =  
60mm standard.

**Footrest Heights:** Specify (mm) | A = Specify (mm) | B = Specify (mm) | K =  
'B' only required for angle adjustable footrests.

## General

Permobil manufactures Panthera wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes


## Reference Diagrams



Backrest Angle



Rear Axle Adjustment



Footrest Height