

## Panthera S3 - Order Form



Dealer Name\*:

Account #\*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference\*:

Client Weight:

Rep Name:

*\*Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

### Chassis

Standard configurations include front frame angle 105°, seat angle 7°, calf strap, standard colour option, backrest, carbon rear axle, front castors, rigid footrest, push to lock brakes and standard rear wheels with aluminium handrims. Seat height with 24" wheels: front 47 cm, rear 43 cm. Seat height can be adjusted with different size rear wheels or rear axle raising kit.

100kg weight limit for seat widths 27-42cm. 150kg weight limit for seat widths 45-50cm and S3 Large chassis.

G548	<b>Panthera S3</b>	Seat width 33-50cm, seat depth 40cm.	\$
G549	<b>Panthera S3 Short</b>	Shorter chassis compared to S3. For seat width 30-33cm, seat depth 27-33cm. For seat width 36-45cm, seat depth 35cm.	\$
G552	<b>Panthera S3 Short Low</b>	Chassis is 5cm shorter and 2.5cm lower compared to S3. Seat width 33-39cm, seat depth 35cm.	\$
G554	<b>Panthera S3 Large</b>	Chassis is 5cm longer and 2.5cm higher compared to S3. Seat width 39-50cm, seat depth 45cm.	\$
G5831	<b>Panthera S3 Long</b>	Chassis is 5cm longer compared to S3. Seat width 33-50cm, seat depth 35, 37.5, 40, 42.5cm. <i>Extended model not kept in stock, longer leadtime.</i>	\$
G5831	<b>Panthera S3 0°</b>	Chassis has flat 0° seat angle, compared to 7° seat angle of S3. Seat width 33-50cm, seat depth 35, 37.5, 40, 42.5cm. 45cm front seat height with 24" rear wheels. <i>Extended model not kept in stock, longer leadtime.</i>	\$
G555	<b>Panthera S3 Short Abducted</b>	Chassis is 5cm shorter compared to S3, with abducted front frame 6cm wider than seat width. Seat width 27-33cm, seat depth 27-33cm.	\$

### Chassis Colour

<b>Standard Colours:</b>	Anodic Black							Standard
<b>Optional Colours:</b>	Blue Classic	Mean Green	Panorange	Brave Bronze	White Pearl	Telemagenta		No Charge
	Grey Metallic	Racing Blue	Ruby Red	Sting Yellow	Horizon Grey	Pine Green		
	Pearl Violet							

*Non-standard colours require additional delivery time of 20 days. Additional information [here](#).*

### Seat Measurements

<b>Seat Width:</b>	27cm	30cm	33cm	36cm	39cm	42cm	45cm	50cm
--------------------	------	------	------	------	------	------	------	------

*Check allowed seat width based on chassis selection.*

<b>Seat Depth:</b>	35cm	37.5cm	40cm	42.5cm	45cm
--------------------	------	--------	------	--------	------

*Check allowed seat depth based on chassis selection.*

### Backrest

<b>Standard Backrest</b>	20-45cm range. 20cm back height not available for 50cm seat width.	Standard
<b>Tapered Backrest 3cm</b>	30-40cm range.	\$
<b>Tapered Backrest 6cm</b>	30-40cm range. <i>Additional delivery time of 20 days.</i>	\$
<b>Flared Backrest 3cm</b>	25-40cm range. <i>Additional delivery time of 20 days.</i>	\$
<b>Curved Backrest</b>	25-40cm range.	\$
<b>Curved Tapered Backrest</b>	35-40cm range. <i>Additional delivery time of 20 days.</i>	\$
<b>Frame for Solid Backrests</b>	30cm. Backrest frame still folds.	\$

Backrest Height: 20cm 25cm 30cm 35cm 40cm 45cm

Check allowed height based on backrest selection.

## Rear Wheels

### Rear Axle

33023BB-C	Carbon Standard	Standard
33023BB-C	Carbon +15mm Use this rear axle when choosing arm supports and/or increased seat height.	\$ 0.00
33035BB-C	Carbon 4° Camber Only for seat widths 27-45cm.	\$
33029BB-C	Carbon 6° Camber Only for seat widths 27-33cm.	\$ 0.00
33033BB-C	Carbon 8° Camber Only for seat widths 27-39cm.	\$
33025HH-C	Increased Seat Height	\$
	10mm 20mm 30mm 40mm 50mm	

### Rear Wheels

Standard Wheels:	22" 24"	Standard
Spinergy X Wheels:	24" 25" 26"	\$
	Black Spokes White Spokes	
One Arm Drive:	22" 24" 26"	\$
	Aluminium Pushrim	Included
	Friction Pushrim	\$

Includes special rear axle. Additional delivery time of 5 days.

### Tyres

Schwalbe Right Run 22"(501), 24"(540), 25" (559), 26"(590)	Standard
Rugged 24" (540)	\$
Schwalbe Marathon Plus 24"(540), 25" (559), 26"(590)	\$
Schwalbe One 24"(540), 25" (559)	\$
PU 24" Not compatible with reTyre skins.	\$

### Handrims

190DD30-C	Aluminium	Standard
190DD11-C	Titanium	\$
19051DDBT-C	Paragrip	\$
19052DDBT-C	Tetragrip	\$
19030DD-C	Maxgrip 22" and 24" only.	\$

## Front Wheels

Castor Selection: 90mm (3.5") 120mm (4.7") 150mm (6")

## Footrest

36001BB-C	V-Shape Footplate	Standard
36111BB-C	V-Shape Footplate, Extended Reinforced Standard for S3 Large.	No charge
3600BBB-C	U-Shape Footplate, Small Heel Bar	No charge
65040BB-C	Carbon Footplate	\$
65140BB-C	Angle-Adjustable Footplate	\$
670BB00-C	Foldable Footplates	No charge
672BB00-C	Foldable Footplates, Extended Longer insert tube.	\$
673BB00-C	Foldable Footplates, Large Bigger footplate, extended 8cm forward.	\$
674BB00-C	Foldable Footplates, Large Extended	\$
36010BB-C	Plastic Plate for Rigid Footrests	\$
	Heel Strap for Foldable Footplates	No charge

## Wheel Locks

4630000-C	Push to Lock	Standard
4630020-C	Push to Lock, Extended Lever	\$
4630150-C	Wheel Lock Blades <i>For solid PU tyres. Not with one hand brakes.</i>	\$
367BBB-C	One Hand Brake, Right Activation	\$
368BBB-C	One Hand Brake, Left Activation	\$

## Side Guards / Armrests

No Side Guards	Standard
Standard Side Guards	\$
Carbon Side Guards, Without Fender	\$
Carbon Side Guards, Small Fender	\$
Carbon Side Guards, Large Fender	\$
Side Guards With Arm Supports	\$

## Push Handles

	No Push Handles		Standard
4252060 VER2-C	Height Adjustable Push Handles		\$
4252063 VER2-C	Height Adjustable Push Handles, +10cm	Additional 10cm height.	\$
4245000-C	Foldable Push Handles		\$
42400BB VER2-C	Push Bar		\$

## Accessories

### Safety

4353000-C	Anti-Tips, Standard Seat Heights	\$
4355000-C	Anti-Tips, Short/Low <i>For short and low chassis.</i>	\$
4353080-C	Anti-Tips, Large <i>For large chassis.</i>	\$

### Cushions

572BB40-C	Seat Cushion, 2.5cm	\$
575BB40-C	Seat Cushion, 5cm	\$
57710BB-C	Wedge Cushion <i>Used to alter seat angle.</i>	\$
990021X-C	Seat Bag	\$
	Mounted On Seat Upholstery      Mounted Under Cushion	

### Spoke Protectors | *Only for 22" and 24" rear wheels.*

4500D11-C	Spoke Protectors, Clear	\$
4500D00-C	Spoke Protectors, with Logo	\$
45003DD-C	Spoke Protectors, Splash	\$
45001DD-C	Spoke Protectors, Car Rim	\$
45004DD-C	Spoke Protectors, Wild Animals	\$
45002DD-C	Spoke Protectors, Smiley	\$

### Torso Supports

6770305-C	Torso Support Adjustable, Left	\$
6770306-C	Torso Support Adjustable, Right	\$
4710012-C	Backrest Wedge <i>To fit behind back upholstery for additional lateral support.</i>	\$

Other Accessories

9900251-C	Tube Protection, Front	\$
9900251-C	Tube Protection, Rear	\$
6810020-C	Tray Table <i>Only together with arm supports.</i>	\$
9900301-C	Crutch Holder	\$
9900261-C	Backpack	\$
reTyre Skins		\$
24"    25"    26"		
260010X-C	FreeWheel Attachment <i>Requires extended reinforced or angle-adjustable footplate.</i>	\$

Set-Up Adjustments

Backrest Angle:	Specify (mm)   A =		
	14-28mm, 20mm standard.		
Rear Axle Adjustment:	Specify (mm)   B =		
	60mm standard.		
Footrest Heights:	Specify (mm)   A =	Specify (mm)   B =	Specify (mm)   K =
	<i>'B' only required for angle adjustable footrests.</i>		

General

Permobil manufactures Panthera wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes


Reference Diagrams



Backrest Angle



Rear Axle Adjustment



Footrest Height