

Panthera U3 - Order Form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

**Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

Chassis

Standard configurations include tapered frame with front frame angle 105°, seat angle 7°, calf strap, standard colour option, backrest, carbon rear axle, front castors, rigid footplate, push to lock brakes and standard rear wheels with aluminium handrims. Seat height with 24" wheels: front 47cm, rear 43cm. Seat height can be adjusted with different size rear wheels or rear axle raising kit. 100kg weight limit for seat widths 33-42cm. 150kg weight limit for 45cm seat width.

G551	U3	Seat width 33-45cm, seat depth 40cm.	\$
G5801	U3 Y Front	Chassis gives full frame width for user's thighs compared to U3. Seat width 33-50cm, seat depth 35, 37.5, 40, 42.5cm. <i>Extended model not kept in stock, longer leadtime.</i>	\$
G5802	U3 Y Front Short	Chassis is 5cm shorter and gives full frame width for user's thighs compared to U3. For seat width 33-45cm, seat depth 35, 37.5cm. For seat width 50cm, seat depth 37.5cm. <i>Extended model not kept in stock, longer leadtime.</i>	\$
G584	U3 Y Front Low	Chassis is 2.5 lower and gives full frame width for user's thighs compared to U3. Seat width 33-50cm, seat depth 35, 37.5, 40, 42.5cm. <i>Extended model not kept in stock, longer leadtime.</i>	\$

Chassis Colour

Standard Colours:	Anodic Black							Standard
Optional Colours:	Blue Classic	Mean Green	Panorange	Brave Bronze	White Pearl	Telemagenta		No Charge
	Grey Metallic	Racing Blue	Ruby Red	Sting Yellow	Horizon Grey	Pine Green		
	Pearl Violet							

Non-standard colours require additional delivery time of 20 days. Additional information [here](#).

Seat Measurements

Seat Width: 33cm 36cm 39cm 42cm 45cm

Check allowed seat width based on chassis selection.

Seat Depth: 35cm 37.5cm 40cm 42.5cm

Check allowed seat depth based on chassis selection.

Backrest

Standard Backrest	20-45cm range. 20cm back height not available for 50cm seat width.	\$ 0.00
Tapered Backrest 3cm	30-40cm range.	\$
Tapered Backrest 6cm	30-40cm range. <i>Additional delivery time of 20 days.</i>	\$
Flared Backrest 3cm	25-40cm range. <i>Additional delivery time of 20 days.</i>	\$
Curved Backrest	25-40cm range.	\$
Curved Tapered Backrest	35-40cm range. <i>Additional delivery time of 20 days.</i>	\$
Frame for Solid Backrests	30cm. Backrest frame still folds.	\$

Backrest Height: 20cm 25cm 30cm 35cm 40cm 45cm

Check allowed height based on backrest selection.

Rear Wheels

Rear Axle

33023BB-C	Carbon Standard	Standard
33023BB-C	Carbon +15mm <i>Use this rear axle when choosing arm supports and/or increased seat height.</i>	\$ 0.00
33035BB-C	Carbon 4° Camber	\$
	Carbon 8° Camber <i>Only for seat widths 33-39cm.</i>	\$
Increased Seat Height		\$
10mm	20mm	30mm
40mm	50mm	

Rear Wheels

Standard Wheels:	22"	24"	Standard
Spinergy X Wheels:	24"	25"	26"
	Black Spokes	White Spokes	
One Arm Drive:	22"	24"	26"
	Aluminium Pushrim		Included
	Friction Pushrim		\$

Includes special rear axle. Additional delivery time of 5 days.

Tyres

Schwalbe Right Run 22"(501), 24"(540), 25" (559), 26"(590)	No charge
Rugged 24" (540)	\$
Schwalbe Marathon Plus 24"(540), 25" (559), 26"(590)	\$
Schwalbe One 24"(540), 25" (559)	\$
PU 24" <i>Not compatible with reTyre skins.</i>	\$

Handrims

190DD30-C	Aluminium	Standard
190DD11-C	Titanium	\$
19051DDBT-C	Paragrip	\$
19052DDBT-C	Tetragrip	\$
19030DD-C	Maxgrip 22" and 24" only.	\$

Front Wheels

Castor Selection:	90mm (3.5")	120mm (4.7")
-------------------	-------------	--------------

Footrest

36915BB-C	U-Shape Footplate	Standard
36006BB-C	U-Shape Footplate, Extended	\$ 0.00
36916BB-C	U-Shape Footplate, Small Heel Bar	\$ 0.00
65050BB-C	Carbon Footplate	\$
65150BB-C	Angle-Adjustable Footplate	\$
36010BB-C	Plastic Plate for Rigid Footrest	\$

Wheel Locks

4630000-C	Push to Lock	Standard
4630020-C	Push to Lock, Extended Lever	\$
4630150-C	Wheel Lock Blades <i>For solid PU tyres. Not with one hand brakes.</i>	\$
367BBBB-C	One Hand Brake, Right Activation	\$
368BBBB-C	One Hand Brake, Left Activation	\$

Side Guards / Armrests

No Side Guards	Standard
Standard Side Guards	\$
Carbon Side Guards, Without Fender	\$
Carbon Side Guards, Small Fender	\$
Carbon Side Guards, Large Fender	\$
Side Guards With Arm Supports	\$

Push Handles

No Push Handles	Standard
4252060 VER2-C Height Adjustable Push Handles	\$
4252063 VER2-C Height Adjustable Push Handles, +10cm Additional 10cm height.	\$
4245000-C Foldable Push Handles	\$
42400BB VER2-C Push Bar	\$

Accessories

Safety

4353000-C Anti-Tips, Standard Seat Heights	\$
4355000-C Anti-Tips, Short/Low For short and low chassis.	\$

Cushions

572BB40-C Seat Cushion, 2.5cm	\$
575BB40-C Seat Cushion, 5cm	\$
57710BB-C Wedge Cushion Used to alter seat angle.	\$
990021X-C Seat Bag	\$
Mounted On Seat Upholstery Mounted Under Cushion	

Spoke Protectors | Only for 22" and 24" rear wheels.

4500D11-C Spoke Protectors, Clear	\$
4500D00-C Spoke Protectors, with Logo	\$
45003DD-C Spoke Protectors, Splash	\$
45001DD-C Spoke Protectors, Car Rim	\$
45004DD-C Spoke Protectors, Wild Animals	\$
45002DD-C Spoke Protectors, Smiley	\$

Torso Supports

6770305-C Torso Support Adjustable, Left	\$
6770306-C Torso Support Adjustable, Right	\$
4710012-C Backrest Wedge To fit behind back upholstery for additional lateral support.	\$

Other Accessories

9900251-C Tube Protection, Front	\$
9900251-C Tube Protection, Rear	\$
6810020-C Tray Table Only together with arm supports.	\$
9900301-C Crutch Holder	\$
9900261-C Backpack	\$
reTyre Skins	\$
24" 25" 26"	
260010X-C FreeWheel Attachment Required angle-adjustable footplate.	\$

Set-Up Adjustments

Backrest Angle: Specify (mm) | A =
14-28mm, 20mm standard.

Rear Axle Adjustment: Specify (mm) | B =
60mm standard.

Footrest Heights: Specify (mm) | A = Specify (mm) | B = Specify (mm) | K =
'B' only required for angle adjustable footrests.

General

Permobil manufactures Panthera wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes


Reference Diagrams



Backrest Angle



Rear Axle Adjustment



Footrest Height