

## Panthera X3 - Order Form



Dealer Name\*:

Account #\*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference\*:

Client Weight:

Rep Name:

*\*Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

### Chassis

Standard configuration includes tapered frame with 100° front angle, 35-46cm depth-adjustable upholstery with 7° seat angle, backrest, front castors, rigid footplate, push to lock brakes, Spinergy rear wheels with 2° camber and aluminium handrims. Seat height A1 with 24" rear wheels: front 47cm, rear 43cm. Seat height can be adjusted with different size rear wheels. 125kg weight limit.

G355

Panthera X3

\$

**Balance Position** | Available in two different balance and height positions. The position cannot be changed later. Seat height 'B' is 15mm lower.

**Standard Seat Height:** A1 | Most Active

A2 | Less Active

**Low Seat Height:** B1 | Most Active

B2 | Less Active

### Seat Measurements

**Seat Width:** 33cm 36cm 39cm 42cm 45cm

### Backrest

**Standard Backrest, Carbon Fibre** Ovalised back tubes.

Standard

**Backrest Frame for Solid Backrests** U3L style, round back tubes. 30cm height. *Custom request, additional leadtime of 5 days.*

No charge

**Backrest Height:** 22cm 23cm 24cm 25cm 26cm 27cm 28cm 29cm 30cm 31cm 32cm 33cm 34cm 35cm

*Only for standard carbon fibre backrest.*

**Backrest Position:** P1 | Most Active

P2 | Active

P3 | Less Active

P4 | Least Active

*Backrest position affects the balance. The position can be easily adjusted later.*

### Rear Wheels

**Spinergy X Wheels:** 24" 25" 26"

Standard

Black Spokes White Spokes

1100023-C

**Quick Release Axles, Steel** Required for user weight >100kg.

No charge

### Tyres

Schwalbe Right Run 24"(540), 25" (559), 26"(590)

Standard

Schwalbe Marathon Plus 24"(540), 25" (559), 26"(590)

\$

Schwalbe One 24"(540), 25" (559)

\$

Rugged 24" (540)

\$

PU Solid 24" (540)

\$

## Handrims

|             |                                          |          |
|-------------|------------------------------------------|----------|
| 190DD30-C   | Aluminium                                | Standard |
| 190DD11-C   | Titanium                                 | \$       |
| 19044DD-C   | Aluminium X                              | \$       |
| 19051DDBT-C | Paragrip                                 | \$       |
| 19052DDBT-C | Tetragrip                                | \$       |
| 19030DD-C   | Maxgrip 24" wheels only. Contains latex. | \$       |

## Front Castors

|           |                                                                                        |           |
|-----------|----------------------------------------------------------------------------------------|-----------|
| 2100016-C | X Style 87mm (3")                                                                      | No charge |
| 2100204-C | S3 Style, 90mm (3.5")                                                                  | No charge |
| 2100205-C | S3 Style, 120mm (4.7") Not compatible with 24" rear wheels and B1/B2 balance position. | No charge |

## Footrest

|           |                                                      |           |
|-----------|------------------------------------------------------|-----------|
| 36915BB-C | U-Shape                                              | Standard  |
| 36006BB-C | U-Shape, Extended Reinforced                         | No charge |
| 36916BB-C | U-Shape, Small Heel Bar                              | No charge |
| 36004BB-C | Aluminium U-Shape                                    | \$        |
| 65060BB-C | Carbon Footplate                                     | \$        |
| 36010BB-C | Plastic Plate for Footrest Not with Carbon Footrest. | \$        |

## Wheel Locks

|           |                                  |           |
|-----------|----------------------------------|-----------|
| 4630000-C | Push to Lock S3                  | Standard  |
| 462BB02-C | One Hand Brake, Right Activation | No charge |
| 462BB01-C | One Hand Brake, Left Activation  | No charge |

## Side Guards

|                                    |    |
|------------------------------------|----|
| Standard Side Guards, S3 Style     | \$ |
| Carbon Side Guards, Without Fender | \$ |
| Carbon Side Guards, Small Fender   | \$ |
| Carbon Side Guards, Large Fender   | \$ |

## Push Handles

|           |                                                        |    |
|-----------|--------------------------------------------------------|----|
| 4200016-C | Arm Hooks, X3                                          | \$ |
| 4200017-C | Foldable Push Handles, X3 For standard carbon backrest | \$ |
| 4245000-C | Foldable Push Handles For U3L style backrest           | \$ |

## Accessories

### Seat Cushion

|           |                                                       |    |
|-----------|-------------------------------------------------------|----|
| 575BB40-C | Seat Cushion 5cm                                      | \$ |
| 575BB35-C | Seat Cushion 5cm, Short                               | \$ |
| 572BB40-C | Seat Cushion 2.5cm                                    | \$ |
| 575BB35-C | Seat Cushion 2.5cm, Short                             | \$ |
| 57710BB-C | Wedge Cushion Used to alter seat angle.               | \$ |
| 990021X-C | Seat Bag                                              | \$ |
|           | Mounted On Seat Upholstery      Mounted Under Cushion |    |

### Safety

|           |                 |    |
|-----------|-----------------|----|
| 4354104-C | Anti-Tips, Pair | \$ |
|-----------|-----------------|----|

## Other

|           |                                                                    |    |
|-----------|--------------------------------------------------------------------|----|
| 9024009-C | Transport Protection, Castor Housing                               | \$ |
| 9024020-C | Transport Protection, Rear                                         | \$ |
| 9900256-C | Tube Protection, Front Chassis                                     | \$ |
| 9900255-C | Tube Protection, Backrest Bar                                      | \$ |
| 9900261-C | Backpack                                                           | \$ |
| 110030D-C | reTyre Skins                                                       | \$ |
|           | 24"    25"    26"                                                  |    |
| 260010X-C | FreeWheel Attachment <i>Requires extended reinforced footrest.</i> | \$ |
|           | Height 6-11.5cm    Height 11.5-16.5cm                              |    |
| 905005X-C | Mototronik Quick Connector <i>Select seat width.</i>               | \$ |
|           | 33-36cm Width    39-42cm Width    45cm Width                       |    |

## Set-Up Adjustments

**Backrest Angle:**      Specify (mm) | A =  
15mm standard.


**Footrest Height:**      Specify (mm) | K =

## General

Permobil manufactures Panthera wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

## Notes

## Reference Diagrams




Balance Position



Backrest Placement



Backrest Angle Adjustment



Footrest Height