

Panthera S3 Swing - Order form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

*Required field

Special considerations that need to be addressed here (e.g., diagnosis):

Chassis

Standard configurations include front frame angle 105°, seat angle 4°, calf strap, standard colour option, backrest, carbon rear axle, front castors, swing-away footrests with foldable footplates, push to lock brakes and standard rear wheels with aluminium handrims. Seat height with 24" wheels: front 45cm, rear 43cm. Seat height can be adjusted with different size rear wheel or rear axle raising kit. 100kg weight limit for seat widths 33-42cm. 150kg weight limit for seat widths 45-50cm.

G547	Panthera S3 Swing	Seat width 36-50cm, seat depth 40cm.	\$
G553	Panthera S3 Swing Short	Chassis is 5cm shorter compared to S3 Swing. Seat width 33-39cm, seat depth 35cm.	\$
G5837	Panthera S3 Swing Low 0°	Chassis has flat 0° seat angle and lower seat height compared to S3 Swing. For seat width 36-45cm, seat depth 35, 37.5, 40. For seat width 50cm, seat depth 37.5, 40cm. 39cm front seat height with 24" rear wheels. <i>Extended model not kept in stock, longer leadtime.</i>	\$
G5818	Panthera S3 Swing Long	Chassis is 5cm longer compared to S3 Swing. Seat width 36-50cm, seat depth 37.5, 40, 42.5, 45cm. <i>Extended model not kept in stock, longer leadtime.</i>	\$

Chassis Colour

Standard Colours:	Anodic Black						Standard
Optional Colours:	Blue Classic	Mean Green	Panorange	Brave Bronze	White Pearl	Telemagenta	
	Grey Metallic	Racing Blue	Ruby Red	Sting Yellow	Horizon Grey	Pine Green	
	Pearl Violet						

Non-standard colours require additional delivery time of 20 days. Additional information [here](#).

Seat Measurements

Seat Width: 33cm 36cm 39cm 42cm 45cm 50cm

Check allowed seat width based on chassis selection.

Seat Depth: 35cm 37.5cm 40cm 42.5cm 45cm

Check allowed seat depth based on chassis selection.

Backrest

Standard Backrest	20-45cm range. 20cm back height not available for 50cm seat width.	Standard
Tapered Backrest 3cm	30-40cm range.	\$
Tapered Backrest 6cm	30-40cm range. <i>Additional delivery time of 20 days.</i>	\$
Flared Backrest 3cm	25-40cm range. <i>Additional delivery time of 20 days.</i>	\$
Curved Backrest	25-40cm range.	\$
Curved Tapered Backrest	35-40cm range. <i>Additional delivery time of 20 days.</i>	\$
Frame for Solid Backrests	35cm. Backrest frame still folds.	\$
Backrest Height:	20cm 25cm 30cm 35cm 40cm 45cm	
	<i>Check allowed height based on backrest selection.</i>	

Rear Wheels

Rear Axle

33023BB-C	Carbon Standard	Standard
33023BB-C	Carbon +15mm <i>Use this rear axle when choosing arm supports and/or increased seat height.</i>	No charge
33035BB-C	Carbon 4° Camber <i>Only for seat widths 33-45cm.</i>	\$
	Carbon 8° Camber <i>Only for seat widths 33-39cm.</i>	\$
Increased Seat Height		\$
10mm	20mm	30mm
40mm	50mm	

Rear Wheels

Standard Wheels:	22"	24"	Standard
Spinergy X Wheels:	24"	25"	26"
	Black Spokes	White Spokes	
One Arm Drive:	22"	24"	26"
	Aluminium Pushrim		Included
	Friction Pushrim		\$

Includes special rear axle. Additional delivery time of 5 days.

Tyres

Schwalbe Right Run 22"(501), 24"(540), 25" (559), 26"(590)	Standard
Rugged 24" (540)	\$
Schwalbe Marathon Plus 24"(540), 25" (559), 26"(590)	\$
Schwalbe One 24"(540), 25" (559)	\$
PU 24" <i>Not compatible with reTyre skins.</i>	\$

Handrims

190DD30-C	Aluminium	Standard
190DD11-C	Titanium	\$
19051DDBT-C	Paragrip	\$
19052DDBT-C	Tetragrip	\$
19030DD-C	Maxgrip 22" and 24" only.	\$

Front Wheels

Castor Selection:	90mm (3.5")	120mm (4.7")	150mm (6")
-------------------	-------------	--------------	------------

Footrest

670BB00-C	Foldable Footplates	Standard
672BB00-C	Foldable Footplates, Extended	\$
673BB00-C	Foldable Footplates, Large	\$
674BB00-C	Foldable Footplates, Large Extended	\$
	Heel Strap for Footplates	No charge

Wheel Locks

4630000-C	Push to Lock	Standard
4630020-C	Push to Lock, Extended Lever	\$
4630150-C	Wheel Lock Blades <i>For solid PU tyres. Not with one hand brakes.</i>	\$
367BBBB-C	One Hand Brake, Right Activation	\$
368BBBB-C	One Hand Brake, Left Activation	\$

Side Guards / Armrests

No Side Guards	Standard
Standard Side Guards	\$
Carbon Side Guards, Without Fender	\$
Carbon Side Guards, Small Fender	\$
Carbon Side Guards, Large Fender	\$
Side Guards With Arm Supports	\$

Push Handles

No Push Handles	Standard
4252060 VER2-C Height Adjustable Push Handles	\$
4252063 VER2-C Height Adjustable Push Handles, +10cm Additional 10cm height.	\$
4245000-C Foldable Push Handles	\$
42400BB VER2-C Push Bar	\$

Accessories

Anti-Tips

4353000-C Anti-Tips, Standard Seat Heights	\$
4355000-C Anti-Tips, Short/Low For short and low chassis.	\$

Cushions

572BB40-C Seat Cushion, 2.5cm	\$
575BB40-C Seat Cushion, 5cm	\$
57710BB-C Wedge Cushion Used to alter seat angle.	\$
990021X-C Seat Bag	\$
Mounted On Seat Upholstery Mounted Under Cushion	

Spoke Protectors | Only for 22" and 24" rear wheels.

4500D11-C Spoke Protectors, Clear	\$
4500D00-C Spoke Protectors, with Logo	\$
45003DD-C Spoke Protectors, Splash	\$
45001DD-C Spoke Protectors, Car Rim	\$
45004DD-C Spoke Protectors, Wild Animals	\$
45002DD-C Spoke Protectors, Smiley	\$

Torso Supports

6770305-C Torso Support Adjustable, Left	\$
6770306-C Torso Support Adjustable, Right	\$
4710012-C Backrest Wedge To fit behind back upholstery for additional lateral support.	\$

Other Accessories

9900251-C Tube Protection, Front	\$
9900251-C Tube Protection, Rear	\$
6810020-C Tray Table Only together with arm supports.	\$
9900301-C Crutch Holder S3	\$
9900261-C Backpack	\$
reTyre Skins	\$
24" 25" 26"	

Set-Up Adjustments

Backrest Angle: Specify (mm) | A =
14-28mm, 20mm standard.

Rear Axle Adjustment: Specify (mm) | B =
60mm standard.

Footrest Heights: Specify (mm) | A = Specify (mm) | B = Specify (mm) | K =
'B' only required for angle adjustable footrests.

Notes


Reference Diagrams



Backrest Angle



Rear Axle Adjustment



Footrest Height