



FUZE M50

2024 Order Form - NZ

Customer and Order Information

Purchase Order *	Tag	
Dealer/Provider Name *		
Purchaser Name and Contact Information *		
Shipping Address *		
User Information (if available)	Weight _____ kgs	Height _____ cm * required

Order Notes/Instructions

Base Chair

Fuze M50 50° Manual Tilt Wheelchair with Mid-Wheel Drive

FUZE-M50 Fuze 50° Manual Tilt Wheelchair Includes articulating rear casters and places drive wheels in the middle of the chair at the centre of gravity. \$

Frame Type & Weight Capacity

Fixed Width Frame Select a Fixed Seat Width below.

02010 136kg/300 lbs Weight Capacity

02011 181kg/400 lbs Weight Capacity Not compatible with Transit Tie Down System. Requires Flat-free Tyres or Heavy Duty Spoke Drive Wheels.

02025 **Adjustable Width Frame** 113kg/250 lbs Weight Capacity. Select an Adjustable Seat Width below. \$

Frame Options

02300 **Transit Tie Down System (TTDS)**

TTDS includes 4 factory-installed transport brackets. TTDS has been tested according to ANSI/RESNA WC 19. See Owners Manual for complete details. Tested to max. Weight Capacity of 136kg/300 lbs.

Seat Width Overall chair width is +8.5".

Fixed Seat Widths Order with Fixed Frame.

<input type="checkbox"/> 03015 15" Seat Width	<input type="checkbox"/> 03016 16" Seat Width	<input type="checkbox"/> 03017 17" Seat Width
<input type="checkbox"/> 03018 18" Seat Width	<input type="checkbox"/> 03019 19" Seat Width	
<input type="checkbox"/> 03020 20" Seat Width	<input type="checkbox"/> 03021 21" Seat Width	<input type="checkbox"/> 03022 22" Seat Width
<input type="checkbox"/> 03023 23" Seat Width		\$
<input type="checkbox"/> 03024 24" Seat Width		\$

Adjustable Seat Widths Order with Adjustable Frame.

<input type="checkbox"/> 03215 15"-19" Range - Set at 15"	<input type="checkbox"/> 03216 15"-19" Range - Set at 16"	<input type="checkbox"/> 03217 15"-19" Range - Set at 17"
<input type="checkbox"/> 03218 15"-19" Range - Set at 18"	<input type="checkbox"/> 03219 15"-19" Range - Set at 19"	

Seat Depth

Standard Seat Depths (Adjustable)

<input type="checkbox"/> 04015 15" Seat Depth	<input type="checkbox"/> 04016 16" Seat Depth	<input type="checkbox"/> 04017 17" Seat Depth
<input type="checkbox"/> 04018 18" Seat Depth	<input type="checkbox"/> 04019 19" Seat Depth	
<input type="checkbox"/> 04020 20" Seat Depth		

Custom Seat Depths Requires Stroller Handle Back - Straight. Not compatible with Reclining or Fold-down Back.

<input type="checkbox"/> 04021 21" Seat Depth	\$
<input type="checkbox"/> 04022 22" Seat Depth	\$

Frame Colour

PDG Premium Colours View colour swatches on PDG's website. pdgmobility.com/colours

- | | | |
|---|---|---|
| <input type="checkbox"/> 05001 Black Shimmer | <input type="checkbox"/> 05015 Sky Blue | <input type="checkbox"/> 05016 Poppy Red |
| <input type="checkbox"/> 05003 Fine Texture Black | <input type="checkbox"/> 05010 Lollypop Blue | <input type="checkbox"/> 05011 Lollypop Red |
| <input type="checkbox"/> 05017 Arctic Grey | <input type="checkbox"/> 05018 Tropical Teal | <input type="checkbox"/> 05019 Tangerine Orange |
| <input type="checkbox"/> 05020 Olive Green | <input type="checkbox"/> 05021 Lush Lilac | <input type="checkbox"/> 05012 Bubble Gum Pink |
| <input type="checkbox"/> 05022 Leaf Green | <input type="checkbox"/> 05023 Boysenberry Purple | <input type="checkbox"/> 05014 Daffodil Yellow |

05200 **Custom Colour** Contact Permobil for specifications.

\$

Armrests

06000 **No Armrests**

Height Adjustable Armrests (Single Post)

\$

06114 10"-14" Range 06116 12"-16" Range

06100 **Fixed Height Armrests (12")** Full Length Armpad only.



Armrest Height If provided, Armrests will be set to this height by PDG. If left blank, the default setting will be used.

Seat Pan to top of Armrest: _____ Does not include height of Armpad



Armrest Position Indicate preference for Armrest position. If left blank, the default setting will be used.

In the interest of maintaining chair stability and avoiding interference with other chair functions, PDG may not be able to accommodate stated preference.

06401 Rear Default Position 06403 Front

Locking Flip-up Armrests Not recommended for users with excessive strength or tone.

\$

Not compatible with Fold-down Back. Only available with Regular Full Length and Desk Length Armpads.

06310 Set at 10" 06311 Set at 11" 06312 Set at 12"

06313 Set at 13" 06314 Set at 14"

Armpads

Regular Armpads

07100 Desk Length - 2" x 10" 07101 Full Length - 2" x 14"

Waterfall Armpads

\$

07105 Desk Length - 2" x 10" 07106 Full Length - 2" x 14"

Wide Armpads May interfere with wheels in some tilt configurations.

\$

07110 Desk Length - 3" x 10" 07111 Full Length - 3" x 12"

Gel Armpads May interfere with wheels in some tilt configurations.

\$

07120 Full Length - 2" x 12" 07121 Full Length - 3.5" x 12"

Armrest Options

07210 **Armrest Wideners** Move Armrests outward by 0.5" on each side for 1" total additional width. Seat size does not change.

\$

Seat Pan

08000 **No Seat Pan**

08001 **Solid Aluminium Seat Pan**

\$

Seat Height

Don't forget **Select desired Seat-to-Floor Height** STF Height is adjustable within the range of the chosen Axle Plate (Low/Standard). RC refers to the size of the Rear Casters.

Chairs with Swing-Away Front Rigging are set up with up to 5° of pre-tilt to avoid interference. This increases the actual STF Height at the front of the seat by up to 1" compared to the selected value. Easy-Out is not affected.

5" Front Casters									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				5"
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	Std. STF			6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>					6"

6" Front Casters & 6" Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				5"
22"		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			Std. STF	6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"

6" Front Casters & 8" Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>			Std. STF	5"
22"			<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"

6" Front Casters & Frog Legs Suspension Forks									
RW	13"	14"	15"	16"	17"	18"	19"	20"	RC
20"				<input type="checkbox"/>	<input type="checkbox"/>				5"
22"				<input type="checkbox"/>	<input type="checkbox"/>			Std. STF	6"
24"	Low STF	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>				6"

Front Casters

Urethane Casters

 Select Front Caster size here. Rear Caster size as per STF tables.

<input type="checkbox"/> 11051	5" x 1.25" Caster size	<input type="checkbox"/> 11061	6" x 1.25" Caster size	
<input type="checkbox"/> 11063	6" x 1.50" Caster size			\$
<input type="checkbox"/> 11260	Soft-Roll Casters	6" x 1.50" Caster size. Urethane tire & Silver Anodized Aluminium hub. Available as Front Casters only.		\$
<input type="checkbox"/> 11400	Frog Legs Suspension Forks	Requires 6" x 1.50" Soft-Roll Casters. Suspension Forks used for Front Casters only.		\$

Drive Wheels

Composite Wheels

<input type="checkbox"/> 12020	20" Wheel size	<input type="checkbox"/> 12022	22" Wheel size	<input type="checkbox"/> 12024	24" Wheel size	
Spoke Wheels						\$
<input type="checkbox"/> 12120	20" Wheel size	<input type="checkbox"/> 12122	22" Wheel size	<input type="checkbox"/> 12124	24" Wheel size	
<input type="checkbox"/> 12224	24" Heavy Duty Spoke Wheels					\$

Axle Type

<input type="checkbox"/> 13000	Quick Release Drive Wheel Axles	
<input type="checkbox"/> 13010	Fixed Drive Wheel Axles	

Drive Wheel Tires

<input type="checkbox"/> 14000	Urethane Tires		
<input type="checkbox"/> 14010	Pneumatic Tires	May result in a finished Seat Height 1/2" greater than stated.	\$
<input type="checkbox"/> 14020	Pneumatic Tires with Flat-Free Inserts	May result in a finished Seat Height 1/2" greater than stated. Not available with Spoke Wheels.	\$

Handrims

<input type="checkbox"/> 15000	No Handrims		
<input type="checkbox"/> 15010	Standard Aluminium Handrims		
<input type="checkbox"/> 15020	Plastic Coated Handrims	\$	
<input type="checkbox"/> 15030	Projection Handrims	Composite Wheels only. 10 Projections, 30°, 1-1.8" Length.	\$

Wheel Locks

User Wheel Locks

<input type="checkbox"/> 16010	Push to Lock	<input type="checkbox"/> 16020	Pull to Lock		
User-activated Unilateral Wheel Lock				Pull to Lock. Available up to 22" Seat Width only.	\$
<input type="checkbox"/> 16310	Left Activation	<input type="checkbox"/> 16311	Right Activation		

Wheel Lock Options

<input type="checkbox"/> 16100	Wheel Lock Extensions	\$
--------------------------------	------------------------------	----

Back Posts

Back Post Style

<input type="checkbox"/> 18100	Stroller Handle Back	<input type="checkbox"/> 18105	Cane Handle Back	Not compatible with Headrest Mounting Plate.	
Back Post Type					
<input type="checkbox"/> 18200	Straight	<input type="checkbox"/> 18205	8° Bend		
Back Post Height Measured from the Seat Pan to the Handle.					
<input type="checkbox"/> 18420	Short (17")	<input type="checkbox"/> 18425	Tall (23")		
Back Post Folding Fold-down Back is not compatible with Reclining Backrest.					
<input type="checkbox"/> 19220	Fold-down	<input type="checkbox"/> 19221	Non-folding		
Seat to Back Angle					
<input type="checkbox"/> 19090	Set at 90°	<input type="checkbox"/> 19096	Set at 98°	<input type="checkbox"/> 19102	Set at 106°
<input type="checkbox"/> 19108	Set at 114°	<input type="checkbox"/> 19114	Set at 122°		

Back Options

Angle/Height Adjustable Stroller Handle Not compatible with Headrest Mounting Plate or Cane Handle Back. \$

19010 Button on Outside 2819011 Button on Inside

Endomorph Back Backrest will be narrower than Seat by total of 2x the Offset selected below. Minimum width of Backrest is 12".
Not compatible with Flip-up Armrests, Reclining Backrest. \$

18102 1" Offset per side 15118103 1.5" Offset per side 2818104 2" Offset per side

19212 **Headrest Mounting Plate** Not compatible with Angle/Height Adjustable Stroller Handle, Cane Handle Back. Not required if mounting to a solid Back. \$

Back Upholstery

20000 **No Back Upholstery** \$

Standard Upholstery Not available with Cane Handle Back. \$

20418 Short (18" Height) 2820421 Regular (21" Height)

20450 **Tension Adjustable Upholstery** Height of Upholstery is 16". Not available with Cane Handle Back. \$

Position of Upholstery Select the desired position of the top of the Upholstery relative to the Seat Pan.

17" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>								

23" Back Post									
Style	17"	18"	19"	20"	21"	22"	23"	24"	25"
Short	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Regular					<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
Tension Adj				<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

Tilt Activation

Attendant Tilt only Select trigger mounting location below.

23011 Right Back Post 23012 Left Back Post

Tilt Limit Optional. If no tilt limit is selected, the full tilt range will be available.

23220 Set at 20° 23230 Set at 30° 23240 Set at 40°

Tip Configuration (Front Caster Lift)

Attendant Tip only Select trigger mounting location below.

23451 Right Back Post 23452 Left Back Post 23453 Omit Trigger No Front Caster Lift

Tip Clearance Limit (Front Caster Lift) Actual tip may vary by up to 0.5". No selection required if Trigger is omitted.

23455 4" Caster Lift Default 23456 2" Caster Lift

Recline

23310 **Reclining Backrest (Manual)** \$

30° maximum recline. Not compatible with Cane Handle Back, Fold-down Back. Not available with 21" & 22" Seat Depth.
Finished Back Post Height is 3" taller than selected Back Post Height.

Front Rigging

25000 **No Front Rigging** 25001 **Sector Blocks only**

Easy-Out Front Rigging Effortless insertion and removal. No Hanger Bracket remains on frame. pdgmobility.com/easy-out

25410 Fixed 70° \$

25420 Universal Angle Adjustable (0°-110°) Set at angle: _____ (use increments of 10°) \$
Adjust angle by loosening one bolt. Tool-free leg-length adjustment. Front Rigging can interfere with casters in the 90-110° settings.
Optional for Universal AA: 25421 AA High-mount Option For shorter Lower Leg Lengths of 8"-13".

Swing Away Front Rigging Legacy Design. Not compatible with Easy-Out, Push to Lock Wheel Locks. Requires pre-tilt of up to 5°.

25130 Swing Away 60° 25140 Swing Away 70°

25200 Swing Away Elevating Legrests Includes Calf Pads. \$

Footplates

26010 **Composite Footplates**

26020 **Aluminium Angle Adjustable Footplates** \$


26040 **One Piece Flip-up Footplate** Requires selection of Composite or Aluminium Angle Adjustable. \$

26050 **Interlocking Angle Adjustable Footplates** Available for Seat Widths up to 21" only. \$

Footrest/Legrest Accessories

Calf Pad Not compatible with Gel Padding.		
<input type="checkbox"/>	27012 For Easy-Out Front Rigging Includes reinforced premium Calf Pad mount. Requires pre-tilt of up to 5°.	\$
<input type="checkbox"/>	27010 For Swing Away Front Rigging	\$
<input type="checkbox"/>	27020 Gel Padding Not compatible with Calf Pad or Easy-Out Front Rigging.	\$
Gel Knee Pad		
<input type="checkbox"/>	27032 For Easy-Out Angle Adjustable Front Rigging	\$
<input type="checkbox"/>	27030 For Swing Away ELR Front Rigging	\$
<input type="checkbox"/>	27040 Length Adjustable Heel Loops	\$
<input type="checkbox"/>	27050 Bodypoint Padded Calf Strap Sizes greater than 22" are not padded.	\$

Front Rigging Extension Tubes

Regular Extension Tubes No selection required for Easy-Out AA High-mount Option.		
<input type="checkbox"/>	28010 Short <u>14"-18"</u> for Easy-Out Front Rigging <u>13"-17"</u> for Swing Away Front Rigging	
<input type="checkbox"/>	28020 Long <u>18"-24"</u> for Easy-Out Front Rigging <u>15"-19"</u> for Swing Away Front Rigging	\$
	Lower Leg Length If provided, Front Rigging will be set to this length by PDG. If left blank, the default setting will be used. Seat Pan to Foot Plates: _____ Measured along the length of the Front Rigging. Consider the height of your Seat Cushion.	

Miscellaneous Options

Spoke Guards Available for 20"-24" Composite Drive Wheels. Made from durable and lightweight aluminium.		\$	
<input type="checkbox"/>	30020 Black	<input type="checkbox"/>	30021 Frame Colour
Stealth Flip Down Abductor Includes wedge-style pad. Requires Solid Seat Pan.		\$	
<input type="checkbox"/>	30060 Small Pad 3"x4"x3"x2" size	<input type="checkbox"/>	30062 Large Pad 5"x5"x4"x3" size
<input type="checkbox"/>	30080 Gel Hip Guides Reduces available Seat Width by 2".	\$	
<input type="checkbox"/>	30090 Hip Belt Mount Provides improved mounting options for Bodypoint Hip Belt or other third-party belts.	\$	
<input type="checkbox"/>	30152 Ventilator Bracket Compatible with LTV1 (laptop) ventilator.	\$	
<input type="checkbox"/>	30160 Clear Lap Tray (Slide On) Available for Seat Widths of 17"-22" only. Will interfere with Drive Wheel when in tilt.	\$	
<input type="checkbox"/>	30200 Q'Strain Lap Belt Model S5-6329-11. For use with WC-19 Transit Tie Down System. 24"-51" Belt Length.	\$	
<input type="checkbox"/>	30222 Medical Necessities Backpack Black with Permobil Logo.	\$	
<input type="checkbox"/>	30170 Adhesive Back Hook Strip, 2" x 30" For Cushion mounting.	\$	

Order Notes