

TiLite T Series - Order form



Dealer Name*:		Account #*:
City, Post Code, Country:		
Phone:	Email:	P.O. Number
Clinician:	Client Reference*:	Client Weight:
Rep Name:		
<div>*Required field</div> Special considerations that need to be addressed here (e.g., diagnosis):		

Frame Style

ATFS1	Aero T Aluminum Frame Full adjustability. 265lb (120kg) weight limit.	\$
TAFS1	TRA Titanium Frame Full adjustability. 265lb (120kg) weight limit.	\$
TRFS1	TR Titanium Frame Minimal adjustability. 265lb (120kg) weight limit.	\$

Frame Options

THDR1	Heavy Duty Frame Not with TR titanium Frame. 350lb (158kg) weight limit. Includes tension adjustable by straps back upholstery. Requires folding aluminum adjustable height backrest, billet aluminum front wheels, Shadow or Spinergy SPOX rear wheels. LX rear wheels only with 0° or 2° camber. Not with Slipstream fork, velcro adjustable upholstery or Glide™ forks.	\$
TQST1	Q'Straint Adaptable Frame Not with TR titanium Frame. 350lb (158kg) weight limit. Includes tension adjustable by straps back upholstery. Requires folding aluminum adjustable height backrest, billet aluminum front wheels, Shadow or Spinergy SPOX rear wheels. LX rear wheels only with 0° or 2° camber. Not with Slipstream fork, velcro adjustable upholstery, open loop footrest, V-front taper or Glide™ forks.	\$
TPWR1	Power Adaptable Frame Requires LiteSpeed billet aluminum front wheel and 0° camber. Not with single-sided forks or titanium camber tubes. Not required for SmartDrive.	\$
TRFR1	Reinforced Frame	\$

Frame Finishes

TFF1	Satin Finish Available with titanium frames only.	No charge
SCBT1	Buffing Pads for Satin Frame	\$
TFF3	Tattoos Tattoos on front frame tubes. Titanium frame finished in satin. Tattoos will be trimmed and adjusted to fit your frame.	\$
	TFF45 Need for Speed TFF44 Flower Garden TFF43 Mosaic	
TFF	Powder Coated Frame Select below.	Standard

Frame Color:	TFF8 High Gloss Black	TFF10 Red Pepper	TFF29 Neon Pink	TFF14 Ocean Blue Metallic
	TFF15 Tangerine Metallic	TFF30 Acid Green Metallic	TFF17 Midnight Blue Pearl	TFF18 White Pearl
	TFF19 Black Lava	TFF27 Titanium Grey	TFF59 Pomegranate Red	TFF62 Blue Mist
	TFF57 Indigo Purple	TFF40 Hyper Blue	TFF41 Flame Orange	TFF42 OD Green
	TFF60 Silver Streak	TFF63 Lavendar Dusk	TFF58 Black Galaxy	TFF61 Seafoam Green
	TFF64 Daffodil Yellow			



Need for Speed



Flower Garden



Mosaic

Color Anodize Package

Standard:	Black	Standard
Anodized Package:	Gold Green Purple Orange Blue Red	\$

Includes caster mount, caster fork and caster caps. If no anodize colour is selected, default will be black.

TiFit: Tailored Individualized Frame

Measurements shown below are our most popular selections. Not seeing what you want? Contact customer service or view our Creative Mobility Division (CMD) catalog [here](#) to view additional options. Please refer to Reference Diagrams at end of form for more information.

Rear Seat Width 'A' | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post.

TSW1:	12"	13"	14"	15"	16"	17"	18"	19"	
ZSW2:	20"								\$

>20" sizes available through CMD.

Front Seat Width (Seat Taper) 'T' | Measure from inside of front frame tube to inside of opposite front frame tube. Required with V-front end, not with standard front end.

TFSW1:	No Taper		Standard
TFSW2:	1" Taper	3" narrower than Rear Seat Width 'A'. Includes tension adjustable by straps seat upholstery. Bolt on adjustable sling available through CMD #20.	\$
TFSW4:	2" Taper	4" narrower than Rear Seat Width 'A'. Includes tension adjustable by straps seat upholstery. Bolt on adjustable sling available through CMD #20.	\$

Seat Depth 'B' | Measure from front of back post to the front edge of seat sling. Frame bend will start 0.5" in front of seat sling unless Custom Frame Depth is specified or the combination of specifications affect the rolling dynamics of the chair.

TSD1:	12"	13"	14"	15"	16"	17"	18"	19"	
TSD2:	20"								\$
TSD3 Custom Frame Depth:	Specify:								No charge

Measure from front of seat sling to beginning of frame bend. Ranges from 0.5"-2" and can be given as specific as 0.25" increments. A custom frame depth outside the specified range may be available through our Creative Mobility Division. A custom frame depth of 20" will have an additional charge.

Ergonomic Seat 'U' | Measure from front of back post to beginning of upward bend of seat tube.

TERG1 TR Frame:	5"	6"	7"	8"	9"	\$
-----------------	----	----	----	----	----	----

Tension adjustable by strap comes standard. Bolt on adjustable sling available through CMD #20. Solid seat available through CMD #16

TERG2 Aero T/TRA Frame:	6"	7"	8"	\$
-------------------------	----	----	----	----

*Not with V-front or seat taper, seat back that locks when folded, desk armrests, depth adjustable backrests, fold down rigid sideguards, fendered adjustable sideguards or Curve cushion. **Bolt on adjustable sling available through CMD #20. Solid seat available through CMD #16.** Minimum seat depth 12" for 6"-7". Minimum seat depth 13" for 8".*

Front Seat Height 'C' | Measure from floor to top of seat tube at front edge of seat sling.

[illegible]

Rear Seat Height 'D' | Measure from floor to top of seat tube at back post.

TRSH1:	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"
	18.5"	19"	19.5"	20"								

Rear Wheel	26"	15.5"-20"
	25"	15"-20"
	24"	14.5"-20" (15" on TRA)
	22"	13.5"-20"
	20"	12.5"-20"

Rear Seat Height Range

Front Angle or Overall Frame Length | *Select one. If Overall Frame Length is selected, front angle will be bent to achieve requested length.*

TFA1	Frame Angle 'L' <i>Measure from frame front to floor. 85° not with 6" front wheels. 90° not with 5" or 6" front wheels.</i>					Standard
	70°	75°	80°	85°	90°	

TFA2	Overall Frame Length 'M'	<i>Measure from front of backrest to front edge of footrest.</i>	No charge
	Specify:		

Centre of Gravity 'H' | Measure from front of back post at the seat tube to the centre of the rear axle. If no centre of gravity is selected the chair will ship with the centre of gravity in the most rearward position. Ergo seat selection may limit available centre of gravity.

TCOG1:	1"	1.25"	1.5"	1.75"	2"	2.25"	2.5"	2.75"	3"	3.25"	3.5"	3.75"
	4"	4.25"	4.5"	4.75"	5"							

Rear Wheel Spacing 'J' | Measure from outside of back post to inside of rear tire. TRA 0.75" rear wheel spacing not with 0° or 2° camber.

TRWS1	Fixed Rear Wheel Spacing	Standard
	0.75" 1" 1.25"	
TCBR12	Adjustable Rear Wheel Spacing Camber Plugs	No charge
TCBR9	Additional Spacers to Increase Rear Wheel Spacing	No charge
TCBR11	Zero Play Aluminum Camber Plugs Not with power adaptable frame or more than 6° camber.	No charge

Camber 'K'

TCBR1:	0° 2° 4° 6°	
	8° camber may be available through our Creative Mobility Division	
TCBR15	Aluminum Camber Tube	Standard
TCBR3	Titanium Camber Tube Not with power adaptable frame.	\$

Footrest**Footrest Width** | Select standard front end 'I' if no seat taper 'T' has been selected. Select V-front end 'V' if a seat taper 'T' has been selected.

TFE1	No V-Front Taper 'I' Find your rear seat width 'A' in the chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube.	Standard
------	--	----------

Rear Seat Width 'A'	12"	10"	8"	
	13"	11"	9"	
	14"	12"	10"	8"
	15"	13"	11"	9"
	16"	14"	12"	10"
	17"	15"	13"	11"
	18"	16"	14"	12"
	19"	17"	15"	13"
	20"	18"	16"	14"
		No Taper	2" Taper	4" Taper
Footrest Width 'I'				

TFE2	V-Front Taper 'V' Measure from inside of front frame tube to inside of opposite front frame tube 2.5" above footrest. Titanium footrests have 1" of adjustment up and down. Not with ergonomic seat frames.	\$
	9" 10" 11" 12" 13" 14" 15"	
CMD008	Custom Abducted Front Frame Width Adapted footrest and/or front seat width. Please specify desired front seat width and desired footrest width to the nearest 0.5" in the Notes section. Some limitations may apply.	\$

Standard Footrests

Seat to Footrest 'E':	13.5" 14" 14.5" 15" 15.5" 16" 16.5" 17" 17.5" 18" 18.5" 19"	
	Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height 'C'.	
TFTR3	Angle Adjustable Footrest	Standard
TFTR1	Titanium Footrest w/Flat ABS Cover	\$
TFTR2	Titanium Open Loop Footrest	\$
TFTR12	Angle Adjustable Flip-Back Footrest Not with 10" or 11" rear seat width 'A', heavy duty frame or ABS cover with 2" sides.	\$

High Mount Footrests

Seat to Footrest 'E':	7"-9" 9"-11" 11"-13" Set at	
	Measure from front edge of seat sling to top rear of footrest.	
TFTR4	Angle Adjustable High Mount Footrest	\$
CMD040	Custom Two-Piece Flip-Up Footplates Please visit our CMD catalog and specify footrest type and rigidizer bar placement in Notes section of this order form. Some limitations may apply.	\$

Residual Limb Support

TRSL1	Residual Limb Support, 8"W x 10"L	\$
	Left Right	
TRSL2	Residual Limb Support, 8"W x 10"L, Pair	\$

Calf Straps

TCLF2	Velcro Adjustable				\$
	Small (9"-12" footrest widths)	Medium (11"-15" footrest widths)	Large (14"-17" footrest widths)		
	Xtra Large (16"-19" footrest widths)				
SP102x	Bodypoint® Padded Velcro Adjustable				\$
	Medium (11"-15" footrest widths)	Large (16"-19" footrest widths)			
TCLF1	Buckle Adjustable, Black Nylon Webbing				\$

Backrest

Back Type

TSB1	Folding Aluminum Adjustable Height Back	<i>Required for heavy duty frame.</i>	Standard
TSB2	Folding Titanium Adjustable Height Back, Rounded Rigidizer Bar	<i>Titanium back will match frame finish.</i>	\$
TSB15	Non-Folding Fixed Height Titanium Back, Rounded Rigidizer Bar	<i>Not with Aero T, TRA or push handles. Titanium back will match frame finish.</i>	\$

Back Options

TSB11	Backrest Locks When Folded	<i>Available on all backs. Not with ergonomic seat.</i>	Standard
TSB12	4" Deep Backrest Rigidizer Bar	<i>Not a rounded style rigidizer bar. Not with non-folding, fixed height titanium back type.</i>	No charge
TSB3	TiShaft™ Back Release Bar, Aluminum	<i>Not with non-folding or fixed height titanium backs.</i>	Standard
TSB5	TiShaft™ Back Release Bar, 2x2 Twill Weave Carbon Fiber	<i>Not with non-folding or fixed height titanium backs.</i>	\$

Push Handles

TSB6	Integrated Push Handles		No charge
TSB7	Integrated Fold-Down Push Handles	<i>Requires tension adjustable by straps back upholstery.</i>	\$
TSB8	Bolt-On Push Handles	<i>Includes 4" deep backrest rigidizer bar.</i>	\$
TSB9	Bolt-On Height-Adjustable Push Handles	<i>Includes 4" deep backrest rigidizer bar.</i>	\$

Folding Back Height 'F' | TSB1/TSB2. Set at 0.5" increments.

Push Handles	None / Bolt-On	8.5"-11.5"	12"-16.5"	16"-20.5"	Set at:
	Integrated	10.5"-12.5"	13.5"-16.5"	17.5"-20.5"	Set at:
	Fold-Down	13"-15.5"	14"-18.5"	18.5"-21"	Set at:
		Short	Medium	Tall	Specify
Folding Back Height					

Fixed Back Height 'F'

Select height:	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"	9"	9.5"	10"	10.5"
	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"
	17"	17.5"	18"	18.5"	19"	19.5"	20"					

Required for fixed side guards. Not with swing-away tubular or swing-away flip-back tubular armrests, integrated push handles, integrated fold-down push handles or bolt-on height adjustable push handles.

Back Angle 'G' | TSBA1 Measure angle from front of back post to floor.

Folding Back:	80°	83°	86°	89°	92°	95°	98°	101°						
Fixed Back:	80°	81°	82°	83°	84°	85°	86°	87°	88°	89°	90°	91°	92°	93°
	94°	95°	96°	97°	98°	99°	100°	101°						

Back Supports

Back Upholstery

TBUP2	Velcro Adjustable	<i>Not with fold-down push handles.</i>	Standard
TBUP3	Tension Adjustable By Straps, Padded Nylon	<i>Required for fold-down push handles.</i>	\$
TBUP5	Tension Adjustable By Straps, AIR Breathable	<i>Required for fold-down push handles.</i>	\$
TUPC3	Ultrasuede Backrest Stripe, Black	<i>Not with AIR breathable upholstery.</i>	\$
TBUP0	Omit Back Upholstery		No charge

ROHO® Agility® Active Aluminum Contour Back Support

AG2-xxF-2A	Agility Active Aluminum, Min Contour with Foam	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$
AG3-xxF-2A	Agility Active Aluminum, Mid Contour with Foam	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$

Shell Height: 8" 10"

ROHO® Agility® Active Carbon Contour Back Support | 2x2 Twill Weave Carbon Fiber

AG2-xxA-2CB	Agility Carbon Min Contour Back with Foam/Air Insert	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$
AG3-xxA-2CB	Agility Carbon Mid Contour Back with Foam/Air Insert	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$

Shell Height: 10" 13" 16"

Comfort Acta-Embrace® Back Support

AE-BSxWxL	Acta-Embrace Back with 2-Point Quick Release Compass Hardware		\$
-----------	---	--	----

Shell Height: 10" 12" 14" 16"

Accessories: Privacy Shield \$

Comfort Acta-Relief Back Support

AR-BSxWxL	Acta-Relief Back with 2-Point Quick Release Compass Hardware		\$
-----------	--	--	----

Shell Height: 10" 12" 14" 16"

Accessories: Privacy Shield \$

Seat Options

Seat Upholstery

ZSUP1	Tension Adjustable Bolt-On Upholstery, Black	Not with ergonomic seat or tapered seat.	Standard
ZSUP2	Tension Adjustable By Straps, Black	Standard with ergonomic seat and seat taper of 1"-3".	\$
ZSUP8	Aluminum Solid Seat	Ergonomic seat requires CMD016. Not with 10" or 11" rear seat width.	\$
ZSUP0	Omit Seat Upholstery		No charge
CMD018	Custom iSleeve	Please provide sleeve width and depth (inches) in Notes section at the end of this order form. Add 1" to device depth and 0.25" to width. Some limitations may apply.	\$

Seat Cushion

CU-SV-xx	Comfort Curve Cushion		\$
M2-FQA-xx	Comfort M2 Anti-Thrust Cushion		\$
EM-FQA-xx	Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion		\$
ALxxHD	ROHO® AirLite Cushion		\$
1Rxx	ROHO® Single Compartment Cushion	Fits seat widths 14"-19". Low Profile (LPC) High Profile(H-CA-SR)	\$
QSxx	ROHO® QUADTRO SELECT® Cushion	Fits seat widths 14"-19". Low Profile (LPC) High Profile (C)	\$
1RHExxC-SR	ROHO® Hybrid Elite Single Valve Cushion	Fits seat widths 14"-18", and 20"	\$

Precision Lock Front Forks

TFK1	TiLite Standard Fork		Standard
TFK2	TiLite Slipstream Single-Sided Fork	Not with power adaptable frame.	\$
TGL1	Glide™ Suspension Fork	Black anodize only. Select rider weight. 125lbs 125-180lbs 180-265lbs	\$
TTFS1	Titanium Fork Stem		\$

Front Wheels

TFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$
	Silver Hub Black Hub	
TFW14	4"x0.75" Plastic Wheel w/Poly Tire <i>Not with Slipstream fork or power adaptable frame.</i>	No charge
TFW3	4"x0.75" Light-Up Micro Caster <i>Not with Slipstream fork or power adaptable frame.</i>	\$
TFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with power assist, power add-on system or heavy duty frame.</i>	\$
TFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$
	Silver Hub Black Hub	
TFW6	5"x1" Plastic Wheel w/Poly Tire <i>Not with 90° front angle, Slipstream fork or power adaptable frame.</i>	No charge
TFW28	5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with 90° front angle or power adaptable frame.</i>	\$
TFW30	5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not with 90° front angle.</i>	\$
	Silver Hub Black Hub	
TFW9	6"x1" Plastic Wheel w/Poly Tire <i>Not with 85° or 90° front angle, Slipstream fork or power adaptable frame.</i>	No charge
TFW32	6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with 85° or 90° front angle, Slipstream fork or power adaptable frame.</i>	\$
TFW33	6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not with 85° or 90° front angle or Slipstream fork.</i>	\$
	Silver Hub Black Hub	

Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

TRW14	TiLite SHADOW™ <i>25" and 26" wheels require 15" minimum seat depth.</i>	Standard
Wheel Size:	20" (451) 22" (501) 24" (540) 25" (559) 26" (590)	
TRW7	Spinergy SPOX <i>25" and 26" wheels require 15" minimum seat depth. Orange, pink, purple and green spokes only available in 20" size.</i>	\$
Wheel Size:	20" (451) 22" (501) 24" (540) 25" (559) 26" (590)	
Spoke Color:	Black Red Blue White Orange Pink Purple Green	
TRW8	Spinergy LX <i>25" and 26" wheels require 15" minimum seat depth.</i>	\$
Wheel Size:	22" (501) 24" (540) 25" (559) 26" (590)	
Spoke Color:	Black Red Blue White Orange Pink Purple Green	
TRW9	Spinergy Design Lab LX Wheels <i>Black spokes only. If same color is selected, anodized hubs and rims match colour anodized package. 25" and 26" wheels require 15" minimum seat depth.</i>	\$
Wheel Size:	22" (501) 24" (540) 25" (559) 26" (590)	
Anodize Color:	Gold Green Blue Red Purple Orange	
TRW16	Spinergy Carbon Blade <i>2x2 Twill Weave Carbon Fiber. Not with heavy duty frame. 25" wheels require 15" minimum seat depth.</i>	\$
Wheel Size:	24" (540) 25" (559)	
TRW17	X-Core Mag <i>Not with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i>	\$
Wheel Size:	20" (451) 22" (501) 24" (540)	

Axles

TAXL1	Stainless Quick Release	Standard
TAXL2	Titanium Quick Release	\$
TAXL3	Quad Release <i>Not with Mag or 26" (590) Spinergy Sprox or LX rear wheels.</i>	\$

Pneumatic Rear Tires

TRTR1	Treaded 22"(501), 24"(540), 26"(590)	No charge
TRTR2	Airless Inserts <i>Not with Spinergy Carbon Blade rear wheels.</i>	\$
TRTR31	MBL SpeedLite 20"(451), 22"(501), 24"(540), 25"(559), 26"(590)	No charge
TRTR30	Primo Sentinel 24"(540), 25"(559), 26"(590)	\$
TRTR8	Schwalbe® Marathon Plus Evolution 24"(540), 25"(559), 26"(590)	\$
TRTR32	Primo C-764 Knobby 24"(540)	\$

Solid Rear Tires

TRTR22	Primo Orion 20"(451), 22"(501), 24"(540), 25"(559) <i>Not with Spinergy Carbon Blade rear wheels.</i>	\$
TRTR23	Primo Active 24"(540), 25"(559) <i>Not with Spinergy Carbon Blade rear wheels.</i>	\$

Handrims

TRIM1	Aluminum, Silver Anodized	Standard
TRIM21	Aluminum, Black Anodized	No charge
TRIM12	Titanium, Satin Finish <i>Not with 20"(451) or 22"(501) rear wheels or 24"(540) knobby tires.</i>	\$
TRIM20	Q-Grip <i>Exclusive high-friction, anti-slip coating. Not with 20"(451) rear wheels.</i>	\$
TRIM5	Plastic Coated, Black	\$
TRIM10	The Surge <i>Large profile. All-In-One oval with Gription™ strip.</i>	\$
TRIM11	The Surge LT <i>Small profile. All-In-One oval with Gription™ strip.</i>	\$
TRIM14	The Natural-Fit®, Standard Grip <i>Original, large profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
TRIM15	The Natural-Fit®, Super Grip <i>Original, large profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
TRIM16	The Natural-Fit®, Without Thumb Grip <i>Original, large profile.</i>	\$
TRIM17	The Natural-Fit® LT, Standard Grip <i>Reduced weight, small profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
TRIM18	The Natural-Fit® LT, Super Grip <i>Reduced weight, small profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
TRIM19	The Natural-Fit® LT, Without Thumb Grip <i>Reduced weight, small profile.</i>	\$
TRIM22	Omit Handrims	No charge

Handrim Tab Length

TRIM7	Long Tabs	Standard
TRIM8	Short Tabs <i>Not with Mag rear wheels or The Natural-Fit® handrims with thumb grip.</i>	\$

Wheel Locks

TWLK7	Compact Composite Push to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
TWLK8	Compact Composite Pull to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
TWLK10	Extension Handles for Compact Composite Push/Pull to Lock	\$
TWLK11	MBL Composite Scissor Lock	\$
TWLK3	Uni-Lock	\$
	Push Pull	
TWLK4	Extension Handles for Uni-Lock	\$
TWLK15	Grade Aids Push to Lock <i>Not with solid rear tires.</i>	\$

Side Guards

TSDG2	Aluminum Rigid Removable	\$
	Small Standard	
TSDG3	Aluminum Quad Removable, with Hole	\$
	Small Standard	
TSDG4	Aluminum Fold Down Rigid <i>Not with ergonomic seats.</i>	\$
	Small Standard	
TSDG7	Carbon Rigid Removable	\$
	Small Standard	
TSDG8	Carbon Quad Removable, with Hole	\$
	Small Standard	
TSDG9	Carbon Fold Down Rigid <i>Not with ergonomic seats.</i>	\$
	Small Standard	
TSDG6	Aluminum Fendered Adjustable <i>Available with 24"-25" rear wheels. Requires rear seat height 14.5"-18.5"</i>	\$
TSDG12	Carbon Fendered Removable	\$
TSDG13	Side Guards Custom Cut to Profile of Rear Wheels	\$
CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 0.5" and side guards type (rigid or fold-down) in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>	\$

Armrests

TARM1	Swing Away Tubular	<i>Not with non-folding, fixed height titanium seat back or height adjustable push handles.</i>	\$
Range:	6.5"-9.5"	9.5"-12.5"	
TARM4	Swing Away Flip Back Tubular 8.5"-11.5"	<i>Not with non-folding, fixed height titanium seat back.</i>	\$
TARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard	<i>Requires 14" minimum seat depth. Not with ergonomic seat.</i>	\$
Range:	9"-12"	12"-15"	
Pad Size:	Desk Length (10")	Full Length (14")	
TARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop	<i>Requires 15" minimum seat depth. Not with ergonomic seat.</i>	\$
Range:	9"-12"	12"-15"	
Pad Size:	Desk Length (10")	Full Length (14")	
TARM5	Armrests Covered in Black Ultrasuede	<i>Not with desk arms.</i>	\$

Positioning Belts

TBLT1	Velcro Adjustable, Black Nylon Webbing	\$
HB205-xxx-B1	Bodypoint® Two Point Push-Button Buckle (Auto Style)	\$
	Medium Large	
EB205-xxx	Bodypoint® EvoFlex® Pelvic Stabilizer, Push Button	\$
	Medium Large	
HW320-25-2	Bodypoint® Optional Band Clamp Mounting Hardware	\$

Anti-Tips

TTIP1	Rear Aluminum, Black	\$
TTIP2	User Friendly Flip-Up	\$

Accessories

TGRD1	Clear Lexan Spoke Guards with TiMotion Logo	Not with Mag rear wheels.	\$
TBCK3	Permobil Medical Necessities Backpack		\$
TPCH3	Permobil Medical Necessities Seat Pouch		\$
SP001	Bodypoint® Mobility Bag		\$
TUEST1	Upper Extremity Support Tray	Requires armrests.	\$
	12"-14" seat widths	13"-15" seat widths	
TLUG1	Luggage Carriers		\$
TMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe		\$
TCH1	Crutch Holder		\$

SmartDrive Power Assist System








MX2-850-xxx	Permobil SmartDrive MX2+	<i>Please provide SmartDrive MX2+ order form. Not available with rear seat heights < 13.5" or 20" rear wheels.</i>	\$
-------------	---------------------------------	---	----

General



Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams



Aero T and TRA
will always have
a rounded rear
frame.



TR rear frame
will be square
when the center
of gravity
measurement is
less than 3.5"

