



TiLite Z Series - Order form



Dealer Name*:

Account #*:

City, Post Code, Country:

Phone:

Email:

P.O. Number

Clinician:

Client Reference*:

Client Weight:

Rep Name:

**Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

Frame Style

AZFS1	Aero Z Aluminum Frame	Full adjustability. 265lb (120kg) weight limit. 150lb (68kg) weight limit for seat sizes 10"x10" to 12"x12".	\$
Z2FS2	ZRA Titanium Frame	Full adjustability. 265lb (120kg) weight limit. 150lb (68kg) weight limit for seat sizes 10"x10" to 12"x12".	\$
Z2FS1	ZR Titanium Frame	Minimal adjustability. 265lb (120kg) weight limit. 150lb (68kg) weight limit for seat sizes 10"x10" to 12"x12".	\$

Frame Upgrades

ZHDR1	Heavy Duty Frame	Not with Aero Z. 300lb (136kg) weight limit. Includes tension adjustable by straps back upholstery. Requires folding aluminum adjustable height backrest, billet aluminum front wheels, Shadow or Spinergy SPOX rear wheels. LX rear wheels only with 0° or 2° camber. Not with Slipstream fork, velcro adjustable upholstery, swing away tubular armrests or swing away flip back tubular armrests.	\$
ZPWR1	Power Adaptable Frame	Requires LiteSpeed billet aluminum front wheel and 0° camber. Not with single-sided forks or titanium camber tubes. Not required for SmartDrive.	\$
ZRFR1	Reinforced Frame		\$

Transit Options

ZTTD03	Transit Tie-Downs	Not with ZR. Tested to WC-19. Consists of four tie-down brackets, two pelvic belt mounting brackets and reinforced frame. Please refer to the Aero Z and ZRA Transit Tie-Down Option Brochure for important information regarding transit tie-down options.	\$
ZWTB	Transport Brackets	Consists of four tie-down brackets. Unoccupied use only. Please refer to the ZR Transport Bracket Supplement for important information regarding transit options.	\$

Frame Finishes

ZFF1	Satin Finish	Available with titanium frames only.	No charge
ZSCBT1	Buffing Pads for Satin Frame		\$
ZFF3	Tattoos	Tattoos on front frame tubes. Titanium frame finished in satin. Tattoos will be trimmed and adjusted to fit your frame.	\$
		ZFF45 Need for Speed ZFF44 Flower Garden ZFF43 Mosaic	
ZFF	Powder Coated Frame	Select below.	Standard

Frame Color:

ZFF8 High Gloss Black	ZFF10 Red Pepper	ZFF29 Neon Pink	ZFF14 Ocean Blue Metallic
ZFF15 Tangerine Metallic	ZFF30 Acid Green Metallic	ZFF17 Midnight Blue Pearl	ZFF18 White Pearl
ZFF19 Black Lava	ZFF27 Titanium Grey	ZFF59 Pomegranate Red	ZFF62 Blue Mist
ZFF57 Indigo Purple	ZFF40 Hyper Blue	ZFF41 Flame Orange	ZFF42 OD Green
ZFF60 Silver Streak	ZFF63 Lavendar Dusk	ZFF58 Black Galaxy	ZFF61 Seafoam Green
ZFF64 Daffodil Yellow			



Need for Speed



Flower Garden



Mosaic

Color Anodize Package

Standard:	Black						Standard
Anodized Package:	Gold	Green	Purple	Orange	Blue	Red	\$

Includes caster mount, caster fork and caster caps. If no anodize colour is selected, default will be black.

Centre of Gravity 'H' | Measure from front of back post at the seat tube to the centre of the rear axle. If no centre of gravity is selected the chair will ship with the centre of gravity in the most rearward position. Ergo seat selection may limit available centre of gravity.

ZCOG1:	0"	0.25"	0.5"	0.75"	1"	1.25"	1.5"	1.75"	2"	2.25"	2.5"	2.75"
	3"	3.25"	3.5"	3.75"	4"	4.25"	4.5"	4.75"	5"			

Rear Wheel Spacing 'J' | Measure from outside of back post to inside of rear tire. ZRA 0.75" rear wheel spacing not available with 0° or 2° camber. Ergonomic seat selection may limit available centre of gravity.

ZRWS1	Fixed Rear Wheel Spacing	Standard
	0.75" 1" 1.25"	
ZCBR12	Adjustable Rear Wheel Spacing Camber Plugs	No charge
ZCBR9	Additional Spacers to Increase Rear Wheel Spacing	No charge
ZCBR11	Zero Play Aluminum Camber Plugs Not with power adaptable frame or more than 6° camber.	No charge

Camber 'K'

ZCBR1:	0° 2° 4° 6°	
	More than 0° camber not with power adaptable frame. 8° camber may be available through our Creative Mobility Division.	
ZCBR15	Aluminum Camber Tube	Standard
ZCBR3	Titanium Camber Tube Not with power adaptable frame.	\$

Footrest

Footrest Width | Select standard front end 'I' if no seat taper 'T' has been selected. Select V-front end 'V' if a seat taper 'T' has been selected.

ZFE1	No V-Front Taper 'I' Find your rear seat width 'A' in the chart below and select one of the options to the right. Measure from inside of front frame tube to inside of opposite front frame tube.	No charge
------	--	-----------

Rear Seat Width 'A'	10"	7.5"		
	11"	8.5"		
	12"	9.5"	8"	
	13"	10.5"	9"	
	14"	11.5"	10"	
	15"	12.5"	11"	9"
	16"	13.5"	12"	10"
	17"	14.5"	13"	11"
	18"	15.5"	14"	12"
	19"	16.5"	15"	13"
	20"	18"	16"	14"
		No Taper	1.5" Taper	3.5" Taper
Footrest Width 'I'				

ZFE2	V-Front Taper 'V' Measure from inside of front frame tube to inside of opposite front frame tube 2.5" above footrest. Titanium footrests have 1" of adjustment up and down. Not with ergonomic seat frames.	\$
	9" 10" 11" 12" 13" 14" 15"	

Standard Footrests

Seat to Footrest 'E':	13.5"	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"
-----------------------	-------	-----	-------	-----	-------	-----	-------	-----	-------	-----	-------	-----

Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height 'C'.

ZFTR3	Angle Adjustable Footrest	Standard
ZFTR1	Titanium Footrest w/Flat ABS Cover	\$
ZFTR2	Titanium Open Loop Footrest	\$
ZFTR12	Angle Adjustable Flip-Back Footrest Not with 10" or 11" rear seat width, heavy duty frame or ABS cover with 2" sides.	\$

High Mount Footrests

Seat to Footrest 'E':	7"-9"	9"-11"	11"-13"	Set at
-----------------------	-------	--------	---------	--------

Measure from front edge of seat sling to top rear of footrest.

ZFTR4	Angle Adjustable High Mount Footrest	\$
-------	---	----

Residual Limb Support

ZRSL1	Residual Limb Support, 8"W x 10"L, Single	\$
	Left Right	
ZRSL2	Residual Limb Support, 8"W x 10"L, Pair	\$

Calf Straps

ZCLF2	Velcro Adjustable				\$
	Small (9"-12" footrest widths)	Medium (11"-15" footrest widths)	Large (14"-17" footrest widths)		
	Xtra Large (16"-19" footrest widths)				
SP102x	Bodypoint® Padded Velcro Adjustable				\$
	Small (8"-11" footrest widths)	Medium (11"-15" footrest widths)	Large (16"-19" footrest widths)		
ZCLF1	Buckle Adjustable, Black Nylon Webbing				\$

Backrest

Back Type

ZSB1	Folding Aluminum Adjustable Height Back	<i>Required for heavy duty frame.</i>	Standard
ZSB2	Folding Titanium Adjustable Height Back, Rounded Rigidizer Bar	<i>Titanium back will match frame finish.</i>	\$
ZSB15	Non-Folding Fixed Height Titanium Back, Rounded Rigidizer Bar	<i>Not with Aero Z, ZRA, transport brackets or push handles. Titanium back will match frame finish.</i>	\$

Back Options

ZSB11	Backrest Locks When Folded	<i>Available on all backs. Not with ergonomic seat.</i>	Standard
ZSB12	4" Deep Backrest Rigidizer Bar	<i>Not a rounded style rigidizer bar. Required for heavy duty frame. Not with non-folding, fixed height titanium back type.</i>	No charge
ZSB16	Depth Adjustable Back Frame	<i>Seat depths of 10"-18" allow 2" of rearward adjustment. Height adjustable push handles have a limited height range. Not with 19" or 20" seat depths, heavy duty frame, tension adjustable by straps seat upholstery or carbon fiber side guards.</i>	No charge
ZSB3	TiShaft™ Back Release Bar, Aluminum	<i>Not with non-folding fixed height titanium backs.</i>	Standard
ZSB5	TiShaft™ Back Release Bar, 2x2 Twill Weave Carbon Fiber	<i>Not with non-folding fixed height titanium backs.</i>	\$

Push Handles

ZSB6	Integrated Push Handles		No charge
ZSB7	Integrated Fold-Down Push Handles	<i>Requires tension adjustable by straps back upholstery.</i>	\$
ZSB8	Bolt-On Push Handles	<i>Includes 4" deep backrest rigidizer bar. Not with titanium back.</i>	\$
ZSB9	Bolt-On Height-Adjustable Push Handles	<i>Includes 4" deep backrest rigidizer bar.</i>	\$

Folding Back Height 'F' | ZSB1/ZSB2. Set at 0.5" increments.

Push Handles	None / Bolt-On	8.5"-11.5"	12"-16.5"	16"-20.5"	Set at:
	Integrated	10.5"-12.5"	13.5"-16.5"	17.5"-20.5"	Set at:
	Fold-Down	13"-15.5"	14"-18.5"	18.5"-21"	Set at:
		Short	Medium	Tall	Specify
Folding Back Height					

Fixed Back Height 'F'

Select height:	8"	8.5"	9"	9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"
	14"	14.5"	15"	15.5"	16"	16.5"	17"	17.5"	18"	18.5"	19"	19.5"
	20"											

Required for fixed side guards. Not with Aero Z, ZRA, swing-away tubular or swing-away flip-back tubular armrests, integrated push handles, integrated fold-down push handles or bolt-on height adjustable push handles.

Back Angle 'G' | ZSBA1 Measure angle from front of back post to floor.

Folding Back:	80°	83°	86°	89°	92°	95°	98°	101°					No charge
Fixed Back:	80°	81°	82°	83°	84°	85°	86°	87°	88°	89°	90°	91°	No charge
	92°	93°	94°	95°	96°	97°	98°	99°	100°	101°			

Back Supports

Back Upholstery

ZBUP2	Velcro Adjustable	<i>Not with fold-down push handles.</i>	Standard
ZBUP3	Tension Adjustable By Straps, Padded Nylon	<i>Required for fold-down push handles.</i>	\$
ZBUP5	Tension Adjustable By Straps, AIR Breathable	<i>Required for fold-down push handles.</i>	\$
ZUPC3	Ultrasuede Backrest Stripe, Black	<i>Not with AIR breathable upholstery.</i>	\$
ZBUP0	Omit Back Upholstery		No charge

ROHO® Agility® Active Aluminum Contour Back Support

AG2-xxF-2A	Agility Active Aluminum, Min Contour with Foam	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$
AG3-xxF-2A	Agility Active Aluminum, Mid Contour with Foam	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$

Shell Height: 8" 10"

ROHO® Agility® Active Carbon Contour Back Support | 2x2 Twill Weave Carbon Fiber

AG2-xxA-2CB	Agility Carbon Min Contour Back with Foam/Air Insert	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$
AG3-xxA-2CB	Agility Carbon Mid Contour Back with Foam/Air Insert	Available on 14" - 18" seat widths. Angle adjustable fixed hardware included.	\$

Shell Height: 10" 13" 16"

Comfort Acta-Embrace® Back Support

AE-BSxWxL	Acta-Embrace Back with 2-Point Quick Release Compass Hardware		\$
-----------	---	--	----

Shell Height: 10" 12" 14" 16"

Accessories: Privacy Shield \$

Comfort Acta-Relief Back Support

AR-BSxWxL	Acta-Relief Back with 2-Point Quick Release Compass Hardware		\$
-----------	--	--	----

Shell Height: 10" 12" 14" 16"

Accessories: Privacy Shield \$

Seat Options

Seat Upholstery

ZSUP1	Tension Adjustable Bolt-On Upholstery, Black	Not with ergonomic seat or tapered seat.	Standard
ZSUP2	Tension Adjustable By Straps, Black	Standard with ergonomic seat and seat taper of 1"-3".	\$
ZSUP8	Aluminum Solid Seat	Ergonomic seat requires CMD #16. Not with 10" or 11" rear seat width.	\$
ZSUP0	Omit Seat Upholstery		No charge
CMD018	Custom iSleeve	Please provide sleeve width and depth (inches) in Notes section at the end of this order form. Add 1" to device depth and 0.25" to width. Some limitations may apply.	\$

Seat Cushion

CU-SV-xx	Comfort Curve Cushion		\$
M2-FQA-xx	Comfort M2 Anti-Thrust Cushion		\$
EM-FQA-xx	Comfort Acta-Embrace™ Comfort-Tek Zero Elevation Cushion		\$
ALxxHD	ROHO® AirLite Cushion		\$
1Rxx	ROHO® Single Compartment Cushion	Fits seat widths 14"-19". Low Profile (LPC) High Profile(H-CA-SR)	\$
QSxx	ROHO® QUADTRO SELECT® Cushion	Fits seat widths 14"-19". Low Profile (LPC) High Profile (C)	\$
1RHExxC-SR	ROHO® Hybrid Elite Single Valve Cushion	Fits seat widths 14"-18", and 20"	\$

Precision Lock Front Forks

ZFK1	TiLite Standard Fork		Standard
ZFK2	TiLite Slipstream Single-Sided Fork	Not with heavy duty or power adaptable frame.	\$
ZGL1	Glide™ Suspension Fork	Black anodize only. Not with heavy duty frame. Select rider weight. 125lbs 125-180lbs 180-265lbs	\$
ZTFS1	Titanium Fork Stem		\$

Front Wheels

ZFW24	3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$
	Silver Hub Black Hub	
ZFW14	4"x0.75" Plastic Wheel w/Poly Tire <i>Not with Slipstream fork or power adaptable frame.</i>	No charge
ZFW3	4"x0.75" Light-Up Micro Caster <i>Not with Slipstream fork or power adaptable frame.</i>	\$
ZFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with power assist, power add-on system or heavy duty frame.</i>	\$
ZFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire	\$
	Silver Hub Black Hub	
ZFW6	5"x1" Plastic Wheel w/Poly Tire <i>Not with 90° front angle, Slipstream fork or power adaptable frame.</i>	No charge
ZFW28	5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with 90° front angle or power adaptable frame.</i>	\$
ZFW30	5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not with 90° front angle.</i>	\$
	Silver Hub Black Hub	
ZFW9	6"x1" Plastic Wheel w/Poly Tire <i>Not with 85° or 90° front angle, Slipstream fork or power adaptable frame.</i>	No charge
ZFW32	6"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with 85° or 90° front angle, Slipstream fork or power adaptable frame.</i>	\$
ZFW33	6"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not with 85° or 90° front angle.</i>	\$
	Silver Hub Black Hub	

Rear Wheels

See rear seat height ranges chart for rear seat height 'D' limitations.

ZRW14	TiLite SHADOW™ <i>25" and 26" wheels require 15" minimum seat depth.</i>	Standard
Wheel Size:	20" (451) 22" (501) 24" (540) 25" (559) 26" (590)	
ZRW7	Spinergy SPOX <i>25" and 26" wheels require 15" minimum seat depth. Orange, pink, purple and green spokes only available in 20" size.</i>	\$
Wheel Size:	20" (451) 22" (501) 24" (540) 25" (559) 26" (590)	
Spoke Color:	Black Red Blue White Orange Pink Purple Green	
ZRW8	Spinergy LX <i>25" and 26" wheels require 15" minimum seat depth.</i>	\$
Wheel Size:	22" (501) 24" (540) 25" (559) 26" (590)	
Spoke Color:	Black Red Blue White Orange Pink Purple Green	
ZRW9	Spinergy Design Lab LX Wheels <i>Black spokes only. If same color is selected, anodized hubs and rims match colour anodized package. 25" and 26" wheels require 15" minimum seat depth.</i>	\$
Wheel Size:	22" (501) 24" (540) 25" (559) 26" (590)	
Anodize Color:	Gold Green Blue Red Purple Orange	
ZRW16	Spinergy Carbon Blade <i>2x2 Twill Weave Carbon Fiber. Not with heavy duty frame. 25" wheels require 15" minimum seat depth.</i>	\$
Wheel Size:	24" (540) 25" (559)	
ZRW17	X-Core Mag <i>Not with heavy duty frame, short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i>	\$
Wheel Size:	20" (451) 22" (501) 24" (540)	

Axles

ZAXL1	Stainless Quick Release	Standard
ZAXL2	Titanium Quick Release	\$
ZAXL3	Quad Release <i>Not with Mag or 26" (590) Spinergy Spox or LX rear wheels.</i>	\$

Pneumatic Rear Tires

ZRTR1	Treaded 22"(501), 24"(540), 26"(590)	Standard
ZRTR2	Airless Inserts <i>Not with Spinergy Carbon Blade rear wheels.</i>	\$
ZRTR31	MBL SpeedLite 20"(451), 22"(501), 24"(540), 25"(559), 26"(590)	No charge
ZRTR30	Primo Sentinel 24"(540), 25"(559), 26"(590)	\$
ZRTR8	Schwalbe® Marathon Plus Evolution 24"(540), 25"(559), 26"(590)	\$
ZRTR32	Primo C-764 Knobby 24"(540)	\$

Solid Rear Tires

ZRTR22	Primo Orion 20"(451), 22"(501), 24"(540), 25"(559) <i>Not with Spinergy Carbon Blade rear wheels.</i>	\$
ZRTR23	Primo Active 24"(540), 25"(559) <i>Not with Spinergy Carbon Blade rear wheels.</i>	\$

Handrims

ZRIM1	Aluminum, Silver Anodized	Standard
ZRIM21	Aluminum, Black Anodized	No charge
ZRIM12	Titanium, Satin Finish <i>Not with 20"(451) or 22"(501) rear wheels or 24"(540) knobby tires.</i>	\$
ZRIM20	Q-Grip <i>Exclusive high-friction, anti-slip coating. Not with 20"(451) rear wheels.</i>	\$
ZRIM5	Plastic Coated, Black	\$
ZRIM10	The Surge <i>Large profile. All-In-One oval with Gription™ strip.</i>	\$
ZRIM11	The Surge LT <i>Small profile. All-In-One oval with Gription™ strip.</i>	\$
ZRIM14	The Natural-Fit®, Standard Grip <i>Original, large profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM15	The Natural-Fit®, Super Grip <i>Original, large profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM16	The Natural-Fit®, Without Thumb Grip <i>Original, large profile.</i>	\$
ZRIM17	The Natural-Fit® LT, Standard Grip <i>Reduced weight, small profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM18	The Natural-Fit® LT, Super Grip <i>Reduced weight, small profile. Not with short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
ZRIM19	The Natural-Fit® LT, Without Thumb Grip <i>Reduced weight, small profile.</i>	\$
ZRIM22	Omit Handrims	No charge

Handrim Tab Length

ZRIM7	Long Tabs	Standard
ZRIM8	Short Tabs <i>Not with Mag rear wheels or The Natural-Fit® handrims with thumb grip.</i>	\$

Wheel Locks

ZWLK7	Compact Composite Push to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
ZWLK8	Compact Composite Pull to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i>	No charge
ZWLK10	Extension Handles for Compact Composite Push/Pull to Lock	\$
ZWLK11	MBL Composite Scissor Lock	\$
ZWLK3	Uni-Lock	\$
	Push Pull	
ZWLK4	Extension Handles for Uni-Lock	\$
ZWLK15	Grade Aids Push to Lock <i>Not with solid rear tires.</i>	\$

Side Guards

ZSDG2	Aluminum Rigid Removable	\$
	Small Standard	
ZSDG3	Aluminum Quad Removable, with Hole	\$
	Small Standard	
ZSDG4	Aluminum Fold Down Rigid <i>Not with ergonomic seats.</i>	\$
	Small Standard	
ZSDG7	Carbon Rigid Removable	\$
	Small Standard	
ZSDG8	Carbon Quad Removable, with Hole	\$
	Small Standard	
ZSDG9	Carbon Fold Down Rigid <i>Not with ergonomic seats.</i>	\$
	Small Standard	
ZSDG6	Aluminum Fendered Adjustable <i>Available with 24"-25" rear wheels. Requires rear seat height 14.5"-18.5"</i>	\$
ZSDG12	Carbon Fendered Removable	\$
ZSDG13	Side Guards Custom Cut to Profile of Rear Wheels	\$
CMD007	Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 0.5" and side guards type (rigid or fold-down) in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i>	\$

Armrests

ZARM1	Swing Away Tubular	<i>Not with non-folding, fixed height titanium seat back or height adjustable push handles.</i>	\$
Range:	6.5"-9.5"	9.5"-12.5"	
ZARM6	Swing Away Fixed Height Flip Back Tubular	<i>Not with heavy duty frame.</i>	No charge
ZARM4	Swing Away Flip Back Tubular 8.5"-11.5"	<i>Not with heavy duty frame.</i>	\$
ZARM2	Removable Height Adjustable Desk Arm with Rigid Side Guard	<i>Requires 14" minimum seat depth. Not with ergonomic seat.</i>	\$
Range:	9"-12"	12"-15"	
Pad Size:	Desk Length (10")	Full Length (14")	
ZARM3	Removable Height Adjustable Desk Arm with Rigid Side Guard and Transfer Loop	<i>Requires 15" minimum seat depth. Not with ergonomic seat.</i>	\$
Range:	9"-12"	12"-15"	
Pad Size:	Desk Length (10")	Full Length (14")	
ZARM5	Armrests Covered in Black Ultrasuede	<i>Not with desk arms.</i>	\$

Positioning Belts

ZBLT1	Velcro Adjustable, Black Nylon Webbing		\$
HB205-xxx-B1	Bodypoint® Two Point Push-Button Buckle (Auto Style)		\$
	Medium	Large	
EB205-xxx	Bodypoint® EvoFlex® Pelvic Stabilizer, Push Button		\$
	Medium	Large	
HW320-25-2	Bodypoint® Optional Band Clamp Mounting Hardware		\$
ZBLT8	Q'Straint® WC-19 Wheelchair Anchored Pelvic Belt Recommended for use with transit tie down. Not with ZR.		\$

Anti-Tips

ZTIP1	Rear Aluminum, Black	\$
ZTIP2	User Friendly Flip-Up	\$

Accessories

ZGRD1	Clear Lexan Spoke Guards with TiMotion Logo	Not with Mag rear wheels.	\$
ZBCK3	Permobil Medical Necessities Backpack		\$
ZPCH3	Permobil Medical Necessities Seat Pouch		\$
SP001	Bodypoint® Mobility Bag		\$
ZUEST1	Upper Extremity Support Tray	Requires armrests.	\$
	12"-14" seat widths	13"-15" seat widths	
ZLUG1	Luggage Carriers		\$
ZMPCT2	Black Neoprene Impact Guards with Ultrasuede Stripe		\$
ZCH1	Crutch Holder		\$

SmartDrive Power Assist System

MX2-850-xxx	Permobil SmartDrive MX2+	<i>Please provide SmartDrive MX2+ order form. Not available with rear seat heights < 13.5" or 20" rear wheels.</i>	\$
-------------	---------------------------------	---	----

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams

