


Permobil Power Wheelchairs Battery Compliance Certificate for Transport

Manufacturer:	Permobil AB
Device type:	Powered wheelchair (non-foldable) WCBD Wheelchair - Battery Dry
Batteries type:	Sealed Non-Spillable , Valve Regulated Lead Acid (VRLA) – Gel/AGM
UN Classification:	UN2800 – Non-spillable Battery



Validity of document

This document is valid only for the Permobil power wheelchairs equipped with original Permobil batteries, marked with Permobil logo:

Part number	Type	Description	Ah	Voltage	Wh
1836054	Non spillable	Permobil VRLA Lead Acid (Gel/AGM)	85	12	1020
1836529	Non spillable	Permobil VRLA Lead Acid (Gel/AGM)	65	12	780

Permobil Battery Transport Compliance *

We hereby confirm that our batteries are sealed non-spillable **VRLA (UN2800)**, compliant with IATA DGR Special Provision A67, **not** lithium-based, therefore **exempt** from **Dangerous Goods Regulations** and the Watt-hour regulations applicable to lithium-based batteries.

Regulation / Organization	Permobil VRLA (Non-Spillable) Battery Compliance
D.O.T. (U.S. Department of Transportation)	✓ Approved — Permobil batteries conform to 49 CFR 173.159 and are classified as <i>non-spillable</i>
IMO / IMDG (International Maritime Organization / International Maritime Dangerous Goods)	✓ Approved — Permobil batteries meet SP 238 (Special Provision 238 : Must pass vibration & pressure differential tests to be exempt) and are <i>not restricted</i> for vessel/maritime transport
ICAO / IATA (International Civil Aviation Organization / International Air Transport Association)	✓ Approved — Permobil batteries fully comply with A67 (Special Provision A67 : Must pass vibration & pressure differential tests to be exempt from Dangerous Goods Regulations) and are <i>not restricted</i> for air transport
UN Classification	✓ UN2800 — Non spillable battery category

*Marked directly on Permobil original batteries | Valid for transport by Air, Water, Road or Rail

Standard transport handling process for WCBD (Wheelchair – Battery Dry) device

- No DG (Dangerous Goods) declaration required – **WCBD** Device is not treated as dangerous goods
- No removal of the batteries during transport is required
- Visual inspection shows no obvious defect or damage on batteries
- Power isolation / accidental activation prevention is required – Permobil wheelchair needs to be switched off with **main circuit breaker** (refer to the User manual for location of circuit breaker on your device)
- Removed batteries or spare batteries need to be placed in a strong rigid packaging with the packaging marked with “Non-Spillable Battery”.

On behalf of Permobil AB,
Date, Place: Sundsvall, 19/05/2026

Niclas Solberger, Product Owner PWC Evergreen