



TiLite Pilot - Order form



Dealer Name*:

Account #*:

City, State, Post Code:

Phone:

Email:

Rep Name:

Clinician:

Client Name | PO Number

Client Weight:

**Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

Frame Style

TPFS1	TiLite Pilot INCLUDING Growth Kit Full adjustability. 120lb (54kg) weight limit. Seat width 'A' grows 2" and seat depth 'B' grows 3". Price includes one-time growth kit and Transit Tie-Down option. Additional kits sold separately.	\$
TPFS2	TiLite Pilot WITHOUT Growth Kit Full adjustability. 120lb (54kg) weight limit. Seat width 'A' grows 2" and seat depth 'B' grows 3".	\$

Transit Options

TPWTB	Transport Brackets Consists of four tie-down brackets.	No charge
TPTDO	Transit Tie-Downs Tested to WC-19. Refer to the Transit Tie-Down Option Brochure for important information regarding transit tie-down options. WC-19 tie-down brackets require the TiLite back and head support. 95lb (43kg) weight limit. Not available for user weight less than 51lbs (23kg).	\$

Frame Finishes

Frame Color:	TPFF8 High Gloss Black	TPFF10 Red Pepper	TPFF29 Neon Pink
	TPFF14 Ocean Blue Metallic	TPFF15 Tangerine Metallic	TPFF30 Acid Green Metallic
	TPFF17 Midnight Blue Pearl	TPFF18 White Pearl	TPFF19 Black Lava
	TPFF27 Titanium Grey	TPFF59 Pomegranate Red	TPFF62 Blue Mist
	TPFF57 Indigo Purple	TPFF40 Hyper Blue	TPFF41 Flame Orange
	TPFF42 OD Green	TPFF60 Silver Streak	TPFF63 Lavendar Dusk
	TPFF58 Black Galaxy	TPFF61 Seafoam Green	TPFF64 Daffodil Yellow

Color Anodize Package

Standard:	ZAN01 Black						Standard
Anodized Package:	ZAN12 Gold	ZAN15 Green	ZAN13 Purple	ZAN14 Orange	ZAN10 Blue	ZAN11 Red	No charge

Includes caster mount, caster fork and caster caps. If no anodize colour is selected, default will be black.

Frame Measurements

Please refer to Reference Diagrams at end of form for more information.

Seat Width 'A' | Measure inside width from side panel to side panel. Seat width grows 2" if selected width is between 8" and 12". 13" width will grow a maximum of 1", 14" will not grow any wider.

TPSW1:	8"	9"	10"	11"	12"	13"	14"
--------	----	----	-----	-----	-----	-----	-----

Seat Depth 'B' | The seat depth can be grown within the selected range.

Seat Depth Range:	9"-12"	11"-15"					
TPSD1:	9"	10"	11"	12"	13"	14"	15"

Front Seat Height 'C' | Measure from the floor to the top edge of the front of the seat pan. Adjustable in 0.5" increments.

TPFSH1:	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"
---------	-----	-------	-----	-------	-----	-------	-----	-------	-----	-------

Rear Seat Height 'D' | Measure from the floor to the top edge of the seat pan at the intersection of the front of the back support pad and seat pan. Adjustable in 0.5" increments.

Seat Depth 'B'	TPRSH1:	12"	12.5"	13"	13.5"	14"	14.5"	15"	15.5"	16"	16.5"
		20"			22"			24"			
	9"	12.5"-15"	12"-13.5"	0.5"-1.5"							
	10"	12"-15"	12"-13.5"	0.5"-2.25"							
	11"	12"-15"	12"-13.5"	0.5"-3"	13.5"-16.5"	13"-15"	0.5"-2.25"	14"-16.5"	13.5"-16"	0.5"-2.25"	
	12"	12"-15"	12"-13.5"	0.5"-3.5"	13.5"-16"	13"-15"	0.5"-3.25"	14"-16.5"	13.5"-16"	0.5"-3.25"	
	13"				13"-16"	13"-15"	0.5"-3.5"	13.5"-16.5"	13.5"-16"	0.5"-3.5"	
	14"				13"-16"	13"-15"	0.5"-3.5"	13.5"-16.5"	13.5"-16"	0.5"-3.5"	
	15"				13"-16"	13"-15"	1.5"-3.5"	13.5"-16.5"	13.5"-15.5"	1.5"-3.5"	
		Front Height	Rear Height	COG	Front Height	Rear Height	COG	Front Height	Rear Height	COG	
Measurements											

Rear Wheel Spacing | Measure from outside of side panel to inside of rear tire at top of wheel. Not all rear wheel spacing is available with all accessories.

TPRWS1	Fixed Rear Wheel Spacing										Standard
		0.75"	1"	1.25"	1.5"	1.75"	2"				
TPCBBR12	Adjustable Rear Wheel Spacing Camber Plugs										Standard
TPCBBR9	Additional Spacers to Increase Rear Wheel Spacing										No charge

Camber | Black aluminum camber tube standard.

Camber Angle:	2° (Minimum 1.5" rear wheel spacing)	4° (Minimum 1.25" rear wheel spacing)
	6° (Minimum 1" rear wheel spacing)	8°
	12°	

Footrest

Seat to Footrest 'E':	3"	3.5"	4"	4.5"	5"	5.5"	6"	6.5"	7"	7.5"	8"	8.5"
	9"	9.5"	10"	10.5"	11"	11.5"	12"	12.5"	13"	13.5"	14"	14.5"
Measure from the top of the footplate to the front of the seat pan. For angle adjustable footrest, measurement must be at least 2" shorter than front seat height. For flip-back footrest, measurement must be at least 3" shorter than front seat height.												
TPFTR3	Angle Adjustable Footrest											Standard
TPFTR24	Flip-Back Depth and Angle Adjustable Footrest Not available on 8", 9" or 14" seat widths.											\$

Backrest

Back Type

TPSB2	Contoured Non-Folding Adjustable Height Rigid Back with Pad 11"-13" height range.	No charge
TPSB3	Aluminum Canes for Aftermarket Back Supports Not available with height adjustable fenders.	No charge

Push Handles | Height adjustable range is 33"-40" measured from the floor with 24" wheels.

TPSB6	Bolt-On Height Adjustable Removable Push Handle, Single	\$
TPSB7	Bolt-On Height Adjustable Removable Push Handle, Dual	\$

Back Supports

Comfort Acta-Embrace® Back Support

AE-BSxWxL	Acta-Embrace® Back Support Available on 10"-14" seat widths. Includes compass can mount hardware.	\$
-----------	---	----

Shell Height:	10"	12"	14"
----------------------	-----	-----	-----

Comfort Back Options

TPLAT1	BodiLink® 3"L x 4"D Lateral Trunk Support, Left Swing away, GlideWear cover, non-padded.	\$
TPLAT1	BodiLink® 3"L x 4"D Lateral Trunk Support, Right Swing away, GlideWear cover, non-padded.	\$
TPHDR1	Comfort BodiLink Head Support w/GlideWear® Cover Uses piano taper midline hardware. Available on Acta-Back, Acta-Embrace and TiLite rigid back supports. The Comfort Head Support is WC-20 certified on the Acta-Back and Acta-Embrace, WC-19 certified when used with the TiLite backrest.	\$
	6" Width 8" Width	

Seat Options

TPSP1	Solid Seat Pan <i>Required.</i>	\$
Seat Cushion <i>Refer to reference charts at end of form for more information.</i>		
CU-SV-xx	Comfort Curve Cushion <i>Stretch-Air cover.</i>	\$
EM-FQA-xx	Comfort Embrace Comfort-Tek Zero Elevation Cushion	\$
1RxxLPC	ROHO® Low Profile® Single Compartment Cushion	\$

Front Wheels

TPFW14	4"x0.75" Plastic Wheel w/Poly Tire	Standard
TPFW3	4"x0.75" Light-Up Micro Caster	\$
TPFW26	4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire	\$
TPFW27	4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire, Black Hub	\$

Rear Wheels

20" rear wheels only available with seat depth range 9"-12". 22" and 24" rear wheels only available with 11"-15" seat depth range.

TPRW14	TiLite SHADOW™	No charge
Wheel Size:	20" (451) 22" (501) 24" (540)	
TPRW7	Spinergy SPOX <i>Orange, pink, purple and green spokes only available in 20" size.</i>	\$
Wheel Size:	20" (451) 22" (501) 24" (540)	
Spoke Color:	Black Red Blue White Orange Pink Purple Green	
TPRW8	Spinergy LX	\$
Wheel Size:	22" (501) 24" (540)	
Spoke Color:	Black Red Blue White Orange Pink Purple Green	
TPRW9	Spinergy Design Lab LX Wheels <i>Black spokes only. If same color is selected, anodized hubs and rims match colour anodized package.</i>	\$
Wheel Size:	22" (501) 24" (540)	
Anodize Color:	Gold Green Blue Red Purple Orange	
TPRW16	Spinergy Carbon Blade <i>2x2 Twill Weave Carbon Fiber.</i>	\$
Wheel Size:	24" (540)	
TPRW17	X-Core Mag <i>Not with short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i>	\$
Wheel Size:	20" (451) 22" (501) 24" (540)	

Axles

TPAXL1	Stainless Quick Release	Standard
TPAXL3	Quad Release <i>Not with Mag rear wheels.</i>	\$

Tires

TPRTR1	Treaded 20"(451), 22"(501), 24"(540)	Standard
TPRTR2	Airless Inserts	\$
TPRTR4	MBL SpeedLite 20"(451), 22"(501), 24"(540)	No charge
TPRTR30	Primo Sentinel 20"(451), 22"(501), 24"(540)	\$
TPRTR8	Schwalbe® Marathon Plus Evolution 24"(540)	\$

Handrims

TPRIM1	Aluminum, Silver Anodized	Standard
TPRIM21	Aluminum, Black Anodized <i>Not with 20" rear wheels.</i>	No charge
TPRIM20	Q-Grip <i>Exclusive high-friction, anti-slip coating. Not with 20" rear wheels.</i>	\$
TPRIM5	Plastic Coated, Black <i>Not with Mag rear wheels.</i>	\$
TPRIM11	The Surge LT <i>Small profile. All-In-One oval with Gription™ strip. Not with 20" rear wheels.</i>	\$
TPRIM17	The Natural-Fit® LT, Standard Grip <i>Reduced weight, small profile. Not with 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
TPRIM18	The Natural-Fit® LT, Super Grip <i>Reduced weight, small profile. Not with 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i>	\$
TPRIM19	The Natural-Fit® LT, Without Thumb Grip <i>Reduced weight, small profile. Not with 20" rear wheels.</i>	\$
TPRIM22	Omit Handrims	No charge

Handrim Tab Length

TPRIM7	Long Tabs	No charge
TPRIM8	Short Tabs <i>Not with Mag rear wheels or The Natural-Fit® handrims with thumb grip.</i>	\$

Wheel Locks

TPWLK30	Compact Push to Lock <i>Redesigned with more clamping force.</i>	No charge
TPWLK31	Compact Pull to Lock <i>Redesigned with more clamping force.</i>	No charge

Mounting Location:	Standard Location 'F'	Forward Mount Location 'G'	No charge
---------------------------	-----------------------	----------------------------	-----------

Side Guards

TPSDG30	Built In Side Guard with Height Adjustable Fender <i>Not with 20" rear wheels.</i>	\$
TPSDG31	Built In Side Guard without Fender	Standard

Armrests

TPARM2	Height Adjustable Desk Arm <i>10" armpad. Works with fendered side guard attachment option.</i>	\$
TPARM31	Height Adjustable Armrest <i>Attaches to built in side guards.</i>	Standard

Positioning

Upper Body Positioning

SH295x-B5	Bodypoint® PivotFit™ Dynamic Shoulder Harness <i>Cam buckle mounts, rear pull.</i>	\$
	Extra Small Small	
SH120x-B3	Bodypoint® Monoflex™ Chest Support	\$
	Extra Small Small	

Pelvic Positioning

TPBLT1	Velcro Adjustable Belt, Black Nylon Webbing	\$
TPBLT8	Q'Straint® WC-19 Wheelchair Anchored Pelvic Belt <i>Recommended for use with transit tie down.</i>	\$
HB275-S38-B1	Bodypoint® 1" Two-Point Padded Hip Belt <i>Swivel buckle.</i>	\$
EB275-S38	Bodypoint® 1" Evoflex® Belt <i>Small push-button buckle.</i>	\$

Lower Body Positioning

FT240x	Bodypoint® Ankle Huggers	\$
	Extra Small Small	

Anti-Tips

TPTIP1	Rear Aluminum, Black	\$
TPTIP7	Single Anti-Tip	\$
	Left Side Right Side	

Accessories

TPGRD30	Clear Polycarbonate Spoke Guards w/TiLite Pilot Sticker Pack <i>Not with Mag rear wheels.</i>	\$
TPBCK3	Permobil Medical Necessities Backpack	\$
TPCH1	Crutch Holder	\$

SmartDrive Power Assist System

MX2-850-xxx	Permobil SmartDrive MX2+ <i>Please provide SmartDrive MX2+ order form. Not available with 8" seat width, rear seat heights < 13.5" or 20" rear wheels.</i>	\$
-------------	--	----

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams

