



TiLite Twist - Order form



Dealer Name*:

Account #*:

City, State, Post Code:

Phone:

Email:

Rep Name:

Clinician:

Client Name | PO Number

Client Weight:

**Required field*

Special considerations that need to be addressed here (e.g., diagnosis):

Frame Style

| | | |
|--------|---|----|
| TWFS2 | TiLite Twist Full adjustability. 165lb (75kg) weight limit. Seat width 'A' grows 2" and seat depth 'B' grows 3". | \$ |
| TWFSW7 | Abduction Frame, 2" Outward Flare Requires a No Taper Footrest. | \$ |

Transit Options

| | | |
|--------|---|----|
| TWTTDO | Transit Tie-Downs Tested to WC-19. Consists of four tie-down brackets. Transit option is not available with user weights less than 51lb (23kg). Please refer to the Twist Transit Tie-Down Option Brochure for important information regarding transit tie-down options. | \$ |
|--------|---|----|

Frame Finishes

| | | | |
|-------------------|----------------------------|---------------------------|----------------------------|
| Frame Color: | TWFF8 High Gloss Black | TWFF10 Red Pepper | TWFF49 Neon Pink |
| | TWFF14 Ocean Blue Metallic | TWFF15 Tangerine Metallic | TWFF30 Acid Green Metallic |
| Caster Arm Color: | TWFF17 Midnight Blue Pearl | TWFF18 White Pearl | TWFF19 Black Lava |
| | TWFF27 Titanium Grey | TWFF59 Pomegranate Red | TWFF62 Blue Mist |
| | TWFF57 Indigo Purple | TWFF40 Hyper Blue | TWFF41 Flame Orange |
| | TWFF42 OD Green | TWFF60 Silver Streak | TWFF63 Lavendar Dusk |
| | TWFF58 Black Galaxy | TWFF61 Seafoam Green | TWFF64 Daffodil Yellow |
| | TWFF8 High Gloss Black | TWFF10 Red Pepper | TWFF49 Neon Pink |
| | TWFF14 Ocean Blue Metallic | TWFF15 Tangerine Metallic | TWFF30 Acid Green Metallic |
| | TWFF17 Midnight Blue Pearl | TWFF18 White Pearl | TWFF19 Black Lava |
| | TWFF27 Titanium Grey | TWFF45 Pomegranate Red | TWFF48 Blue Mist |
| | TWFF43 Indigo Purple | TWFF40 Hyper Blue | TWFF41 Flame Orange |
| | TWFF42 OD Green | TWFF46 Silver Streak | TWFF49 Lavendar Dusk |
| | TWFF44 Black Galaxy | TWFF47 Seafoam Green | TWFF50 Daffodil Yellow |

Color Anodize Package

| | | | | | | | |
|-------------------|-------------|-------------|--------------|--------------|------------|-----------|----------|
| Standard: | ZAN01 Black | | | | | | Standard |
| Anodized Package: | ZAN12 Gold | ZAN15 Green | ZAN13 Purple | ZAN14 Orange | ZAN10 Blue | ZAN11 Red | \$ |

Includes caster mount, caster fork and caster caps. If no anodize colour is selected, default will be black.

Frame Measurements

Seat Width 'A' | Measure from outside of seat tube at back post to the outside of the opposite seat tube at back post. Seat width grows 2". All upholstery and components (exception cushion) needed to grow the seat width are included. TiShaft back release bar available at no charge upon request.

| | | | | | | |
|--------|--------|--------|---------|---------|---------|---------|
| TWSW5: | 8"-10" | 9"-11" | 10"-12" | 11"-13" | 12"-14" | 13"-15" |
| TWSW5: | Set at | | | | | |

Default will be smallest seat width for selected range. Please specify for other settings (1" increments).

Seat Depth 'B' | Measure from front of back post to the front edge of seat sling. Seat depth grows 3". Front edge of seat sling will be 1.5" from the beginning of the bend on no-taper and tapered frames, and 2" on abduction frames. All upholstery and components needed to grow the seat depth are included (except cushion).

TWSD4: 8"-11" 9"-12" 10"-13" 11"-14" 12"-15" 13"-16"

TWSD4: Set at

8"-11" range requires 18" or 20" rear wheels. Default will be smallest seat depth for selected range. Please specify for other settings (1" increments).

Front Seat Height 'C' | Measure from floor to top of seat tube at front edge of seat sling. Twist adjusts up and down in 1/2" increments. 2" taper footrest width 'I' available with 12"-15" or 13"-16" seat depths only. For 2" taper footrest width 'I' add 1.5" to minimum front seat height.

TWFSH1: 13" 13.5" 14" 14.5" 15" 15.5" 16" 16.5" 17" 17.5" 18" No charge

| Front Wheel | 3" | 13"-14.5" | 13.5"-15" | 13"-16" | 14"-17" | 14.5"-17.5" | 15.5"-18.5" |
|-------------------------|----------|--------------|--------------|---------------|---------------|---------------|---------------|
| | 3" Glide | 14.5"-15.5" | 15"-16.5" | 14"-17" | 15"-18" | 16"-19" | 16.5"-19.5" |
| | 4" | 13.5"-15" | 14"-15.5" | 13.5"-16.5" | 14"-17.5" | 15"-18" | 16"-18.5" |
| | 4" Glide | 14.5"-16" | 15.5"-17" | 14.5"-17.5" | 15.5"-18.5" | 16"-19" | 17"-20" |
| | 5" | | | 14"-17" | 14.5"-17.5" | 15.5"-18.5" | 16.5"-19.5" |
| | | 8"-11" Depth | 9"-12" Depth | 10"-13" Depth | 11"-14" Depth | 12"-15" Depth | 13"-16" Depth |
| Front Seat Height Range | | | | | | | |

Rear Seat Height 'D' | Measure from floor to top of seat tube at back post. Twist adjusts up and down in 1/2" increments.

TWRSH1: 12" 12.5" 13" 13.5" 14" 14.5" 15" 15.5" 16" 16.5" 17" No charge

| Rear Wheel | 18" | 12"-17" |
|------------|-----|-----------|
| | 20" | 13"-17.5" |
| | 22" | 14"-18.5" |
| | 24" | 15"-19.5" |
| | | |

Rear Seat Height Range

Front Angle 'L' | Measure from frame front to floor.

TWFA1: 70° 80° 90°

No charge

Rear Wheel Spacing 'J' | Measure from outside of back post to inside of rear tire. 3/4" rear wheel spacing not available with 0° or 2° camber. Not all rear wheel spacing is available with all accessories.

| | | | | | | | |
|---------|---|----|-------|------|-------|----|-----------|
| TWRWS1 | Fixed Rear Wheel Spacing | | | | | | Standard |
| | 0.75" | 1" | 1.25" | 1.5" | 1.75" | 2" | |
| TWCBR12 | Adjustable Rear Wheel Spacing Camber Plugs | | | | | | No charge |
| TWCBR9 | Additional Spacers to Increase Rear Wheel Spacing | | | | | | No charge |

Camber 'K'

| | | | | | | | |
|---------|----------------------|--|--|--|--|--|----------|
| TWCBR15 | Aluminum Camber Tube | | | | | | Standard |
|---------|----------------------|--|--|--|--|--|----------|

TWCBR1: 0° 2° 4° 6° 8° 12°

Footrest

Footrest Width 'I' | Measure from inside of front frame tube to inside of opposite front frame tube.

| | | | |
|--------|----------|---|-----------|
| TWFSW1 | No Taper | 2" narrower than seat width. | Standard |
| TWFSW4 | 2" Taper | 4" narrower than seat width. Available with 12"-15" or 13"-16" seat depths only. Add 1.5" to minimum seat to footrest 'E' ranges. | No charge |

| Seat Width 'A' | 8" | 6" | |
|----------------|----------------|----------|----------|
| | 9" | 7" | |
| | 10" | 8" | |
| | 11" | 9" | |
| | 12" | 10" | 8" |
| | 13" | 11" | 9" |
| | 14" | 12" | 10" |
| | 15" | 13" | 11" |
| | | No Taper | 2" Taper |
| | Footrest Width | | |

Standard Footrests

| | | | | | | | | | | | | |
|-----------------------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|-----|-------|
| Seat to Footrest 'E': | 9" | 9.5" | 10" | 10.5" | 11" | 11.5" | 12" | 12.5" | 13" | 13.5" | 14" | 14.5" |
| | 15" | 15.5" | 16" | 16.5" | | | | | | | | |

Measure from front edge of seat sling to top rear of footrest. Measurement must be at least 2" shorter than Front Seat Height 'C'.

| | | |
|--------|---|--|
| TWFT3 | Angle Adjustable Footrest | Standard |
| TWFT24 | Flip-Back Depth and Angle Adjustable Footrest | Not when both 13"-15" seat width and abduction frame are ordered. \$ |

High Mount Footrests

| | | | | | | | | | | | | |
|-----------------------|-----|-------|-----|-------|-----|------|----|------|----|------|-----|-------|
| Seat to Footrest 'E': | 5" | 5.5" | 6" | 6.5" | 7" | 7.5" | 8" | 8.5" | 9" | 9.5" | 10" | 10.5" |
| | 11" | 11.5" | 12" | 12.5" | 13" | | | | | | | |

Measure from front edge of seat sling to top rear of footrest.

| | | |
|--------|--|--|
| TWFT4 | Angle Adjustable High Mount Footrest | \$ |
| TWFT27 | Flip-Back Depth and Angle Adjustable High Mount Footrest | Not when both 13"-15" seat width and abduction frame are ordered. \$ |

| Seat Depth 'B' | 3" Front Wheel | | 4" Front Wheel | | 5" Front Wheel | |
|-----------------------|----------------|-----------|----------------|-----------|----------------|-----------|
| | 8"-11" | 9"-12" | 5"-8" | 9"-12" | 5"-8" | |
| | 9"-12" | 10"-13" | 6"-9" | 10"-13" | 6"-9" | |
| | 10"-13" | 9"-13.5" | 5"-9.5" | 9"-13.5" | 5"-9.5" | 9"-14" |
| | 11"-14" | 10"-14.5" | 6"-10.5" | 10"-14.5" | 7"-11.5" | 10"-15" |
| | 12"-15" | 11"-15" | 7"-11.5" | 11"-15" | 8"-12.5" | 11"-15.5" |
| | 13"-16" | 12"-16" | 8"-12.5" | 12"-16" | 7.5"-12" | 12"-16.5" |
| | | Standard | High-Mount | Standard | High-Mount | Standard |
| Footrest Height Range | | | | | | |

Residual Limb Support

| | | |
|-------|---|----|
| TWRS1 | Residual Limb Support, 8"W x 10"L, Single | |
| | Right Left | |
| TWRS2 | Residual Limb Support, 8"W x 10"L, Pair | \$ |
| TWRS1 | Residual Limb Support, 8"W x 10"L, Single | \$ |
| | Right Left | |

Calf Straps

| | | |
|--------|---|----|
| TWCLF2 | Velcro Adjustable | \$ |
| | Small (9"-12" footrest widths) Medium (11"-15" footrest widths) | |
| SP102x | Bodypoint® Padded Velcro Adjustable | \$ |
| | Medium (11"-15" footrest widths) | |
| TWCLF1 | Buckle Adjustable, Black Nylon Webbing | \$ |

Backrest

Back Type

| | | |
|--------|--|--|
| TWSB1 | Folding Aluminum Adjustable Height Back | With TiShaft aluminum back release bar. Standard |
| TWSB11 | Backrest Locks When Folded | Standard |
| TWSB22 | Height Adjustable Backrest Rigidizer Bar | \$ |

Push Handles

| | | |
|--------|--|---|
| TWSB6 | Integrated Push Handles | No charge |
| TWSB7 | Integrated Fold-Down Push Handles | Requires tension adjustable by straps back upholstery. \$ |
| TWSB8 | Bolt-On Push Handles | Requires 12"-14" or 13"-15" seat width range. \$ |
| TWSB9 | Bolt-On Height-Adjustable Push Handles | Requires 12"-14" or 13"-15" seat width range. \$ |
| TWSB20 | Bolt-On Center Mount Height Adjustable Push Handle | \$ |
| | Small (9.25") Standard (14.25") Tall (19.25") | |
| TWSB21 | Stroller Handle | \$ |

Back Height 'F' | Set at 1/2" increments.

| | | | | | |
|---------------------|----------------|-------------|-------------|-------------|---------|
| Push Handles | None / Bolt-On | 8.5"-11.5" | 12"-16.5" | 16.5"-20.5" | Set at: |
| | Integrated | 10.5"-12.5" | 13.5"-16.5" | 17.5"-20.5" | Set at: |
| | Fold-Down | 13"-15.5" | 14.5"-19" | 18.5"-21" | Set at: |
| | | Short | Medium | Tall | Specify |
| Folding Back Height | | | | | |

Back Angle 'G' | Measure angle from front of back post to floor. 90° will be perpendicular to the floor. To approximate the angle from the back post to the seat, subtract 2° from measurement G for every 1/2" of seat dump.

TSWBA1: 80° 83° 86° 89° 92° 95° 98° 101° No charge

Back Upholstery

| | | | | | | |
|--------|---|--------------------------------------|------------|--------------|-------------|---------------------|
| TWBUP2 | Velcro Adjustable | Not with fold-down push handles. | | | | Standard |
| TWBUP3 | Tension Adjustable By Straps, Padded Nylon | Required for fold-down push handles. | | | | \$ |
| TWBUP5 | Tension Adjustable By Straps, AIR Breathable | Required for fold-down push handles. | | | | \$ |
| TWBUP5 | Ultrasuede Backrest Stripe, Black | Not with AIR breathable upholstery. | | | | \$ |
| | Black | Sapphire | Red Pepper | Burnt Orange | Ultraviolet | NiteLite Reflective |
| TWBUP0 | Omit Back Upholstery | | | | | No charge |

Back Supports**Comfort Acta-Embrace® Back Support**

AE-BSxWxL Acta-Embrace® Back Support Available on 10"-14" seat widths. Includes compass can mount hardware. \$

Shell Height: 10" 12" 14"

Comfort Back Options

| | | | |
|--------|---|--|----|
| TPLAT1 | BodiLink® 3"L x 4"D Lateral Trunk Support, Left | Swing away, GlideWear cover, non-padded. | \$ |
| TPLAT1 | BodiLink® 3"L x 4"D Lateral Trunk Support, Right | Swing away, GlideWear cover, non-padded. | \$ |
| TPHDR1 | Comfort BodiLink Head Support w/GlideWear® Cover | Uses piano taper midline hardware. Available on Acta-Back, Acta-Embrace and TiLite rigid back supports. The Comfort Head Support is WC-20 certified on the Acta-Back and Acta-Embrace, WC-19 certified when used with the TiLite backrest. | \$ |
| | 6" Width | 8" Width | |

Seat Options**Seat Upholstery**

| | | |
|--------|---|-----------|
| TWSUP1 | Tension Adjustable Bolt-On Upholstery, Black | Standard |
| TWSUP0 | Omit Seat Upholstery | No charge |

Seat Cushion

| | | |
|-----------|--|----|
| CU-SV-xx | Comfort Curve Cushion | \$ |
| EM-FQA-xx | Comfort Embrace Zero Elevation Cushion | \$ |
| 1RxxLPC | ROHO® Low Profile® Single Compartment Cushion | \$ |

Precision Lock Front Forks

| | | | | |
|-------|---|------------|------------|-----------|
| TWFK1 | TiLite Standard Fork | | | No charge |
| TWFK2 | TiLite Slipstream Single-Sided Fork | | | \$ |
| TWGL1 | Glide™ Suspension Fork <i>Black anodize only. Not with 5" front wheels. Select rider weight.</i> | | | \$ |
| | 125lbs | 125-180lbs | 180-265lbs | |
| TWTS1 | Titanium Fork Stem | | | \$ |

Front Wheels

See front seat height ranges chart for limitations.

| | | |
|--------|--|-----------|
| TWFW24 | 3"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire | \$ |
| | Silver Hub Black Hub | |
| TWFW14 | 4"x0.75" Plastic Wheel w/Poly Tire <i>Not with Slipstream fork.</i> | No charge |
| TWFW3 | 4"x0.75" Light-Up Micro Caster <i>Not with Slipstream fork.</i> | \$ |
| TWFW26 | 4"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire | \$ |
| TWFW27 | 4"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire | \$ |
| | Silver Hub Black Hub | |
| TWFW6 | 5"x1" Plastic Wheel w/Poly Tire <i>Not with 90° front angle or Slipstream fork.</i> | No charge |
| TWFW28 | 5"x1.4" LiteSpeed-7 Plastic Wheel w/Soft Roll Tire <i>Not with 90° front angle.</i> | \$ |
| TWFW30 | 5"x1.4" LiteSpeed-7 Billet Aluminum Wheel w/Soft Roll Tire <i>Not with 90° front angle.</i> | \$ |
| | Silver Hub Black Hub | |

Rear Wheels

See rear seat height ranges chart for limitations.

| | | |
|-----------------------|--|----------|
| TWRW14 | TiLite SHADOW™ | Standard |
| Wheel Size: | 18" (400) 20" (451) 22" (501) 24" (540) | |
| TWRW7 | Spinergy SPOX <i>Orange, pink, purple and green spokes only available in 20" size.</i> | \$ |
| Wheel Size: | 20" (451) 22" (501) 24" (540) | |
| Spoke Color: | Black Red Blue White Orange Pink Purple Green | |
| TWRW8 | Spinergy LX | \$ |
| Wheel Size: | 22" (501) 24" (540) | |
| Spoke Color: | Black Red Blue White Orange Pink Purple Green | |
| TWRW9 | Spinergy Design Lab LX Wheels <i>Black spokes only. If same color is selected, anodized hubs and rims match colour anodized package.</i> | \$ |
| Wheel Size: | 22" (501) 24" (540) | |
| Anodize Color: | Gold Green Blue Red Purple Orange | |
| TWRW16 | Spinergy Carbon Blade <i>2x2 Twill Weave Carbon Fiber.</i> | \$ |
| Wheel Size: | 24" (540) | |
| TWRW17 | X-Core Mag <i>Not with short tabs, quad release axles, spoke guards or omit handrims options. Requires aluminum silver anodized handrims.</i> | \$ |
| Wheel Size: | 20" (451) 22" (501) 24" (540) | |
| Axles | | |
| TWAXL1 | Stainless Quick Release | Standard |
| TWAXL2 | Titanium Quick Release | \$ |
| TWAXL3 | Quad Release <i>Not with Mag rear wheels.</i> | \$ |

Pneumatic Rear Tires

| | | |
|---------|---|-----------|
| TWRTR1 | Treaded | Standard |
| TWRTR2 | Airless Inserts <i>Not with 18" rear wheels.</i> | \$ |
| TWRTR4 | MBL SpeedLite 20"(451), 22"(501), 24"(540) | No charge |
| TWRTR30 | Primo Sentinel 24"(540) | \$ |
| TWRTR8 | Schwalbe® Marathon Plus Evolution 24"(540) | \$ |
| TWRTR32 | Primo C-764 Knobby 24"(540) | \$ |

Solid Rear Tires

| | | |
|---------|---|----|
| TWRTR22 | Primo Orion 18" (400), 20"(451), 22"(501), 24"(540) <i>Not with Spinergy Carbon Blade rear wheels.</i> | \$ |
| TWRTR23 | Primo Active 24"(540) <i>Not with Spinergy Carbon Blade rear wheels.</i> | \$ |

Handrims

| | | |
|---------|---|-----------|
| TWRIM1 | Aluminium, Silver Anodized | Standard |
| TWRIM21 | Aluminium, Black Anodized <i>Not with 18" rear wheels.</i> | No charge |
| TWRIM12 | Titanium, Satin Finish <i>Not with 18", 20" or 22" rear wheels, or 24" knobby tires.</i> | \$ |
| TWRIM20 | Q-Grip <i>Exclusive high-friction, anti-slip coating. Not with 18" or 20" rear wheels.</i> | \$ |
| TWRIM5 | Plastic Coated, Black <i>Not with 18" rear wheels.</i> | \$ |
| TWRIM10 | The Surge <i>Large profile. All-In-One oval with Gription™ strip. Not with 18" or 20" rear wheels.</i> | \$ |
| TWRIM11 | The Surge LT <i>Small profile. All-In-One oval with Gription™ strip. Not with 18" or 20" rear wheels.</i> | \$ |
| TWRIM14 | The Natural-Fit®, Standard Grip <i>Original, large profile. Not with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i> | \$ |
| TWRIM15 | The Natural-Fit®, Super Grip <i>Original, large profile. Not with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i> | \$ |
| TWRIM16 | The Natural-Fit®, Without Thumb Grip <i>Original, large profile.</i> | \$ |
| TWRIM17 | The Natural-Fit® LT, Standard Grip <i>Reduced weight, small profile. Not with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i> | \$ |
| TWRIM18 | The Natural-Fit® LT, Super Grip <i>Reduced weight, small profile. Not with 18" or 20" rear wheels, short tabs or Spinergy Carbon Blade rear wheels.</i> | \$ |
| TWRIM19 | The Natural-Fit® LT, Without Thumb Grip <i>Reduced weight, small profile.</i> | \$ |
| TWRIM22 | Omit Handrims | No charge |

Handrim Tab Length

| | | |
|--------|---|----------|
| TWRIM7 | Long Tabs | Standard |
| TWRIM8 | Short Tabs <i>Not with Mag rear wheels or The Natural-Fit® handrims with thumb grip.</i> | \$ |

Wheel Locks

| | | |
|---------|---|-----------|
| TWWLK7 | Compact Composite Push to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i> | No charge |
| TWWLK8 | Compact Composite Pull to Lock <i>Available with 0.75"-1.5" rear wheel spacing only.</i> | No charge |
| TWWLK10 | Extension Handles for Compact Composite Push/Pull to Lock | \$ |
| TWWLK11 | MBL Composite Scissor Lock | \$ |
| TWWLK3 | Uni-Lock | \$ |
| | Push Pull | |
| TWWLK4 | Extension Handles for Uni-Lock | \$ |

Side Guards

| | | |
|---------|---|----|
| TWSDG2 | Aluminum Rigid Removable | \$ |
| | Extra Small Small Standard | |
| TWSDG3 | Aluminum Quad Removable, with Hole | \$ |
| | Small Standard | |
| TWSDG7 | Carbon Rigid Removable | \$ |
| | Small Standard | |
| TWSDG8 | Carbon Quad Removable, with Hole | \$ |
| | Small Standard | |
| TWSDG12 | Carbon Fendered Removable | \$ |
| TWSDG13 | Side Guards Custom Cut to Profile of Rear Wheels <i>Available on aluminum and carbon fiber non-fendered side guards.</i> | \$ |
| CMD007 | Custom Offset Profiled Side Guards (Removable or Fold-Down) <i>Please specify desired offset to nearest 1/2" and side guards type in the Notes section of this order form. Side guard material must be selected. Some limitations may apply.</i> | \$ |

Armrests

| | | |
|--------|---|----|
| TWARM1 | Swing Away Tubular | \$ |
| TWARM4 | Flip Back Tubular <i>Not with short seat back height or velcro adjustable back upholstery.</i> | \$ |
| TWARM2 | Removable Height Adjustable Desk Arm with Rigid Side Guard <i>10" armpad. Requires 10" minimum seat depth.</i> | \$ |
| TWARM5 | Armrests Covered in Black Ultrasuede <i>Not with desk arms.</i> | \$ |

Positioning

Upper Body Positioning

| | | | |
|-----------|---|-------------------------------|----|
| SH295x-B5 | Bodypoint® PivotFit™ Dynamic Shoulder Harness | Cam buckle mounts, rear pull. | \$ |
| | Extra Small | Small | |
| SH120x-B3 | Bodypoint® Monoflex™ Chest Support | | \$ |
| | Extra Small | Small | |

Pelvic Positioning

| | | | |
|--------------|--|--|----|
| TWBLT1 | Velcro Adjustable Belt, Black Nylon Webbing | | \$ |
| HB275-S38-B1 | Bodypoint® Two Point Centre Pull Swivel Buckle, Small | | \$ |
| HB235-S38-B1 | Bodypoint® Two Point Centre Pull XS Push Button, Small | | \$ |
| EB275-S38 | Bodypoint® Evoflex Belt, Swivel Buckle, Small | | \$ |
| EB235-S38 | Bodypoint® Evoflex Belt, Centre Pull XS Push Button, Small | | \$ |
| TWBLT8 | Q'Straint® WC-19 Wheelchair Anchored Pelvic Belt | Recommended for use with transit tie down. | \$ |
| HW320-25-2 | Bodypoint® Optional Band Clamp Mounting Hardware | | \$ |

Lower Body Positioning

| | | | |
|--------|--------------------------|-------|----|
| FT240x | Bodypoint® Ankle Huggers | | \$ |
| | Xtra Small | Small | |

Anti-Tips

| | | | |
|--------|---|---|----|
| TWTIP1 | Rear Aluminum, Black | Seat widths 8"-10", 9"-11" and 10"-12" will ship with Centre Mount Single Non-Swivel Caster Anti-Tip at no additional cost. | \$ |
| TWTIP2 | User Friendly Flip-Up | Requires 12"-14" or 13"-15" seat widths. | \$ |
| TWTIP5 | Center Mount Swivel Caster Anti-Tip | | \$ |
| TWTIP6 | Center Mount Non-Swivel Caster Anti-Tip | | \$ |

Accessories

| | | | |
|---------|---|---------------------------|----|
| TWGRD1 | Clear Lexan Spoke Guards with TiMotion Logo | Not with Mag rear wheels. | \$ |
| TWBCK3 | Permobil Medical Necessities Backpack | | \$ |
| TWPCH3 | Permobil Medical Necessities Seat Pouch | | \$ |
| SP001 | Bodypoint® Mobility Bag | | \$ |
| TWUEST1 | Upper Extremity Support Tray | Requires armrests. | \$ |
| | 12"-14" seat widths | 13"-15" seat widths | |
| TWMPCT2 | Black Neoprene Impact Guards with Ultrasuede Stripe | | \$ |
| TWCH1 | Crutch Holder | | \$ |

SmartDrive Power Assist System

| | | | |
|-------------|--------------------------|--|----|
| MX2-850-xxx | Permobil SmartDrive MX2+ | Please provide SmartDrive MX2+ order form. Not available with 8" seat width, rear seat heights < 13.5" or 20" rear wheels. | \$ |
|-------------|--------------------------|--|----|

General

Permobil manufactures TiLite manual wheelchairs in a wide variety of configurations to meet the varied needs of wheelchair users. The final selection of a particular model, its options and accessories, and its measurements rests solely with the wheelchair user and the advising healthcare professional. These wheelchairs must be fitted by a qualified healthcare professional. Prices, standard features, options and accessories are subject to change without notice. All rights reserved.

Notes

TiFit Reference Diagrams

