ATTENTION!

PLEASE HANDLE WITH CARE

Remove both large wheels by pushing in button at center of wheel.

Ensure that **seat cushion has been removed** and given to owner of the chair.

When lifting and putting down wheelchair, be careful to **not damage frame**.

For questions, always consult the owner of the chair.

CONTACT AND WHEELCHAIR INFO

Contact Name:			
Contact phone number:			
Contact e-mail:			
Chair model:	Colour(s) of chair:	Colour(s) of chair:	
Chair measurements: Weight:	Height:	Width:	

Please refer to your owner's manual for additional safety tips and transportation guidelines.