

Power wheelchairs

M4 Nexa™

permobil


Compact without compromise

The M4 Nexa features a compact footprint and tight turning radius, with full power functions on seat widths as small as 12”. ActiveReach’s anterior tilt combined with a low 15.5” seat-to-floor height makes access and level transfers easy. Move through the doorways, hallways, and tight corners of everyday life with ease, and get closer to the people and moments that matter.

For users who everyone else asks to compromise.

- **Full power functions, including ActiveReach, at a 12” seat width**
- **15.5” seat-to-floor height with anterior tilt**
- **Narrow 21.5” width, with tight turning radius**






Seating designed to move with you

The Nexa Seating System is designed around the body’s natural points of articulation, significantly reducing displacement during power positioning. The new pivoting backrest attachment allows dynamic movement for consistent contact.

The result? Seating that feels secure and stays with you, so you can remain comfortable through every position change throughout the day.

- **Full seat-depth access preserved across seating positions**
- **Validated through extensive testing, including pressure mapping and repeated tilt/recline cycles**





A more precise fit, from the start


The broadest seat depth range in the category, 10” to 22”, in a single seating system. This gives clinicians the flexibility to configure the chair around the individual, rather than asking users to adapt to the chair.

- **Seat depth range, 10” - 22”, in one seating system**
- **Configure around the individual’s dimensions and positioning needs**
- **Grow into your chair, not out of it - adjust as a user’s needs change, without a completely new seating platform**

Work smarter

Less time servicing means more time supporting users with simpler access, adjustment, and configuration, plus the connected confidence of +Connect, QuickConfig, Fleet Management, and MyPermobil. So users can stay focused on their day, not their equipment.

- **Always-on connectivity for battery and range visibility**
- **Fault alerts before the chair stops, not after**
- **Usage history to support a user’s evolving needs**



Shroud colours

Specifications

Maximum user weight	136 kg
Maximum speed	10 km/h
Range	32.2 km*
Base width	21.5” (54.3 cm)
Base length	33.5” (85 cm)
Base turning radius	21” (53.3 cm)
Ground clearance	3”
Recommended battery type	52 Ah | 68 Ah
Weight including batteries	144 kg (52 Ah) 160.5 kg (68 Ah)
Drive electronics	R-net Power Platform 120A
Seat-to-floor height	15.5”, 17”

Seat elevation	Optional 12”
Drive wheel	12.5”
Castors	6” (front) | 7” (rear)
Tilt options	Posterior: 0° – 50° ActiveReach: 10°
Legrest angle - power	85° – 175°
Backrest angle - power recline	90° – 150°
Seat width	12” – 22”
Seat depth	10” - 22”

** Your actual range will depend on battery size, battery condition, and driving environment. | † Measured without cushion.*
Specifications may vary based upon configuration. Please refer to your owners manual for warnings and instructions for the safe operation of your wheelchair in a variety of driving conditions. Chair weights are approximate and include base, batteries, power seat functions, and seating system without accessories.

+Connect. More control. Less complexity.

Connecting users and service technicians through one intuitive platform improves communication and ensures you are better prepared for the unexpected. Power Platform and +Connect form the core digital infrastructure of our power wheelchairs, combining advanced assistive technology with smart connectivity.

Connected from the factory, each chair is online and ready for proactive service, remote insights, and efficient support right from the start. Every connected chair reflects Permobil’s commitment to reliability, transparency, and confidence throughout the life of the chair.



Learn more

Available on the App Store | Google Play

App Store is a trademark of Apple Inc. Google Play is a trademark of Google LLC.

permobil

Permobil Australia | 1300 845 483 | orders.au@permobil.com | www.permobil.com/en-au