## Mid-wheel









| M4 Nexa | M5 Corpus | M Corpus VS | M300 Corpus HD |
|---|---|---|---|
| 136 kg | 150 kg | 120 kg | 204 kg |
| 10 km/h | 10 km/h | 10 km/h | 8 km/h |
| 32.2 km* | 32.2 km* | 32.2 km* | 32.2 km* |
| 21.5” | 25” | 25” | 26” |
| 33.5” | 37.5” | 36.5” | 38” |
| 21” | 21.5” | 22” | 22” |
| 3” | 3” | 3” | 3” |
| 144 kg (52 Ah) / 160.5 kg (68 Ah) | 192 kg (72 Ah) | 193 kg (72 Ah) | 203 kg (72 Ah) |
| 15.5”, 17” | 17.5”, 18.5”, 19.5” | 17.5”, 18.5”, 19.5” | 18”, 19”, 20” (elevator)<br>18”, 19”, 20” (fixed) |
| *Optional* 12” | *Optional* 14” | 12” | *Optional* 8” |
| Posterior: 0° – 50°<br>ActiveReach: 10° | Posterior: 0° – 50°<br>ActiveReach: 5°, 10° & 20° | Posterior: 0° – 50°<br>Standing | Posterior: 0° – 45° |
| 90° – 150° (power) | 85° – 120° (manual)<br>85° – 180° (power) | 85° – 180° (power) | 85° – 120° (manual)<br>85° – 150° (power) |
| 85° – 175° (power) | 85° – 180° (manual & power) | 85° – 180° (power) | 85° – 175° (manual & power) |
| 10” – 22” | 14” – 22” (by 1” incr.) | 14” – 22” (by 1” incr.) | 19” – 25” (by 1” incr.) |
| 12” - 22” | 17”, 19”, 21”, 23” | 17”, 19”, 21”, 23” | 19”, 21”, 23”, 25” |
| 12” - 22” ** | 20”, 23” – 28” (by 1” incr.) | 23” – 28” (by 1” incr.) | 21”, 24” – 29” (by 1” incr.) |

***Physical back shell size. Overall backrest height is easily adjustable.*



*M4 Nexa*



Permobil power wheelchairs equipped with Power Platform and +Connect deliver advanced assistive technology and smart connectivity from the start.

Connected from the factory, these chairs are online and ready for proactive service, remote insights, and efficient support from day one. This connectivity reflects Permobil’s commitment to reliability, transparency, and confidence over the life of every chair.

Available on the App Store | ANDROID APP ON Google Play

*App Store is a trademark of Apple Inc. Google Play is a trademark of Google LLC.*

permobil

Permobil Australia | 1300 845 483 | info.au@permobil.com | www.permobil.com/en-au

BR-AU-PW004.1

Power wheelchairs

# Wheelchair matrix




permobil

| | Front-wheel | | | Mid-wheel | |
|---|---|---|---|---|---|
| |  |  |  |  |  |
| | F3 Corpus | F5 Corpus | F5 Corpus VS | M1 | M3 |
| **Max user weight** | 150 kg | 150 kg | 136 kg | 136 kg | 150 kg |
| **Max speed** | 10 km/h<br>8 km/h *without support wheels* | 10 km/h | 10 km/h | 8 km/h | 10 km/h |
| **Range** | 32.2 km* | 32.2 km* | 32.2 km* | 25.8 km* | 32.2 km* |
| **Base width** | 24” | 25.5” | 25.5” | 24” | 24” |
| **Base length** | 40” *with support wheels*<br>36” *without support wheels* | 42” *with support wheels*<br>38” *without support wheels* | 43” | 35.5” | 35.5” |
| **Base turning radius** | 26.5” | 30” | 30” | 20.5” | 20.5” |
| **Ground clearance** | 3” | 3” | 3” | 3” | 3” |
| **Weight w/ batteries** | 185 kg (72 Ah) | 187 kg (72 Ah) | 196 kg (72 Ah) | 145 kg (60 Ah) | 184 kg (72 Ah) |
| **Seat-to-floor height** | 17.5”, 18.5”, 19.5” | 17.5”, 18.5”, 19.5” | 17.5”, 18.5”, 19.5” | 16”, 17”, 18” | 17.5”, 18.5”, 19.5” |
| **Seat elevator** | *Optional* 12” | *Optional* 14” | 14” | — | *Optional* 12” |
| **Power tilt options** | Posterior: 0° – 50°<br>ActiveReach: 5°, 10°, 20° & 30° | Posterior: 0° – 50°<br>ActiveReach: 5°, 10°, 20° & 45° | Posterior: 0° – 50°<br>Standing | Posterior: 0° – 50° | Posterior: 0° – 50°<br>ActiveReach: 5°, 10° & 20° |
| **Recline options** | 85° – 120° (manual)<br>85° – 180° (power) | 85° – 120° (manual)<br>85° – 180° (power) | 85° – 180° (power) | — | 85° – 120° (manual)<br>85° – 180° (power) |
| **Legrest options** | 85° – 180° (manual & power) | 85° – 180° (manual & power) | 85° – 180° (power) | 85° – 140° (manual)<br>65° – 80° (adjust. swing-away) | 85° – 180° (manual & power) |
| **Seat depth** | 14” – 22” (by 1” incr.) | 14” – 22” (by 1” incr.) | 14” – 22” (by 1” incr.) | 14” – 22” (by 1” incr.) | 14” – 22” (by 1” incr.) |
| **Seat width** | 17”, 19”, 21”, 23” | 17”, 19”, 21”, 23” | 17”, 19”, 21”, 23” | 14” – 22” (by 2” incr.) | 17”, 19”, 21”, 23” |
| **Backrest height** | 20”, 23” – 28” (by 1” incr.) | 20”, 23” – 28” (by 1” incr.) | 23” – 28” (by 1” incr.) | 18” – 29” | 20”, 23” – 28” (by 1” incr.) |

*All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories. | * Your actual range will depend on battery size, battery condition, and driving environment.*



*Explorer Mini*

| | Front-wheel | | Mid-wheel | Rear-wheel | |
|---|---|---|---|---|---|
| |  |  |  |  |  |
| Paediatric Series | Explorer Mini | Koala | K300 PS Jr. | M300 PS Jr. | K450 MX |
| **Intended age range** | 12 months – 3 years | — | — | — | — |
| **Max user weight** | 16 kg | 35 kg | 75 kg | 75 kg | 60 kg |
| **Max speed** | 2.4 km/h | 6.5 km/h | 8 km/h | 10 km/h | 8 km/h |
| **Range** | 5.6 km* | 24.1 km* | 20.9 km* | 20.9 km* | 27.4 km* |
| **Base width** | 19” | 23” | 24” | 24” | 25” |
| **Base length** | 25” | 32” | 40” | 36.5” | 35” |
| **Min turning radius** | — | 21” | 25” | 20” | 28” |
| **Ground clearance** | 1” | 2.5” | 3” | 3” | 2.5” |
| **Weight w/ batteries** | 23.6 kg | 102 kg (40 Ah) | 123 kg (60 Ah) | 142 kg (60 Ah) | 153 kg (60 Ah) |
| **Seat-to-floor height** | 12.1” – 21.7” | 15” (elevator) | 16.5” (elevator)<br>16.5”, 17.25”, 18.5” (fixed) | 17” (elevator)<br>17”, 17.75”, 18.5” (fixed) | 3” – 26” (elevator) |
| **Seat elevator** | — | Up to 8” | *Optional* 8” | *Optional* 8” | 8” |
| **Power tilt options** | — | 0° – 25° | 0° – 45° | 0° – 45° | 0° – 45° |
| **Recline options** | — | 85° – 120° (manual) | 85° – 120° (manual) | 85° – 120° (manual) | 85° – 105° (manual) |
| **Legrest options** | — | 90° – 150° (manual) | 70° (swing-away)<br>90 – 150° (centre mount) | 70° or 90° (swing-away)<br>90° – 150° (centre mount) | 90° – 125° (manual) |
| **Seat depth** | — | 10”, 12”, 14” | 10” – 18” | 10” – 18” | 10” – 18” (by 1” incr.) |
| **Seat width** | — | 10”, 12”, 14” | 11” – 16” (by 1” incr.) | 11” – 16” (by 1” incr.) | 10”, 12”, 14”, 16” |
| **Backrest height** | **Total chair height**<br>29” – 37” | 14”, 16”, 20” | 16.5” – 22.5” (by 1” incr.) | 16.5” – 22.5” (by 1” incr.) | 14”, 16”, 20” (MX)<br>16.5” – 22.5” (PS Jr) |

*All specifications may vary based upon configuration. | Chair weights are approximate and include base, batteries, power seat functions and seating system without accessories.*
**Your actual range will depend on battery size, battery condition, and driving environment..*